

Thompson Retractor

Uncompromised Exposure™



Product Catalog

15TH EDITION

THIS PAGE INTENTIONALLY LEFT BLANK

Table of Contents

New Features	4
Thompson Retractor Frame Components.....	5-23
Thompson Retractor Blades	24-59
Unique System Components.....	60-73
Product Index	74-87
Product Measurement Guide	88-89
Service and Warranty.....	90

NEW Thompson Retractor Features

Faster, Lighter, More Efficient

Thompson's Elite Rail Clamp, unilateral crossbar, and retractor handles have all been updated for improved performance and ease of use.



Thompson is dedicated to providing the most forward thinking, proven, and high quality features.



Elite III Rail Clamps

Sterile, quick tightening knob locks rail clamp to bed securely in half the time, when compared with the Elite II Rail Clamp, while lightweight base design reduces overall weight by 1 pound.



Quick Angle Handles

Versatile, low-profile positioning of blades helps maintain a low-profile operating field and achieve better exposure. Blades can be toed in or out 45° with a T-Handle.



Improved Cam Joints

Now easier to lock and unlock, Cam II Clip-on Handles have been upgraded to include comfort "pinch" joints for easier removal from the frame.




NOTICE:

Elite II Rail Clamps and Angling Handles still available, but may require a longer lead time. Please call for more information.

Retractor Frame Components

Rail Clamps.....	6-7
Rail Attachments.....	7
Crossbars.....	8
Frame Arms.....	9
Lower Abdominal Bars.....	9
Hinged Ring.....	9
Extension Arms.....	10
Straight Arms.....	10
90° Rods.....	11
Upright Arms.....	11
Big-Big Joints.....	12-13
Big-Little Joints.....	14-15
Retractor Handles + Handle Accessories.....	16-19
Illumination.....	20
Accessories.....	21
Instrument Trays.....	22-23



ELITE III RAIL CLAMPS (CURRENT)

REF	ITEM DESCRIPTION	PART #
A	Elite III Rail Clamp 11" (25cm) with 1 Cam Joint	43902ACS
	Elite III Rail Clamp 16" (41cm) with 1 Cam Joint	43902ACL
B	Elite III Rail Clamp 18" (45cm) with 1 Cam Joint	43902AC
	Elite III Rail Clamp 11" (25cm) with 1 (3/8") Cam Joint	43902S
	Elite III Rail Clamp 18" (45cm) with 2 Cam Joints	43903AC
	Elite III Rail Clamp 16" (41cm) with 2 Cam Joints	43903ACL
	Elite III Rail Clamp 22" (55cm) with 1 Cam Joint	43904AC
	Elite III Rail Clamp 22" (55cm) with 2 Cam Joints	43905AC
	Elite III Rail Clamp 26" (65cm) with 1 Cam Joint	43906AC
	Elite III Rail Clamp 26" (65cm) with 2 Cam Joints	43907AC

OLDER GENERATION ELITE RAIL CLAMPS

REF	ITEM DESCRIPTION	PART #
	Elite II Rail Clamp 11" (25cm) with 1 Cam Joint	41902ACS
C	Elite II Rail Clamp 16" (41cm) with 1 Cam Joint	41902ACL
	Elite II Rail Clamp 18" (45cm) with 1 Cam Joint	41902AC
	Elite II Rail Clamp 11" (25cm) with one 3/8" Cam Joint	41902S
	Elite II Rail Clamp 18" (45cm) with 2 Cam Joints	41903AC
	Elite II Rail Clamp 16" (41cm) with 2 Cam Joints	41903ACL
	Elite II Rail Clamp 22" (55cm) with 1 Cam Joint	41904AC
	Elite II Rail Clamp 22" (55cm) with 2 Cam Joints	41905AC
	Elite II Rail Clamp 26" (65cm) with 1 Cam Joint	41906AC
	Elite II Rail Clamp 26" (65cm) with 2 Cam Joints	41907AC




INFINITE HEIGHT RAIL CLAMPS (IHRC) (CURRENT)

REF	ITEM DESCRIPTION	PART #
D	IHRC with 1 Cam Joint 11" - 22"	43902IRC
	IHRC with 2 Cam Joints 11" - 22"	43903IRC

OLDER GENERATION IHRC'S

REF	ITEM DESCRIPTION	PART #
E	IHRC with 1 Cam Joint 11" - 22"	41902IRC
	IHRC with 2 Cam Joints 11" - 22"	41903IRC




OLDER GENERATION RAIL CLAMPS

REF	ITEM DESCRIPTION	PART #
F	Power Rail Clamp	41950




TABLE ADAPTER (CURRENT)

REF	ITEM DESCRIPTION	PART #
G	Jackson Spine Frame Adapter	41927



RAIL EXTENDERS (CURRENT)

REF	ITEM DESCRIPTION	PART #
H	Rail Extender 15" Single Clamp	5844
	Rail Extender 15" Single Clamp (reverse of 5844)	5844B
	Rail Extender 20" Single Clamp	41938
I	Rail Extender 14" Long with 2 ¼" Offset Dual Clamp	41917
	Rail Extender 22" Dual Clamp	41929




UNILATERAL CROSSBARS (CURRENT)

REF	ITEM DESCRIPTION	PART #
A	Crossbar with 2 Cam III Joints 24" (11" x 13")	43900BC
	Crossbar with 2 Cam III Joints 24" (8" x 16")	43900BCW
B	Crossbar with 2 Cam III Joints (7" x 19 1/2") (in Hip System)	43990
C	Crossbar with 1 Cam III Joint Mini-Incision (in Kidney System)	43800

OLDER GENERATION UNILATERAL CROSSBARS

REF	ITEM DESCRIPTION	PART #
D	Crossbar with 2 Cam Joints 24" (11" x 13")	41900BC
	Crossbar with 2 Cam Joints 24" (8" x 16")	41900BCW
	Crossbar with 2 Adjustable Cam Joints 24"	41900BAC
	Crossbar with 2 Cam Joints (7" x 19 1/2") (in Hip System)	41990
	Crossbar with 1 Cam Joint Mini-Incision (in Kidney System)	44800




BILATERAL CROSSBARS (CURRENT)

REF	ITEM DESCRIPTION	PART #
E	Bilateral Crossbar 7 1/4" x 13 1/2" x 7 1/4"	41900F
F	34" Bilateral Crossbar Hinged 11 1/2" x 8 1/2" x 11 1/2"	41900PH
	41 1/2" Bilateral Crossbar Hinged 14 1/2" x 10" x 14 1/2"	41900PHL

OLDER GENERATION BILATERAL CROSSBARS

REF	ITEM DESCRIPTION	PART #
G	34" Bilateral Crossbar Hinged 11" x 8 1/2" x 11"	41900H
	41 1/2" Bilateral Crossbar Hinged 14" x 10" x 14"	41900HL
	Bilateral Crossbar with 2 Joints	41900FC




ANGLED LATERAL ARMS (CURRENT)

REF	ITEM DESCRIPTION	PART #
H	26" (6" x 20") (15cm x 50cm)	44126
I	24" (11" x 13") (28cm x 33cm)	44124W
J	24" (8" x 16") (20cm x 40cm)	44124N
	20" (10" x 10") (25cm x 25cm)	44120
K	18" (8" x 10") (20cm x 25cm)	44118N
	17" (7" x 10") (25cm x 18cm)	44117N
L	12" (6" x 6") (15cm x 15cm)	44014S
	19" (12" x 7") 3/8" Rod	44119
M	17" (7" x 10") 3/8" Rod	44117
	Double Angled Arm 22" (7" x 9" x 7")	44122DB
N	Double Angled Arm 28" (7 1/4" x 13 1/2" x 7 1/4")	44128DB


OLDER GENERATION ANGLED LATERAL ARMS

REF	ITEM DESCRIPTION	PART #
	12" (30cm)	44014
	12" (30cm)	44014N
	24" (61cm)	44124




LOWER ABDOMINAL BARS (CURRENT)

REF	ITEM DESCRIPTION	PART #
O	Lower Abdominal Bar 14" x 6"	44112SP
	Lower Abdominal Bar 10" x 8"	44111SP
	Lower Abdominal Bar 10" x 6"	44110SP




HINGED RING (CURRENT)

REF	ITEM DESCRIPTION	PART #
P	Ring Half Hinged 14 1/2" (2 required for use.)	44218
Q	Ring Extension 3" (2 required for use.)	44220
R	Extension Arm 14" with Cam Latch Mount	44614LC



EXTENSION ARMS (CURRENT)

REF	ITEM DESCRIPTION	PART #
A	Extension Arm 7 1/2" x 8 1/2" x 7 1/2" (in Kidney System)	44200
B	Extension Arm 12" with Cam Latch Joint	44014LC
C	Extension Arm 12" with Hinge Joint	44014HS
D	3/8" Extension Arm 12" with Hinge Joint	44014H
E	Hinged Extension Arm 12" x 12" (30cm x 30cm)	44115
	Angled Closure Arm 9" (4 1/2"x 4 1/2")	44109B
F	Angled Closure Arm 13" (6 1/2"x 6 1/2")	44113B
	Angled Closure Arm 18" (9"x 9")	44118B
	Angled Closure Arm 21" (10 1/2"x 10 1/2")	44121B
G	Extension Arm 14" with Cam Latch Mount	44614LC
H	Extension Arm 14" with Clip Joint	44014C
	Extension Arm 16" with Clip Joint	44016C
	Extension Arm 20" with Clip Joint	44020C



STRAIGHT ARMS (CURRENT)

REF	ITEM DESCRIPTION	PART #
I	Straight Arm 12" (30cm)	44012N
J	Straight Arm 20" (50cm)	44120N
	Straight Arm with 1 Cam Joint 18"	44118C
K	Straight Arm with 1 Cam Joint 20" (50cm)	44120C



90° ROD (CURRENT) *

REF	ITEM DESCRIPTION	PART #
L	90° Collapsible Rod with Grooves	44734




OLDER GENERATION UPRIGHT ARMS

REF	ITEM DESCRIPTION	PART #
M	Upright Arm 20"	44646C
	Upright Arm with Universal Joint 20"	44646N

OLDER GENERATION 90° RODS

REF	ITEM DESCRIPTION	PART #
	90° Collapsible Rod with 1 Cam Joint 14" x 19"	44733C
	90° Collapsible Rod with 1 Cam Joint 17" x 19"	44733AC
	90° Collapsible Rod with 1 Cam Joint 21" x 21"	44733BC
	90° Collapsible Rod 14" x 19"	44733
	90° Collapsible Rod 17" x 19"	44733A
	90° Collapsible Rod 21" x 21"	44733B
	90° Rigid Rod 14" x 19"	40733
	90° Rigid Rod 17" x 19"	40733A
	90° Rigid Rod 21" x 21"	40733B
	90° Rigid Rod 10" x 20"	40732
	Extension Arm 90° 10" x 6"	44643

* Current only as part of the Body Contouring System.



CAM III SERRATED SLIDE-ON JOINT (CURRENT)

ITEM DESCRIPTION	PART #
Cam III Serrated Slide-on Joint ½" x ½" *	43112C




CAM I SLIDE-ON JOINT (CURRENT)

ITEM DESCRIPTION	PART #
Cam I Slide-on Joint ½" x ½"	42110C



OLDER GENERATION CAM II SLIDE-ON JOINT

ITEM DESCRIPTION	PART #
Cam II Slide-on Joint ½" x ½"	42110S



OLDER GEN. CAM I SERRATED SLIDE-ON JOINT

ITEM DESCRIPTION	PART #
Cam I Serrated Slide-on Joint ½" x ½" *	42112C



OLDER GENERATION LATCH CAM JOINT

ITEM DESCRIPTION	PART #
Latch Cam Joint ½" x ½"	42111LC



OLDER GENERATION SLIDE-ON JOINT

ITEM DESCRIPTION	PART #
Slide-on Joint ½" x ½"	42110N

* Joint is too wide to slide over elbows on angled arms.



OLDER GENERATION LATCH JOINT

ITEM DESCRIPTION	PART #
Latch Joint 1/2" x 1/2"	42111LN



OLDER GENERATION ONE PIECE UNIVERSAL JOINT

ITEM DESCRIPTION	PART #
Hinged 1/2" x 1/2"	42111H



OLDER GENERATION SERRATED UNIVERSAL JOINT

ITEM DESCRIPTION	PART #
Slide-on 1/2" x 1/2" *	42112N



CAM II CLIP-ON JOINT (CURRENT)

ITEM DESCRIPTION	PART #
Cam II Low Profile Clip-On Joint ½" x ¼"	42126C



UNIVERSAL CLIP-ON JOINT (CURRENT)

ITEM DESCRIPTION	PART #
Clip-On Joint ½" x ¼"	42126NCP



UNIVERSAL HINGE JOINT (CURRENT)

ITEM DESCRIPTION	PART #
One Piece Universal Hinge Joint ½" x ¼"	42126H



UNIVERSAL SLIDE-ON JOINT (CURRENT)

ITEM DESCRIPTION	PART #
Slide-On Joint ½" x ¼"	42125N



UNIVERSAL LATCH JOINT (CURRENT)

ITEM DESCRIPTION	PART #
Screw Down Universal Latch Joint ½" x ¼"	42126LN

(* Joint is too wide to slide over elbows on angled arms)



UNIVERSAL SLIDE-ON JOINT (CURRENT)

ITEM DESCRIPTION	PART #
Universal Slide-On Joint Reversed ¼" x ½"	42125NM



OLDER GENERATION CAM I JOINTS

ITEM DESCRIPTION	PART #
Cam I Slide-On Joint 1/2" x 1/4"	42125C
Cam I Joint 1/2" x 3/8"	42125OC
Cam I Joint 1/4" x 1/4"	42124C




OLDER GENERATION CAM LATCH JOINT

ITEM DESCRIPTION	PART #
Cam Latch Clip-On Joint 1/2" x 1/4"	42126LC



OLDER GENERATION CAM II SLIDE-ON JOINT

ITEM DESCRIPTION	PART #
Cam II Low Profile Slide-On Joint 1/2" x 1/4"	42125S

**MICRO-ADJUSTABLE II HANDLES (CURRENT)**

REF	ITEM DESCRIPTION	S-LOCK PART #	INTERCH. PART #
A	Clip-on Quick Angle 10"	SL45006CGP	45006CGP
	Clip-on Quick Angle 7 ¾"	SL45007CGP	45007CGP
	Clip-on Quick Angle 15"	SL45006CGRP	45006CGRP
	Clip-on Quick Angle 20"	SL45007PGP	45007PGP
B	Clip-on 10 ¾" (non-angling)	SL45006STCP	45006STCP
	Clip-on 8 ¾" (non-angling)	SL45007STCP	45007STCP
	Slide-on Quick Angle 7 ¾"	SL45007STG	45007STG
	Clip-On 10 ¾" Bariatric (non-angling)	—	45006BTCP
	Clip-On 8 ¾" Bariatric (non-angling)	—	45007BTCP

OLDER GENERATION MICRO-ADJUSTABLE II HANDLES

REF	ITEM DESCRIPTION	S-LOCK PART #	INTERCH. PART #
C	Clip-On Angling 10"	SL45006CA	45006CA
	Clip-On Angling 7 ¾"	SL45007CA	45007CA
	Clip-On Angling 15"	SL45006CAR	45006CAR
	Clip-On Angling 20 ½"	SL45007P	45007P

OLDER GEN. MICRO-ADJ II HANDLES (Continued)

REF	ITEM DESCRIPTION	S-LOCK PART #	INTERCH. PART #
	Clip-On 10 ¾" (non-angling)	SL45006STC	45006STC
	Clip-On 8 ¾" (non-angling)	SL45007STC	45007STC
	Slide-On Angling 10"	SL45006STA	45006STA
	Slide-On Angling 7 ¾"	SL45007STA	45007STA
	Mini Micro Clip-On Angling 8" (cervical use only)	SL45010A	45010A
	Clip-On 15" (non-angling)	—	45006CR
	Clip-On 10 ¾" Bariatric (non-angling)	—	45006BTC
	Clip-On 8 ¾" Bariatric (non-angling)	—	45007BTC
	Slide-On Angling 8"	—	45009CS
	Slide-On 8 ¾" (non-angling)	—	45006ST
	Slide-On 10 ¾" (non-angling)	—	45006CS
	10 ¾" Bariatric (non-angling)	—	45006BS
	Slide-On 10" (non-angling)	—	45006LH
	Slide-On 8 ¾" Bariatric (non-angling)	—	45006BT
	Latch Interchangeable 10"	—	45006CL
	Latch Angling 10"	—	45009CL

OLDER GEN. MICRO-ADJ II ROTATABLE HANDLES

REF	ITEM DESCRIPTION	PART #
	Clip-On Rotatable 10"	45008STCP
D	Clip-On Rotatable 10"	45008STC
	Latch Rotatable 10"	45008C
	Slide-On Rotatable 10"	45008CH
	Slide-On Rotatable 10"	45008CS

OLDER GENERATION MICRO-ADJUSTABLE I HANDLES

ITEM DESCRIPTION	INTERCH. PART #
Universal Slide-on 8 ¾" (non-angling)	45007H
Universal Slide-on Angling 10"	45007AH
Slide-on 10"	45006LHB



CAM II CLIP-ON HANDLES (CURRENT)

REF	ITEM DESCRIPTION	S-LOCK PART #	INTERCH. PART #
E	Clip-on Quick Angle 8"	SL42126WGP	42126WGP
	Clip-on Quick Angle 10"	SL42126WGLP	42126WGLP
F	Clip-on 8" (non-angling)	SL42126WHP	42126WHP
	Clip-on Quick Angle 8" (3/8")	SL42127WGP	42127WGP

OLDER GENERATION CAM II HANDLES

REF	ITEM DESCRIPTION	S-LOCK PART #	INTERCH. PART #
G	Clip-On Angling 8"	SL42126WA	42126WA
	Clip-On Angling 10"	SL42126WAL	42126WAL
	3/8" Clip-On Angling 8"	SL42127WA	42127WA
	Clip-On 8" (non-angling)	SL42126WH	42126WH
	Slide-On Angling 8"	SL42125WA	42125WA
	Slide-On (non-angling) 8"	SL42125WH	42125WH
	Slide-On (non-angling) with 1/2" Handle Rod	—	42110WH



UNIVERSAL JOINTS + HANDLES (CURRENT)

REF	ITEM DESCRIPTION	S-LOCK PART #	INTERCH. PART #
H	Quick Angle Handle 8"	SL45003GS	45003GS
	Quick Angle Handle 10"	SL45003GSL	45003GSL
I	Handle 8" (non-angling)	SL45003L	45003L
	Handle 10" (non-angling)	SL45003LM	45003LM
J	Clip-on Joint 1/2" x 1/4"	42126NCP	42126NCP
K	Slide-on Joint 1/2" x 1/4"	42125N	42125N

OLDER GENERATION HANDLES

REF	ITEM DESCRIPTION	PART #
L	S-Lock Angling Handle 8"	SL45003AS
	S-Lock Angling Handle 10"	SL45003ASL
	Angling Handle 8"	45003AS
	Angling Handle 10"	45003ASL
	S-Lock Handle 8"	SL45003
	Interchangeable Blade Handle 8"	45003
	Angling Retractor Handle	45003A
	Rotatable Handle 6"	45000
M	Rotatable Handle 8"	45001
	Rotatable Handle 10"	45001L




ROTATABLE HANDLES NOTICE:

Replacement Rotatable Knob available. Call to order.

Handles for Non-Thompson Systems

The handles on this page allow for Thompson's features and blades to be used with other brands of retractor systems.



RATCHET RETRACTOR HANDLES (CURRENT)


	ITEM DESCRIPTION	S-LOCK PART #	INTERCH. PART #
A	Quick Angle 10 ¼"	SL45005GL	45005GL
B	10 ¼" (non-angling)	SL45005L	45005L

OLDER GENERATION RATCHET RETRACTOR HANDLES

	ITEM DESCRIPTION	PART #
C	S-Lock Angling 10 ¼"	SL45005AL
	Angling 10 ¼"	45005AL
D	Rotatable	45005

SynFrame® is a Registered Trademark of Depuy Synthes, Inc.® Thompson® Surgical Instruments, Inc. has no affiliation with the SynFrame® or Depuy Synthes, Inc.®

Omni-Tract® is a Registered Trademark of Integra®. Thompson® Surgical Instruments, Inc. has no affiliation with Omni-Tract® or Integra®.




THOMPSON ADAPTER HANDLES FOR SYNFRAME® (CURRENT)

	ITEM DESCRIPTION	S-LOCK PART #	INTERCH. PART #
E	Quick Angle 8 ¾"	SL45014GL	45014GL
F	8 ¾" (non-angling)	SL45014L	45014L

OLDER GENERATION THOMPSON ADAPTER HANDLES FOR SYNFRAME®

	ITEM DESCRIPTION	S-LOCK PART #	INTERCH. PART #
G	Angling 8 ¾"	SL45014AL	45014AL




THOMPSON ADAPTER HANDLES FOR OMNI-TRACT® (CURRENT) *

	ITEM DESCRIPTION	S-LOCK PART #	INTERCH. PART #
H	Quick Angle 8 ¾"	SL45004GL	45004GL
I	8 ¾" (non-angling)	SL45004L	45004L

OLDER GENERATION THOMPSON ADAPTER HANDLES FOR OMNI-TRACT®

	ITEM DESCRIPTION	PART #
J	S-Lock Angling 8 ¾"	SL45004AL
	Angling 8 ¾"	45004AL
	Rotatable	45004

* Also see standard S-Lock Clip-on Handles.



HAND HELD ADAPTERS (CURRENT)

REF	ITEM DESCRIPTION	PART #
K	Quick Angle Hand Held Adapter	SL42128G
L	Manipulator Handle	42228

OLDER GENERATION HAND HELD ADAPTERS

REF	ITEM DESCRIPTION	PART #
M	Angling Hand Held Adapter	SL42128F
N	Locking Cam Hand Held Adapter	42128C
O	Hand Held Adapter	42128



T-HANDLE (CURRENT)

REF	ITEM DESCRIPTION	PART #
P	T-Handle (used with Quick Angle Handles)	60020


OLDER GENERATION WRENCHES

REF	ITEM DESCRIPTION	PART #
Q	Angling Head Wrench (used with Angling Handles)	43032
R	Screw Down Joint Wrench (used with Screw Down Joints)	42127



FLEXIBLE BLADE HOLDER (CURRENT)

REF	ITEM DESCRIPTION	PART #
S	Flexible Blade Holder with Cam Joint (used to hold Nerve Root Retractors)	42140C




LED ILLUMINATION (CURRENT)

	ITEM DESCRIPTION	PART #
A	Mini Retractor Lite LED (cervical use only)	40002MRL
B	Mini Retractor Lite Clip (Radiopaque)	40026MRL
C	PLA MIS Retractor Lite LED (for use with MIS PLA® System only)	40002MIS
D	LED Light Source	LLS-2000-2




XENON ILLUMINATION (CURRENT)

	ITEM DESCRIPTION	PART #
E	Lite Wand Xe	40020X
	Lite Wand Xe with ACMI Cable	40001X
	Lite Wand Xe with Storz Cable	40002X
	Lite Wand Xe with Wolf Cable	40004X
	Lite Wand Xe with Olympus Cable	40019X
F	Retractor Lite Xe	40020RX
	Retractor Lite Xe with ACMI Cable	40001RX
	Retractor Lite Xe with Storz Cable	40002RX
	Retractor Lite Xe with Wolf Cable	40004RX
	Retractor Lite Xe with Olympus Cable	40019RX
G	Lite Clip (Radiopaque)	40026
H	Lite Wand II Olympus Adapter	40017
I	Lite Wand II Wolf Adapter	40018
J	Lite Wand II Storz Adapter	40019



DEPTH GAUGES (CURRENT)

	ITEM DESCRIPTION	PART #
K	Anterior Cervical	46521
L	Anterior Lumbar	51236




HYSTERECTOMY TRAY (CURRENT)

	ITEM DESCRIPTION	PART #
N	Hysterectomy Tray	41933




SALS (CURRENT)

	ITEM DESCRIPTION	PART #
M	Suction for Anterior Lumbar Surgery (SALS)	51234
	Replacement Stylet	51229



INSTRUMENT TRAYS (CURRENT)

REF	ITEM DESCRIPTION	PART #
A	22" x 11" x 3 ½" w/Pin Mat	50000G
	26" x 10" x 5"	50000BL




ORGANIZED INSTRUMENT TRAY - HPB (CURRENT)

REF	ITEM DESCRIPTION	PART #
B	Hepatobiliary	50000HPB




LITE / ACCESSORY INSTRUMENT TRAY (CURRENT)

REF	ITEM DESCRIPTION	PART #
C	6 ¼" x 14" x 2 ½" (plastic)	50000LW



ACCESSORY INSTRUMENT TRAY (CURRENT)

REF	ITEM DESCRIPTION	PART #
D	Accessory - 3 Compartments	50000A3




ORGANIZED INSTRUMENT TRAY -
BOLLING MITRAL VALVE (CURRENT)

REF	ITEM DESCRIPTION	PART #
E	Bolling	50000BOL



INSTRUMENT TRAY - DRAKE HOOK (CURRENT)

REF	ITEM DESCRIPTION	PART #
F	Drake Hook	50000DH





ORGANIZED INSTRUMENT TRAYS - ANTERIOR LUMBAR SPINE (CURRENT)

REF	ITEM DESCRIPTION	PART #
	AL Renal Vein Blades	50000ALR
	Anterior Lumbar	50000ALT
G	Expanded Reverse Lip ASA	50000ERL
	Reverse Lip ASA	50000RLT

ORGANIZED INSTRUMENT TRAYS - ANTERIOR CERVICAL SPINE (CURRENT)

REF	ITEM DESCRIPTION	PART #
	Anterior Cervical	50000ACT
I	Expanded Anterior Cervical	50000EAC
	Cervical Blade	50000UAC
	Cervical Blade Tray (Insert for 50000UAC)	50000UACT



ORGANIZED INSTRUMENT TRAYS - MIS PLA (CURRENT)

REF	ITEM DESCRIPTION	PART #
	Midline Blades	50000ML
H	MIS Spine Frame	50000MSF
	PLA	50000PLA
	PLA One Step Dilator	50000PLAO

ORGANIZED INSTRUMENT TRAYS - CCA (CURRENT)

REF	ITEM DESCRIPTION	PART #
J	CCA Blades	50000CCAB
	CCA Frame	50000CCAF

Blades

General and Vascular Surgery Blades.....25-45

Balfour 25-26

Harrington ("Sweetheart")..... 27

Weinberg 28

Kelly 28

Splanchnich 29

Fence 29

Malleable 30-32

Mini-Open 32-33

Richardson 33-35

St. Marks 35

Deaver 36

Renal Vein 37-38, 49

Fenestrated..... 38

Murphy Rake..... 39

Volkman 39

Aortic Aneurysm 40

Allison Lung 40

Double Angled 41

Heany..... 42

Bladder 42-43

Prostate..... 44

Rectal 45

Spine Surgery Blades46-59

Anterior Lumbar..... 46-49


Cervical 50-57

Posterior Lumbar 58-59



BLADE NOTICE:
Non-S-Lock blades may require a longer lead time.
Many blades also available in aluminum, call for information.

BALFOUR BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
[with lips] 32mm x 32mm (1 ¼" x 1 ¼")	SL46142	46142	45042	45142 *	43142
[with lips] 51mm x 54mm (2" x 2")	SL46142P	46142P	45042P	45142P	43142P
[with lips] 76mm x 51mm (3" x 2")	SL46140P	46140P	45040P	45140P	43140P
[with lips] 70mm x 73mm (2 ¾" x 3")	SO46140	46140	45040	45140	43140 *
[with lips] 83mm x 73mm (3 ¼" x 3")	SO46141	46141	45041	45141	43141 *
[with lips] 83mm x 123mm (3 ¼" x 5")	SO46141D	46141D	45041D	45141D	43141D
[with lips] 83mm x 175mm (3 ¼" x 6 7/8")	SO46141F	46141F	45041F	45141F	43141F
[small lips] 70mm x 73mm (2 ¾" x 3")	SO46140C	46140C	45040C	45140C	43140C
[small lips] 83mm x 73mm (3 ¼" x 3")	SO46141C	46141C	45041C	45141C	43141C
Radiolucent [with lips] 32mm x 32mm (1 ¼" x 1 ¼") [Anodized]	SL46142EB	46142EB	—	—	—
Radiolucent [with lips] 32mm x 32mm (1 ¼" x 1 ¼")	SL46142ET	—	—	—	—
Radiolucent [with lips] 51mm x 51mm (2" x 2") [Anodized]	SL46142PEB	46142PEB	—	—	—
Radiolucent [with lips] 65mm x 72mm (2 ½" x 2 ¾") ANTERIOR LUMBAR USE ONLY	SL46140ET	—	—	—	—
Radiolucent [with lips] 65mm x 72mm (2 ½" x 2 ¾")	SO46140ET	46140ET	45040E	—	—
Radiolucent [with lips] 65mm x 72mm (2 ½" x 2 ¾") [Anodized]	SO46140EB	46140EB	—	—	—
Radiolucent [small lips] 65mm x 72mm (2 ½" x 2 ¾")	SO46140CET	46140CET	45040CE	—	—
Radiolucent [small lips] 65mm x 72mm (2 ½" x 2 ¾") [Anodized]	SO46140CEB	46140CEB	—	—	—
Radiolucent [with lips] 83mm x 70mm (3 ¼" x 2 ¾") ANTERIOR LUMBAR USE ONLY	SL46141ET	—	—	—	—
Radiolucent [with lips] 83mm x 70mm (3 ¼" x 2 ¾")	SO46141ET	46141ET	45041E	—	—
Radiolucent [with lips] 83mm x 70mm (3 ¼" x 2 ¾") [Anodized]	SO46141EB	46141EB	—	—	—
Radiolucent [small lips] 83mm x 70mm (3 ¼" x 2 ¾")	SO46141CET	46141CET	—	—	—
Radiolucent [small lips] 83mm x 70mm (3 ¼" x 2 ¾") [Anodized]	SO46141CEB	46141CEB	—	—	—
Radiolucent [with lips] 83mm x 127mm (3 ¼" x 5")	SO46141DET	—	—	—	—
Radiolucent [with lips] 83mm x 127mm (3 ¼" x 5") [Anodized]	SO46141DEB	46141DEB	—	—	—

* Fixed / Swivel Blades with 6" handle available, add 'S' to front of part number.

**NOTICE:**

Any Balfour blades may be ordered with small lips or without lips. Call your Thompson representative to order.

WIDE BALFOUR (ROCHARD) BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
83mm x 57mm (3 1/4" x 2 1/4")	SO46241A	46241A	—	—	43241A
100mm x 61mm (4" x 2 3/8")	SO46241B	46241B	—	—	43241B
114mm x 63mm (4 1/2" x 2 1/2")	SO46241C	46241C	—	—	43241C
Radiolucent 83mm x 57mm (3 1/4" x 2 1/4") [Anodized Aluminum]	SO46241AEB	46241AEB	—	—	—
Radiolucent 102mm x 57mm (4" x 2 1/4") [Anodized Aluminum]	SO46241BEB	46241BEB	—	—	—
Radiolucent 114mm x 57mm (4 1/2" x 2 1/4") [Anodized Aluminum]	SO46241CEB	46241CEB	—	—	—

NON-SLIP BALFOUR BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Right 77mm x 94mm (3" x 3 3/4")	SO46169C	46169C	—	—	—
Left 77mm x 94mm (3" x 3 3/4")	SO46168C	46168C	—	—	—
Right 77mm x 66mm	SO46169B	46169B	—	—	—
Left 77mm x 66mm	SO46168B	46168B	—	—	—

HARRINGTON (SWEETHEART) BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
38mm x 127mm (1 ½" x 5")	SL46162P	46162P	P45062	45162P	43162P
64mm x 102mm (2 ½" x 4")	SL46165	46165	—	—	43165
64mm x 152mm (2 ½" x 6")	SL46160	46160	45060	45160	43160
64mm x 203mm (2 ½" x 8")	SL46163	46163	45063	45163	43163
64mm x 254mm (2 ½" x 10")	SL46164	46164	45064	45164	43164
MALLEABLE 64mm x 152mm (2 ½" x 6")	—	46160S	—	—	—
114mm x 203mm (4 ½" x 8")	—	46161M	—	—	—
114mm x 254mm (4 ½" x 10")	SL46161	46161	45061	45161	43161
38mm x 203mm (1 ½" x 8")	—	46166	—	—	—
Radiolucent 64mm x 152mm (2 ½" x 6") Titanium	SL46160ET	46160ET	—	—	—
Radiolucent 64mm x 152mm (2 ½" x 6") Titanium	SO46160ET	46160ET	—	—	—
Radiolucent 64mm x 152mm (2 ½" x 6") Aluminum	SL46160E	46160E	45060E	—	—
Radiolucent 64mm x 203mm (2 ½" x 8")	SL46163ET	—	—	—	—
Radiolucent 64mm x 254mm (2 ½" x 10")	SL46164ET	—	—	—	—
Radiolucent 64mm x 102mm (2 ½" x 4")	SL46165ET	—	—	—	—

WEINBERG BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
102mm x 152mm (4" x 6")	—	Obsolete **	—	—	42520 *
83mm x 133mm (3 ¼" x 5 ¼")	—	Obsolete **	—	—	42521 *
102mm x 203mm (4" x 8")	—	Obsolete **	—	—	42522 *

* ½" fixed handle requires ½" x ½" joint for use, see pages 12 - 13. (Joint #43112C recommended.)

** Interchangeable Weinberg blades obsolete due to possible handle damage when under extremely heavy retraction.

KELLY BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
38mm x 51mm (1 ½" x 2")	SL46127	46127	45027	45127 *	43127 *
51mm x 64mm (2" x 2 ½")	SL46128	46128	45028	45128	43128 *
64mm x 76mm (2 ½" x 3")	SO46129	46129	45029	45129	43129 *
76mm x 89mm (3" x 3 ½")	SO46130	46130	45030	45130	43130 *
Radiolucent 38mm x 51mm (1 ½" x 2")	SL46127ET	46127ET	—	—	—
Radiolucent 51mm x 64mm (2" x 2 ½")	SL46128ET	46128ET	—	—	—
Radiolucent 64mm x 76mm (2 ½" x 3")	SO46129ET	46129ET	—	—	—
Radiolucent 76mm x 89mm (3" x 3 ½")	SO46130ET	46130ET	—	—	—

* Fixed / Swivel Blades with 6" handle available, add 'S' to front of part number.

SPLANCHNICH BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
51mm x 114mm (2" x 4 1/2")	SL46145	46145	—	45145	—
51mm x 152mm (2" x 6")	SL46146	46146	—	45146	—
Radiolucent 25mm x 152mm (1" x 6")	SL46148ET	46148ET	45048E	—	—
Radiolucent 25mm x 203mm (1" x 8")	SL46149ET	46149ET	45049E	—	—
Radiolucent 51mm x 114mm (2" x 4 1/2")	SL46145ET	46145ET	45045E	—	—
Radiolucent 51mm x 152mm (2" x 6")	SL46146ET	46146ET	45046E	45145E	—
Radiolucent 51mm x 203mm (2" x 8")	SL46147ET	46147ET	45047E	—	—

FENCE BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
102mm x 102mm (4" x 4")	SO46134	46134	45034	45134	43134
102mm x 127mm (4" x 5")	SO46135	46135	45035	45135	43135
102mm x 152mm (4" x 6")	SO46136	46136	45036	45136	43136
102mm x 203mm (4" x 8")	SO46136A	46136A	45036A	45136A	43136A
Radiolucent 102mm x 102mm (4" x 4")	SO46134ET	46134ET	—	—	—
Radiolucent 102mm x 127mm (4" x 5")	SO46135ET	46135ET	—	—	—
Radiolucent 102mm x 152mm (4" x 6")	SO46136ET	46136ET	45036E	—	—
Radiolucent 102mm x 203mm (4" x 8")	SO46136AET	46136AET	—	—	—

MALLEABLE BLADES




ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
25mm x 110 mm (1" x 4 5/16")	SL45656A	45656A	—	—	—
13mm x 90mm (1/2" x 3 1/2")	SL45657A	45657A	—	—	—
13mm x 64mm (1/2" x 2 1/2")	SL46115	46115	45015	45115	—
13mm x 152mm (1/2" x 6")	SL46116	46116	45016	45116	—
13mm x 203mm (1/2" x 8")	SL46117	46117	45017	45117	—
25mm x 203mm (1" x 8")	SL46118	46118	45018	45118	43118
25mm x 104mm (1" x 4")	SL46118P	46118P	P45018	45118P	—
25mm x 254mm (1" x 10")	SL46118A	46118A	45018A	45118A	—
38mm x 254mm (1 1/2" x 10")	SL46118B	46118B	45018B	45118B	—
25mm x 127mm (1" x 5")	SL46119	46119	45019	45119	—
51mm x 102mm (2" x 4")	SL46122P	46122P	P45022	45122P	—
51mm x 203mm (2" x 8")	SL46122	46122	45022	45122	43122
51mm x 254mm (2" x 10")	SL46122A	46122A	45022A	45122A	43122A
38mm x 203mm (1 1/2" x 8")	SL46123	46123	45023	45123	—
64mm x 254mm (2 1/2" x 10")	SL46124	46124	45024	45124	43124 *
76mm x 254mm (3" x 10")	SO46124A	46124A	45024A	45124A	43124A
102mm x 254mm (4" x 10")	SO46124B	46124B	45024B	45124B	43124B
Radiolucent 25mm x 110mm (1" x 4 1/2")	SL45656AT	45656AT	—	—	—
Radiolucent 13mm x 90mm (1/2" x 3 1/2")	SL45657AT	45657AT	—	—	—
Radiolucent 13mm x 64mm (1/2" x 2 1/2")	SL46115ET	46115ET	—	—	—
Radiolucent 13mm x 152mm (1/2" x 6")	SL46116ET	46116ET	—	—	—
Radiolucent 13mm x 203mm (1/2" x 8")	SL46117ET	46117ET	—	—	—
Radiolucent 25mm x 203mm (1" x 8")	SL46118ET	46118ET	—	—	—
Radiolucent 25mm x 102mm (1" x 4")	SL46118PET	46118PET	—	—	—
Radiolucent 51mm x 254mm (2" x 10")	SL46122AET	46122AET	—	—	—

Continued on next page...

MALLEABLE BLADES (Continued)

ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent 51mm x 203mm (2" x 8")	SL46122ET	46122ET	45022E	—	—
Radiolucent 64mm x 254mm (2 ½" x 10")	SL46124ET	46124ET	—	—	—
Radiolucent 76mm x 254mm (3" x 10")	SO46124AET	46124AET	—	—	—
Radiolucent 102mm x 254mm (4" x 10")	SO46124BET	46124BET	—	—	—
38mm x 102mm (1 ½" x 4")	—	46120	45020	45120	—
102mm x 165mm (4" x 6 ½")	—	46114	—	45114	—
76mm x 165mm (3" x 6 ½")	—	46113	—	45113	—
51mm x 165mm (2" x 6 ½")	—	46112	—	45112	—
25mm x 152mm (1" x 6")	—	46119A	45019A	45119A	—
38mm x 152mm (1 ½" x 6")	—	46120B	45020B	45120B	—
38mm x 127mm (1 ½" x 5")	—	46120A	45020A	45120A	—
25mm x 76mm (1" x 3")	—	46118C	—	—	—


* Fixed Blade with 6" handle available, add 'S' to front of part number.

MALLEABLE FINGER BLADES

ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
76mm (3")	SL46191SP	46191SP	—	—	—
102mm (4")	SL46191P	46191P	—	—	—
152mm (6")	SL46191	46191	45091	45191	43191 *
203mm (8")	SL46191L	46191L	—	—	—
Radiolucent 152mm (6")	SL46191ET	46191ET	—	—	—

* Fixed Blade with 6" handle available, add 'S' to front of part number.

MALLEABLE SLOTTED BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
51mm x 203mm (2" x 8")	SL46122V	46122V	—	—	—
76mm x 203mm (3" x 8")	SO46124AV	46124AV	—	—	—
51mm x 165mm (2" x 6 ½")	—	46202	—	—	—
76mm x 165mm (3" x 6 ½")	—	46203	—	—	—

MINI-OPEN MALLEABLE BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
25mm x 75mm (1" x 3")	SL46220	46220	—	—	—
25mm x 100mm (1" x 4")	SL46221	46221	—	—	—
25mm x 125mm (1" x 5")	SL46222	46222	—	—	—

MINI-OPEN DOUBLE ANGLED BLADES




ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
51mm x 127mm (2" x 5")	SL46224	46224	—	—	—
51mm x 152mm (2" x 6")	SL46225	46225	—	—	—
51mm x 165mm (2" x 6 ½")	SL46226	46226	—	—	—
Deep 51mm x 127mm (2" x 5")	SL46227	46227	—	—	—

MINI-OPEN RICHARDSON BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
51mm x 76mm (2" x 3")	SL46155	46155	45055	—	—
51mm x 76mm (2" x 3")	SO46155	—	—	—	—
51mm x 102mm (2" x 4")	SL46156	46156	—	—	—
51mm x 102mm (2" x 4")	SO46156	—	—	—	—
51mm x 127mm (2" x 5")	SL46157	46157	—	—	—
51mm x 127mm (2" x 5")	SO46157	—	—	—	—
51mm x 152mm (2" x 6")	SL46158	46158	—	45158	—
51mm x 152mm (2" x 6")	SO46158	—	—	—	—

RICHARDSON BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
44mm x 51mm (1 ¾" x 2")	SL46149	46149	—	—	—
51mm x 76mm (2" x 3")	SL46154	46154	—	—	—
51mm x 102mm (2" x 4")	SL46150	46150	45050	45150	43150
51mm x 127mm (2" x 5")	SO46151	46151	45051	45151 *	43151 *
51mm x 152mm (2" x 6")	SO46152	46152	45052	45152	43152
51mm x 178mm (2" x 7")	SO46153	46153	45053	45153	—
Radiolucent 51mm x 76mm (2" x 3")	SL46154ET	—	—	—	—
Radiolucent 51mm x 102mm (2" x 4")	SL46150ET	—	—	—	—
Radiolucent 51mm x 127mm (2" x 5")	SO46151ET	—	—	—	—
Radiolucent 51mm x 152mm (2" x 6")	SO46152ET	—	—	—	—
Radiolucent 51mm x 178mm (2" x 7")	SO46153ET	46153ET	—	—	—

* Fixed / Swivel Blades with 6" handle available, add 'S' to front of part number.

MINI RICHARDSON BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
19mm x 25mm ($\frac{3}{4}$ " x 1")	SL46125	46125	—	45125 *	43125 *
19mm x 51mm ($\frac{3}{4}$ " x 2")	SL46125A	46125A	—	45125A *	43125A
25mm x 32mm (1" x 1 $\frac{1}{4}$ ")	SL46126	46126	—	45126 *	43126 *

* Fixed / Swivel Blades with 6" handle available, add 'S' to front of part number.
Older style Fixed Carotid Blades also available, call for more information.

ST. MARKS BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
[no lip] 64mm x 178mm (2 $\frac{1}{2}$ " x 7")	SO46187A	46187A	45087A	45187A	—
[with lip] 64mm x 178mm (2 $\frac{1}{2}$ " x 7")	SO46188	46188	45088	45188	43188
[with lip] 64mm x 140mm (2 $\frac{1}{2}$ " x 5 $\frac{1}{2}$ ")	SO46189	46189	45089	45189	—
Radiolucent [with lip] 64mm x 140mm (2 $\frac{1}{2}$ " x 5 $\frac{1}{2}$ ")	SO46189ET	46189ET	—	—	—
Radiolucent [with lip] 64mm x 178mm (2 $\frac{1}{2}$ " x 7")	SO46188ET	46188ET	—	—	—


DEAVER BLADES (NON-MALLEABLE)



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
25mm x 127mm (1" x 5")	SL46170	46170	45070	45170	43170 *
38mm x 127mm (1 ½" x 5")	SL46172	46172	45072	45172	43172 *
51mm x 127mm (2" x 5")	SL46174	46174	45074	45174	43174
64mm x 127mm (2 ½" x 5")	SO46175	46175	45075	45175	43175
76mm x 127mm (3" x 5")	SO46176	46176	45076	45176	43176 *
102mm x 127mm (4" x 5")	SO46178	46178	—	—	—
25mm x 102mm (1" x 4")	—	—	P45070	—	—
51mm x 102mm (2" x 4")	—	—	P45074	—	—
Radiolucent 25mm x 127mm (1" x 5")	SL46170ET	46170ET	—	—	—
Radiolucent 38mm x 127mm (1 ½" x 5")	SL46172ET	46172ET	—	—	—
Radiolucent 51mm x 127mm (2" x 5")	SL46174ET	46174ET	—	—	—
Radiolucent 64mm x 127mm (2 ½" x 5")	SO46175ET	46175ET	—	—	—
Radiolucent 76mm x 127mm (3" x 5")	SO46176ET	46176ET	45076E	—	—
Radiolucent 76mm x 127mm (3" x 5") ANTERIOR LUMBAR USE ONLY	SL46176ET	—	—	—	—
Radiolucent 102mm x 127mm (4" x 5")	SO46178ET	46178ET	—	—	—
Radiolucent 102mm x 127mm (4" x 5") ANTERIOR LUMBAR USE ONLY	SL46178ET	—	—	—	—

* Fixed Blades with 6" handle available, add 'S' to front of part number.

CONCAVE RENAL VEIN BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
25mm x 140mm (1" x 5 1/2") Rigid	SL46119B	46119B	45019B	45119B	—
25mm x 191mm (1" x 7 1/2") Rigid	SL46119C	46119C	45019C	—	43119C
Radiolucent 25mm x 140mm (1" x 5 1/2") Rigid	—	—	45019BET	—	—
Radiolucent 25mm x 140mm (1" x 5 1/2") Malleable	SL46119BET	46119BET	45019BE	—	—
Radiolucent 25mm x 190mm (1" x 7 1/2") Malleable	SL46119CET	46119CET	45019CE	—	—

RIGID RENAL VEIN BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent 25mm x 203mm (1" x 8")	SL46118RET	46118RET	45018RET	—	—
Radiolucent 25mm x 140mm (1" x 5 1/2")	SL46119RET	46119RET	45019RE	45119RE	—
Radiolucent 25mm x 191mm (1" x 7 1/2")	SL46120RET	46120RET	—	45120RE	—
Radiolucent 25mm x 178mm (1" x 7")	—	—	45020RE	—	—
Radiolucent 25mm x 254mm (1" x 10")	SL46121RET	46121RET	45021RE	45121RE	—
Radiolucent 25mm x 305mm (1" x 12")	SL46122RET	46122RET	45022RE	—	—
Radiolucent 25mm x 152mm (1" x 6")	—	46118CET	45018CE	—	—
Radiolucent 51mm x 152mm (2" x 6")	—	46118DET	—	—	—
Radiolucent 51mm x 203mm (2" x 8")	—	46118EET	—	—	—

Renal Vein Blades continued on next page, additional Renal Vein Blades shown with Anterior Lumbar Spine Blades.

MALLEABLE RENAL VEIN BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
25mm x 203mm (1" x 8")	SL46118R	46118R	45018R	45118R	43118R


Additional Renal Vein Blades shown on previous page, and with Anterior Lumbar Spine Blades.

FENESTRATED BLADES




ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
44mm x 89mm (1 ¾" x 3 ½")	SL46143	46143	—	45143	43143
48mm x 127mm (1 7⁄8" x 5")	SL46144	46144	—	45144	—

MURPHY RAKE BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Sharp - 3 Prong 19mm x 13mm ($\frac{3}{4}$ " x $\frac{1}{2}$ ")	SL46216	46216	—	45216	—
Sharp - 4 Prong 25mm x 13mm (1" x $\frac{1}{2}$ ")	SL46214	46214	—	45214	43214
Blunt - 3 Prong 19mm x 13mm ($\frac{3}{4}$ " x $\frac{1}{2}$ ")	SL46217	46217	—	45217	—
Blunt - 4 Prong 25mm x 13mm (1" x $\frac{1}{2}$ ")	SL46215	46215	—	45215	—

VOLKMAN BLADES




ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Sharp 25mm x 13mm (1" x $\frac{1}{2}$ ")	SL46210	46210	—	45210	—
Blunt 25mm x 13mm (1" x $\frac{1}{2}$ ")	SL46211	46211	—	45211	—
Sharp 16mm x 13mm ($\frac{5}{8}$ " x $\frac{1}{2}$ ")	SL46212	46212	—	45212	—
Blunt 16mm x 13mm ($\frac{5}{8}$ " x $\frac{1}{2}$ ")	SL46213	46213	—	45213	—

AORTIC ANUERYSM BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Large 457mm x 229mm (18" x 9")	SO46124C	46124C	45024C	45124C	—
Small 305mm x 216mm (12" x 8 ½")	SO46124D	46124D	45024D	45124D	—

ALLISON LUNG RETRACTORS



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
229mm (9")	SL46183	46183	45083	45183	—

DOUBLE ANGLED BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
13mm x 70mm (½" x 2 ¾")	SL46663	46663	—	—	—
13mm x 64mm (½" x 2 ½")	—	—	—	45663Q	—
13mm x 76mm (½" x 3")	SL46662	46662	—	45662Q	—
25mm x 64mm (1" x 2 ½")	SL46661	46661	45661A	45661Q	—
25mm x 76mm (1" x 3")	SL46660	46660	45660	45660Q	—

DOUBLE ANGLED BULB BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
25mm x 64mm (1" x 2 ½")	SL46664	46664	—	45664Q	—

HEANY BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
25mm x 102mm (1 " x 4")	SL46180	46180	45080	45180	43180
25mm x 127mm (1 " x 5")	SL46182	46182	45082	45182	43182 *

* Fixed Blade with 6" handle available, add 'S' to front of part number.

DOUBLE BLADED BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
76mm x 114mm (3" x 4 ½")	SL46190	46190	45090	45190	43190 *


* Fixed Blade with 6" handle available, add 'S' to front of part number.

HOLTGREWE BLADES




ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
60mm x 152mm (2 3/8" x 6") Malleable	SL46137	46137	45037	45137	43137

"Y" BLADDER BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
81mm x 216mm (3 1/8" x 8 1/2")	SL46139	46139	45039	45139	43139

BLADDER BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
35mm x 89mm (1 3/8" x 3 1/2")	—	46192	45092	45192	S43192

HARRIS PROSTATE RETRACTORS



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
8mm (3/8")	SL45753AM	45753AM	—	45665Q	—

MALLEABLE PROSTATE BLADE



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
280mm x 40mm (11" x 1 1/2")	—	46223	—	—	—

HILL-FERGUSON RECTAL BLADES




ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
32mm x 76mm (1 ¼" x 3")	SL46185	46185	—	45185	—
25mm x 64mm (1" x 2 ½")	SL46186	46186	—	45186	—

RECTAL BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
13mm x 117mm (½" x 4 ½")	SL46187	46187	—	45187	—

BRAU REVERSE LIP RIGID BLADES (SHORT NECK)




ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent Rigid 25mm x 100mm (1" x 4")	SL46260RET	46260RET	—	—	—
Radiolucent Rigid 25mm x 110mm (1" x 4 3/8")	SL46261RET	46261RET	—	—	—
Radiolucent Rigid 25mm x 120mm (1" x 4 3/4")	SL46262RET	46262RET	—	—	—
Radiolucent Rigid 25mm x 130mm (1" x 5 1/8")	SL46263RET	46263RET	—	—	—
Radiolucent Rigid 25mm x 140mm (1" x 5 1/2")	SL46264RET	46264RET	—	—	—
Radiolucent Rigid 25mm x 150mm (1" x 5 7/8")	SL46265RET	46265RET	—	—	—
Radiolucent Rigid 25mm x 160mm (1" x 6 1/4")	SL46266RET	46266RET	—	—	—
Radiolucent Rigid 25mm x 170mm (1" x 6 5/8")	SL46267RET	46267RET	—	—	—
Radiolucent Rigid 25mm x 180mm (1" x 7")	SL46268RET	46268RET	—	—	—
Radiolucent Rigid 25mm x 190mm (1" x 7 1/2")	SL46269RET	46269RET	—	—	—
Radiolucent Rigid 25mm x 200mm (1" x 7 7/8")	SL46270RET	46270RET	—	—	—
Radiolucent Rigid 25mm x 210mm (1" x 8 1/4")	SL46271RET	46271RET	—	—	—
Radiolucent Rigid 25mm x 220mm (1" x 8 5/8")	SL46272RET	46272RET	—	—	—
Radiolucent Rigid 25mm x 230mm (1" x 9")	SL46273RET	46273RET	—	—	—
Radiolucent Rigid 25mm x 240mm (1" x 9 1/2")	SL46274RET	46274RET	—	—	—
Radiolucent Rigid 25mm x 250mm (1" x 9 7/8")	SL46275RET	46275RET	—	—	—
Radiolucent Rigid 32mm x 100mm (1 1/4" x 4")	SL46280RET	46280RET	—	—	—
Radiolucent Rigid 32mm x 110mm (1 1/4" x 4 3/8")	SL46281RET	46281RET	—	—	—
Radiolucent Rigid 32mm x 120mm (1 1/4" x 4 3/4")	SL46282RET	46282RET	—	—	—
Radiolucent Rigid 32mm x 130mm (1 1/4" x 5 1/8")	SL46283RET	46283RET	—	—	—
Radiolucent Rigid 32mm x 140mm (1 1/4" x 5 1/2")	SL46284RET	46284RET	—	—	—
Radiolucent Rigid 32mm x 150mm (1 1/4" x 5 7/8")	SL46285RET	46285RET	—	—	—
Radiolucent Rigid 32mm x 160mm (1 1/4" x 6 1/4")	SL46286RET	46286RET	—	—	—
Radiolucent Rigid 32mm x 170mm (1 1/4" x 6 5/8")	SL46287RET	46287RET	—	—	—
Radiolucent Rigid 32mm x 180mm (1 1/4" x 7")	SL46288RET	46288RET	—	—	—

Continued on next page...

BRAU REVERSE LIP RIGID BLADES (SHORT NECK) (Continued)

ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent Rigid 32mm x 190mm (1 ¼" x 7 ½")	SL46289RET	46289RET	—	—	—
Radiolucent Rigid 32mm x 200mm (1 ¼" x 7 ⅞")	SL46290RET	46290RET	—	—	—
Radiolucent Rigid 32mm x 210mm (1 ¼" x 8 ¼")	SL46291RET	46291RET	—	—	—
Radiolucent Rigid 32mm x 220mm (1 ¼" x 8 ⅞")	SL46292RET	46292RET	—	—	—
Radiolucent Rigid 32mm x 230mm (1 ¼" x 9")	SL46293RET	46293RET	—	—	—
Radiolucent Rigid 32mm x 240mm (1 ¼" x 9 ½")	SL46294RET	46294RET	—	—	—
Radiolucent Rigid 32mm x 250mm (1 ¼" x 9 ⅞")	SL46295RET	46295RET	—	—	—

BRAU REVERSE LIP TAPERED BLADES (SHORT NECK)

ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent Tapered 25mm x 100mm (1" x 4")	SL46260TET	46260TET	—	—	—
Radiolucent Tapered 25mm x 150mm (1" x 5 ⅞")	SL46265TET	46265TET	—	—	—
Radiolucent Tapered 25mm x 200mm (1" x 7 ⅞")	SL46270TET	46270TET	—	—	—
Radiolucent Tapered 25mm x 250mm (1" x 9 ⅞")	SL46275TET	46275TET	—	—	—

BRAU REVERSE LIP MALLEABLE BLADES

ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent Malleable 25mm x 200mm (1" x 7 ⅞")	SL46192ET	46192ET	45092E	—	—
Radiolucent Malleable 32mm x 200mm (1 ¼" x 7 ⅞")	SL46193ET	46193ET	45093E	—	—

BRAU REVERSE LIP BLADES (LONG NECK)



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent Rigid 25mm x 100mm (1" x 4")	SL46194RET	46194RET	45094RE	45194RE	—
Radiolucent Rigid 25mm x 150mm (1" x 5 7/8")	SL46195RET	46195RET	45095RE	45195RE	—
Radiolucent Rigid 25mm x 200mm (1" x 7 7/8")	SL46196RET	46196RET	45096RE	45196RE	—
Radiolucent Rigid 32mm x 100mm (1 1/4" x 4")	SL46197RET	46197RET	45097RE	45197RE	—
Radiolucent Rigid 32mm x 150mm (1 1/4" x 5 7/8")	SL46198RET	46198RET	45098RE	45198RE	—
Radiolucent Rigid 32mm x 200mm (1 1/4" x 7 7/8")	SL46199RET	46199RET	45099RE	45199RE	—
Radiolucent Rigid 25mm x 255mm (1" x 10")	SL46200RET	46200RET	—	45200RE	—
Radiolucent Rigid 32mm x 255mm (1 1/4" x 10")	SL46201RET	46201RET	—	45201RE	—
Radiolucent Tapered 25mm x 150mm (1" x 5 7/8")	SL46195TET	46195TET	45095TE	—	—
Radiolucent Tapered 25mm x 200mm (1" x 7 7/8")	SL46196TET	46196TET	45096TE	—	—
Radiolucent Rigid 25mm x 100mm (1" x 4") (flat blade)	SL46194ET	46194ET	45094E	—	—
Radiolucent Rigid 25mm x 150mm (1" x 5 7/8") (flat blade)	SL46195ET	46195ET	45095E	45195E	—
Radiolucent Rigid 25mm x 200mm (1" x 7 7/8") (flat blade)	SL46196ET	46196ET	45096E	—	—
Radiolucent Rigid 32mm x 100mm (1 1/4" x 4") (flat blade)	SL46197ET	46197ET	45097E	—	—
Radiolucent Rigid 32mm x 150mm (1 1/4" x 5 7/8") (flat blade)	SL46198ET	46198ET	45098E	—	—
Radiolucent Rigid 32mm x 200mm (1 1/4" x 7 7/8") (flat blade)	SL46199ET	46199ET	45099E	—	—
Radiolucent Rigid 25mm x 255mm (1" x 10") (flat blade)	SL46200ET	46200ET	—	—	—
Radiolucent Rigid 32mm x 255mm (1 1/4" x 10") (flat blade)	SL46201ET	46201ET	—	—	—
Radiolucent Rigid Concave 25mm x 200mm (1" x 7 7/8")	SL46196CET	46196CET	—	—	—
Radiolucent Rigid Concave 32mm x 200mm (1 1/4" x 7 7/8")	SL46199CET	46199CET	—	—	—
Radiolucent Rigid Concave 25mm x 255mm (1" x 10")	SL46200CET	46200CET	—	—	—
Radiolucent Rigid Concave 32mm x 255mm (1 1/4" x 10")	SL46201CET	46201CET	—	—	—

CONCAVE RENAL VEIN BLADES (SHORT NECK)



S-Lock



Interchangeable

ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent 25mm x 80mm (1" x 3 1/8")	SL46558	46558	—	—	—
Radiolucent 25mm x 90mm (1" x 3 1/2")	SL46559	46559	—	—	—
Radiolucent 25mm x 100mm (1" x 4")	SL46560	46560	—	—	—
Radiolucent 25mm x 110mm (1" x 4 3/8")	SL46561	46561	—	—	—
Radiolucent 25mm x 120mm (1" x 4 3/4")	SL46562	46562	—	—	—
Radiolucent 25mm x 130mm (1" x 5 1/8")	SL46563	46563	—	—	—
Radiolucent 25mm x 140mm (1" x 5 1/2")	SL46564	46564	—	—	—
Radiolucent 25mm x 150mm (1" x 5 7/8")	SL46565	46565	—	—	—
Radiolucent 25mm x 160mm (1" x 6 1/4")	SL46566	46566	—	—	—
Radiolucent 25mm x 170mm (1" x 6 5/8")	SL46567	46567	—	—	—
Radiolucent 25mm x 180mm (1" x 7")	SL46568	46568	—	—	—
Radiolucent 25mm x 190mm (1" x 7 1/2")	SL46569	46569	—	—	—
Radiolucent 25mm x 200mm (1" x 7 7/8")	SL46570	46570	—	—	—
Radiolucent 25mm x 210mm (1" x 8 1/4")	SL46571	46571	—	—	—
Radiolucent 25mm x 220mm (1" x 8 5/8")	SL46572	46572	—	—	—
Radiolucent 25mm x 230mm (1" x 9")	SL46573	46573	—	—	—
Radiolucent 25mm x 240mm (1" x 9 1/2")	SL46574	46574	—	—	—
Radiolucent 25mm x 250mm (1" x 9 7/8")	SL46575	46575	—	—	—

CONCAVE LATERAL BLADES




ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
18mm x 70mm (5/8" x 2 3/4")	SL45640A	45640A	—	S45640Q	—
19mm x 55mm (3/4" x 2 1/4")	SL45150A	45150A	—	—	—
25mm x 30mm (1" x 1 1/4")	SL45156A	45156A	—	—	—
25mm x 35mm (1" x 1 1/2")	SL45157A	45157A	—	—	—
25mm x 50mm (1" x 2")	SL45151A	45151A	—	—	—
25mm x 60mm (1" x 2 3/8")	SL45158A	45158A	—	—	—
25mm x 70mm (1" x 2 3/4")	SL45159A	45159A	—	—	—
25mm x 85mm (1" x 3 3/8")	SL45152A	45152A	—	—	—
Radiolucent 23mm x 35mm (1" x 1 3/8")	SL46312AT	46312AT	—	—	—
Radiolucent 23mm x 40mm (1" x 1 1/2")	SL46314AT	46314AT	—	—	—
Radiolucent 23mm x 45mm (1" x 1 3/4")	SL46316AT	46316AT	—	—	—
Radiolucent 23mm x 50mm (1" x 2")	SL46318AT	46318AT	—	—	—
Radiolucent 23mm x 55mm (1" x 2 1/8")	SL46320AT	46320AT	—	—	—
Radiolucent 23mm x 60mm (1" x 2 3/8")	SL46322AT	46322AT	—	—	—
Radiolucent 23mm x 65mm (1" x 2 1/2")	SL46324AT	46324AT	—	—	—
Radiolucent 23mm x 70mm (1" x 2 3/4")	SL46326AT	46326AT	—	—	—
Radiolucent 23mm x 75mm (1" x 3")	SL46328AT	46328AT	—	—	—
Radiolucent 23mm x 80mm (1" x 3 1/8")	SL46330AT	46330AT	—	—	—
Radiolucent 23mm x 85mm (1" x 3 3/8")	SL46332AT	46332AT	—	—	—
Radiolucent 23mm x 90mm (1" x 3 1/2")	SL46334AT	46334AT	—	—	—
Radiolucent 23mm x 95mm (1" x 3 3/4")	SL46336AT	46336AT	—	—	—
Radiolucent 23mm x 100mm (1" x 4")	SL46338AT	46338AT	—	—	—
Radiolucent 23mm x 105mm (1" x 4 1/8")	SL46340AT	46340AT	—	—	—
Radiolucent 23mm x 110mm (1" x 4 3/8")	SL46342AT	46342AT	—	—	—
Radiolucent 23mm x 115mm (1" x 4 1/2")	SL46344AT	46344AT	—	—	—
Radiolucent 23mm x 120mm (1" x 4 3/4")	SL46346AT	46346AT	—	—	—

Continued on next page...

CONCAVE LATERAL BLADES (Continued)

ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent 22mm x 60mm (7/8" x 2 3/8")	SL45150AT	45150AT	—	—	—
Radiolucent 22mm x 70mm (7/8" x 2 3/4")	SL45151AT	45151AT	—	—	—
Radiolucent 22mm x 80mm (7/8" x 3 1/8")	SL45152AT	45152AT	—	—	—

CONCAVE LATERAL BLADES WITH TEETH



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent with Teeth 23mm x 35mm (1" x 1 3/8")	SL46312ATT	46312ATT	—	—	—
Radiolucent with Teeth 23mm x 40mm (1" x 1 1/2")	SL46314ATT	46314ATT	—	—	—
Radiolucent with Teeth 23mm x 45mm (1" x 1 3/4")	SL46316ATT	46316ATT	—	—	—
Radiolucent with Teeth 23mm x 50mm (1" x 2")	SL46318ATT	46318ATT	—	—	—
Radiolucent with Teeth 23mm x 55mm (1" x 2 1/8")	SL46320ATT	46320ATT	—	—	—
Radiolucent with Teeth 23mm x 60mm (1" x 2 3/8")	SL46322ATT	46322ATT	—	—	—
Radiolucent with Teeth 23mm x 65mm (1" x 2 1/2")	SL46324ATT	46324ATT	—	—	—
Radiolucent with Teeth 23mm x 70mm (1" x 2 3/4")	SL46326ATT	46326ATT	—	—	—
Radiolucent with Teeth 23mm x 75mm (1" x 3")	SL46328ATT	46328ATT	—	—	—
Radiolucent with Teeth 23mm x 80mm (1" x 3 1/8")	SL46330ATT	46330ATT	—	—	—
Radiolucent with Teeth 23mm x 85mm (1" x 3 3/8")	SL46332ATT	46332ATT	—	—	—
Radiolucent with Teeth 23mm x 90mm (1" x 3 1/2")	SL46334ATT	46334ATT	—	—	—
Radiolucent with Teeth 23mm x 95mm (1" x 3 3/4")	SL46336ATT	46336ATT	—	—	—
Radiolucent with Teeth 23mm x 100mm (1" x 4")	SL46338ATT	46338ATT	—	—	—
Radiolucent with Teeth 23mm x 105mm (1" x 4 1/8")	SL46340ATT	46340ATT	—	—	—
Radiolucent with Teeth 23mm x 110mm (1" x 4 3/8")	SL46342ATT	46342ATT	—	—	—
Radiolucent with Teeth 23mm x 115mm (1" x 4 1/2")	SL46344ATT	46344ATT	—	—	—
Radiolucent with Teeth 23mm x 120mm (1" x 4 3/4")	SL46346ATT	46346ATT	—	—	—
Radiolucent with Teeth 22mm x 60mm (7/8" x 2 3/8")	SL45150ATT	45150ATT	—	—	—
Radiolucent with Teeth 22mm x 70mm (7/8" x 2 3/4")	SL45151ATT	45151ATT	—	—	—
Radiolucent with Teeth 22mm x 80mm (7/8" x 3 1/8")	SL45152ATT	45152ATT	—	—	—

CONCAVE TAPERED BLADES




ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
19mm x 35mm ($\frac{3}{4}$ " x $1\frac{3}{8}$ ")	SL45643A	45643A	—	—	—
19mm x 45mm ($\frac{3}{4}$ " x $1\frac{3}{4}$ ")	SL45645A	45645A	—	S45645Q	—
19mm x 50mm ($\frac{3}{4}$ " x 2")	SL45646A	45646A	—	—	—
45mm x 50mm ($1\frac{3}{4}$ " x 2")	SL45641A	45641A	—	—	—
45mm x 65mm ($1\frac{3}{4}$ " x $2\frac{3}{4}$ ")	SL45642A	45642A	—	—	—
Radiolucent 23mm x 35mm ($1"$ x $1\frac{3}{8}$ ")	SL45643AT	45643AT	—	—	—
Radiolucent 23mm x 40mm ($1"$ x $1\frac{1}{2}$ ")	SL45644AT	45644AT	—	—	—
Radiolucent 23mm x 45mm ($1"$ x $1\frac{3}{4}$ ")	SL45645AT	45645AT	—	—	—
Radiolucent 23mm x 50mm ($1"$ x 2")	SL45646AT	45646AT	—	—	—
Radiolucent 23mm x 55mm ($1"$ x $2\frac{1}{8}$ ")	SL45647AT	45647AT	—	—	—
Radiolucent 23mm x 60mm ($1"$ x $2\frac{3}{8}$ ")	SL45648AT	45648AT	—	—	—
Radiolucent 23mm x 65mm ($1"$ x $2\frac{1}{2}$ ")	SL45649AT	45649AT	—	—	—
Radiolucent 23mm x 70mm ($1"$ x $2\frac{3}{4}$ ")	SL45650AT	45650AT	—	—	—
Radiolucent 23mm x 75mm ($1"$ x 3")	SL45651AT	—	—	—	—
Radiolucent 23mm x 80mm ($1"$ x $3\frac{1}{8}$ ")	SL45652AT	—	—	—	—
Radiolucent 51mm x 50mm ($2"$ x 2")	SL45641AT	45641AT	—	—	—
Radiolucent 51mm x 65mm ($2"$ x $2\frac{1}{2}$ ")	SL45642AT	45642AT	—	—	—
Radiolucent with Teeth 23mm x 35mm ($1"$ x $1\frac{3}{8}$ ")	SL45643ATT	45643ATT	—	—	—
Radiolucent with Teeth 23mm x 40mm ($1"$ x $1\frac{1}{2}$ ")	SL45644ATT	45644ATT	—	—	—
Radiolucent with Teeth 23mm x 45mm ($1"$ x $1\frac{3}{4}$ ")	SL45645ATT	45645ATT	—	—	—
Radiolucent with Teeth 23mm x 50mm ($1"$ x 2")	SL45646ATT	45646ATT	—	—	—
Radiolucent with Teeth 23mm x 55mm ($1"$ x $2\frac{1}{8}$ ")	SL45647ATT	45647ATT	—	—	—
Radiolucent with Teeth 23mm x 60mm ($1"$ x $2\frac{3}{8}$ ")	SL45648ATT	45648ATT	—	—	—
Radiolucent with Teeth 23mm x 65mm ($1"$ x $2\frac{1}{2}$ ")	SL45649ATT	45649ATT	—	—	—
Radiolucent with Teeth 23mm x 70mm ($1"$ x $2\frac{3}{4}$ ")	SL45650ATT	45650ATT	—	—	—

Continued on next page...

CONCAVE TAPERED BLADES (Continued)


ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent with Teeth 23mm x 75mm (1" x 3")	SL45651ATT	—	—	—	—
Radiolucent with Teeth 23mm x 80mm (1" x 3 1/8")	SL45652ATT	—	—	—	—
Radiolucent with Teeth 51mm x 50mm (2" x 2")	SL45641ATT	45641ATT	—	—	—
Radiolucent with Teeth 51mm x 65mm (2" x 2 1/2")	SL45642ATT	45642ATT	—	—	—

CONCAVE LONGITUDINAL BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
25mm x 30mm (1" x 1 1/4")	SL45160A	45160A	—	—	—
25mm x 50mm (1" x 2")	SL45161A	45161A	—	—	—
25mm x 70mm (1" x 2 3/4")	SL45162A	45162A	—	—	—
Radiolucent 23mm x 30mm (1" x 1 1/4")	SL45160AT	45160AT	—	—	—
Radiolucent 23mm x 40mm (1" x 2")	SL45161AT	45161AT	—	—	—
Radiolucent 23mm x 50mm (1" x 2 3/4")	SL45162AT	45162AT	—	—	—
Radiolucent 23mm x 60mm (1" x 2 3/8")	SL45163AT	45163AT	—	—	—
Radiolucent 23mm x 70mm (1" x 2 3/4")	SL45164AT	45164AT	—	—	—
Radiolucent 23mm x 80mm (1" x 3 1/8")	SL45165AT	45165AT	—	—	—
Radiolucent 23mm x 90mm (1" x 3 1/2")	SL45166AT	45166AT	—	—	—
Radiolucent 23mm x 100mm (1" x 4")	SL45167AT	45167AT	—	—	—
Radiolucent 23mm x 110mm (1" x 4 3/8")	SL45168AT	45168AT	—	—	—


CONCAVE MICRO BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ITEM DESCRIPTION	S-LOCK #	INTERCH. #
Radiolucent 14mm x 35mm (½" x 1 ⅜")	SL46510AT	46510AT	Radiolucent with Teeth 14mm x 35mm (½" x 1 ⅜")	SL46510ATT	46510ATT
Radiolucent 14mm x 40mm (½" x 1 ½")	SL46512AT	46512AT	Radiolucent with Teeth 14mm x 40mm (½" x 1 ½")	SL46512ATT	46512ATT
Radiolucent 14mm x 45mm (½" x 1 ¾")	SL46514AT	46514AT	Radiolucent with Teeth 14mm x 45mm (½" x 1 ¾")	SL46514ATT	46514ATT
Radiolucent 14mm x 50mm (½" x 2")	SL46516AT	46516AT	Radiolucent with Teeth 14mm x 50mm (½" x 2")	SL46516ATT	46516ATT
Radiolucent 14mm x 55mm (½" x 2 ⅙")	SL46518AT	46518AT	Radiolucent with Teeth 14mm x 55mm (½" x 2 ⅙")	SL46518ATT	46518ATT
Radiolucent 14mm x 60mm (½" x 2 ⅔")	SL46520AT	46520AT	Radiolucent with Teeth 14mm x 60mm (½" x 2 ⅔")	SL46520ATT	46520ATT
Radiolucent 14mm x 65mm (½" x 2 ½")	SL46523AT *	46523AT *	Radiolucent with Teeth 14mm x 65mm (½" x 2 ½")	SL46523ATT *	46523ATT *
Radiolucent 14mm x 70mm (½" x 2 ¾")	SL46524AT *	46524AT *	Radiolucent with Teeth 14mm x 70mm (½" x 2 ¾")	SL46524ATT *	46524ATT *
Radiolucent 14mm x 75mm (½" x 3")	SL46526AT *	46526AT *	Radiolucent with Teeth 14mm x 75mm (½" x 3")	SL46526ATT *	46526ATT *
Radiolucent 14mm x 80mm (½" x 3 ⅙")	SL46528AT *	46528AT *	Radiolucent with Teeth 14mm x 80mm (½" x 3 ⅙")	SL46528ATT *	46528ATT *

* Slightly thicker blade.

ODONTOID BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent 20mm x 90mm (¾" x 3 ½")	SL45666AT	45666AT	—	—	—

MIS TAPERED (SYD) BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent 15mm x 23mm x 35mm ($\frac{5}{8}$ " x 1" x 1 $\frac{3}{8}$ ")	SL45624AT	45624AT	—	—	—
Radiolucent 15mm x 23mm x 40mm ($\frac{5}{8}$ " x 1" x 1 $\frac{1}{2}$ ")	SL45625AT	45625AT	—	—	—
Radiolucent 15mm x 23mm x 45mm ($\frac{5}{8}$ " x 1" x 1 $\frac{3}{4}$ ")	SL45626AT	45626AT	—	—	—
Radiolucent 15mm x 23mm x 50mm ($\frac{5}{8}$ " x 1" x 2")	SL45627AT	45627AT	—	—	—
Radiolucent 15mm x 23mm x 55mm ($\frac{5}{8}$ " x 1" x 2 $\frac{1}{8}$ ")	SL45628AT	45628AT	—	—	—
Radiolucent 15mm x 23mm x 60mm ($\frac{5}{8}$ " x 1" x 2 $\frac{3}{8}$ ")	SL45629AT	45629AT	—	—	—
Radiolucent 15mm x 23mm x 65mm ($\frac{5}{8}$ " x 1" x 2 $\frac{1}{2}$ ")	SL45630AT	45630AT	—	—	—
Radiolucent 15mm x 23mm x 70mm ($\frac{5}{8}$ " x 1" x 2 $\frac{3}{4}$ ")	SL45631AT	45631AT	—	—	—
Radiolucent 15mm x 23mm x 75mm ($\frac{5}{8}$ " x 1" x 3")	SL45632AT	45632AT	—	—	—

MIDLINE BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Right 70mm x 83mm (2 $\frac{3}{4}$ " x 3 $\frac{1}{4}$ ")	SL45750A	45750A	—	—	—
Left 70mm x 83mm (2 $\frac{3}{4}$ " x 3 $\frac{1}{4}$ ")	SL45751A	45751A	—	—	—
Left 57mm x 83mm (2 $\frac{1}{4}$ " x 3 $\frac{1}{4}$ ")	SL45755A	45755A	—	—	—
Right 57mm x 83mm (2 $\frac{1}{4}$ " x 3 $\frac{1}{4}$ ")	SL45754A	45754A	—	—	—
Radiolucent Right 70mm x 83mm (2 $\frac{3}{4}$ " x 3 $\frac{1}{4}$ ")	SL45750AT	45750AT	—	—	—
Radiolucent Left 70mm x 83mm (2 $\frac{3}{4}$ " x 3 $\frac{1}{4}$ ")	SL45751AT	45751AT	—	—	—

CONCAVE MEDIAL BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
13mm x 51mm (½" x 2")	SL45212A	45212A	—	—	—
13mm x 75mm (½" x 3")	SL45214A	45214A	—	—	—
13mm x 100mm (½" x 4")	SL45216A	45216A	—	—	—
19mm x 51mm (¾" x 2")	SL45222A	45222A	—	—	—
19mm x 75mm (¾" x 3")	SL45224A	45224A	—	—	—
19mm x 100mm (¾" x 4")	SL45226A	45226A	—	—	—
25mm x 51mm (1" x 2")	SL45232A	45232A	—	—	—
25mm x 75mm (1" x 3")	SL45234A	45234A	—	—	—
25mm x 100mm (1" x 4")	SL45236A	45236A	—	—	—
32mm x 70mm (1 ¼" x 2 ¾")	SL45153A	45153A	—	—	—
38mm x 70mm (1 ½" x 2 ¾")	SL45154A	45154A	—	—	—
38mm x 85mm (1 ½" x 3 ⅜")	SL45155A	45155A	—	—	—

ARMY-NAVY BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
13mm x 55mm (½" x 2 ⅛")	SL45620A	45620A	—	S45620Q	45620Q
13mm x 65mm (½" x 2 ⅝")	SL45621A	45621A	—	S45621Q	45621Q
13mm x 75mm (½" x 3 ⅙")	SL45622A	45622A	—	S45622Q	45622Q

MEYERDING BLADES




ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
16mm x 51mm ($\frac{5}{8}$ " x 2")	SL45614A	45614A	—	—	—
16mm x 75mm ($\frac{5}{8}$ " x 3")	SL45600A	45600A	—	—	—
16mm x 100mm ($\frac{5}{8}$ " x 4")	SL45601A	45601A	—	S45601Q	—
25mm x 60mm (1" x 2 $\frac{3}{8}$ ")	SL45610A	45610A	—	—	—
25mm x 75mm (1" x 3")	SL45611A	45611A	—	S45611Q	—
25mm x 100mm (1" x 4")	SL45602A	45602A	—	—	—
25mm x 125mm (1" x 5")	SL45603A	45603A	—	—	—
51mm x 60mm (2" x 2 $\frac{3}{8}$ ")	SO45612A	45612A	—	—	—
51mm x 75mm (2" x 3")	SO45613A	45613A	—	—	—
51mm x 100mm (2" x 4")	SO45615A	45615A	—	—	—
51mm x 125mm (2" x 5")	SO45616A	45616A	—	S45616Q	—
Radiolucent 51mm x 60mm (2" x 2 $\frac{3}{8}$ ")	SO45612AET	45612AET	—	—	—
With Teeth 16mm x 51mm ($\frac{5}{8}$ " x 2")	SL45614AT	45614AT	—	—	—
With Teeth 16mm x 75mm ($\frac{5}{8}$ " x 3")	SL45600AT	45600AT	—	—	—
With Teeth 16mm x 100mm ($\frac{5}{8}$ " x 4")	SL45601AT	45601AT	—	—	—
With Teeth 25mm x 60mm (1" x 2 $\frac{3}{8}$ ")	SL45610AT	45610AT	—	—	—
With Teeth 25mm x 75mm (1" x 3")	SL45611AT	45611AT	—	—	—
With Teeth 25mm x 100mm (1" x 4")	SL45602AT	45602AT	—	—	—
With Teeth 25mm x 115mm (1" x 4 $\frac{1}{2}$ ")	SL45618AT	—	—	—	—
With Teeth 25mm x 125mm (1" x 5")	SL45603AT	45603AT	—	—	—
With Teeth 38mm x 165mm (1 $\frac{1}{2}$ " x 6 $\frac{1}{2}$ ")	SL45619AT	—	—	—	—
With Teeth 51mm x 60mm (2" x 2 $\frac{3}{8}$ ")	SO45612AT	45612AT	—	—	—
With Teeth 51mm x 75mm (2" x 3")	SO45613AT	45613AT	—	—	—
With Teeth 51mm x 100mm (2" x 4")	SO45615AT	45615AT	—	—	—
With Teeth 51mm x 125mm (2" x 5")	SO45616AT	45616AT	—	—	—
Radiolucent with Teeth 51mm x 60mm (2" x 2 $\frac{3}{8}$ ")	SO45612AETT	45612AETT	—	—	—

CONCAVE POSTERIOR LUMBAR BLADES



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
Radiolucent 23mm x 80mm ($\frac{7}{8}$ " x 3 $\frac{1}{8}$ ")	SO45670ATTC	45670ATTC	—	—	—
Radiolucent 45mm x 80mm (1 $\frac{3}{4}$ " x 3 $\frac{1}{8}$ ")	SO45671ATTC	45671ATTC	—	—	—
Radiolucent 60mm x 80mm (2 $\frac{3}{8}$ " x 3 $\frac{1}{8}$ ")	SO45672ATTC	45672ATTC	—	—	—
Radiolucent 23mm x 100mm ($\frac{7}{8}$ " x 4")	SO45673ATTC	45673ATTC	—	—	—
Radiolucent 45mm x 100mm (1 $\frac{3}{4}$ " x 4")	SO45674ATTC	45674ATTC	—	—	—
Radiolucent 60mm x 100mm (2 $\frac{3}{8}$ " x 4")	SO45675ATTC	45675ATTC	—	—	—
Radiolucent 23mm x 120mm ($\frac{7}{8}$ " x 4 $\frac{3}{4}$ ")	SO45676ATTC	45676ATTC	—	—	—
Radiolucent 45mm x 120mm (1 $\frac{3}{4}$ " x 4 $\frac{3}{4}$ ")	SO45677ATTC	45677ATTC	—	—	—
Radiolucent 60mm x 120mm (2 $\frac{3}{8}$ " x 4 $\frac{3}{4}$ ")	SO45678ATTC	45678ATTC	—	—	—
Radiolucent 23mm x 150mm ($\frac{7}{8}$ " x 5 $\frac{7}{8}$ ")	SO45679ATTC	45679ATTC	—	—	—
Radiolucent 45mm x 150mm (1 $\frac{3}{4}$ " x 5 $\frac{7}{8}$ ")	SO45680ATTC	45680ATTC	—	—	—
Radiolucent 60mm x 150mm (2 $\frac{3}{8}$ " x 5 $\frac{7}{8}$ ")	SO45681ATTC	45681ATTC	—	—	—
Radiolucent 23mm x 180mm ($\frac{7}{8}$ " x 7")	SO45682ATTC	45682ATTC	—	—	—
Radiolucent 45mm x 180mm (1 $\frac{3}{4}$ " x 7")	SO45683ATTC	45683ATTC	—	—	—
Radiolucent 60mm x 180mm (2 $\frac{3}{8}$ " x 7")	SO45684ATTC	45684ATTC	—	—	—

PRONG RETRACTORS



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
75mm (3")	SL45760A	45760A	—	—	—
100mm (4")	SL45761A	45761A	—	—	—

NERVE ROOT RETRACTORS



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
90° 5mm x 123mm (1/4" x 4 7/8")	—	45752A *	—	—	—
45° 8mm x 145mm (3/8" x 5 3/4")	—	45753A *	—	—	—
9mm x 114mm (3/8" x 4 1/2")	—	45650A	—	—	—
13mm x 114mm (1/2" x 4 1/2")	—	45651A	—	—	—

* Flexible Blade Holder (42140C) required for use.

PERESS BLADES **



ITEM DESCRIPTION	S-LOCK #	INTERCH. #	ROTATABLE #	SWIVEL #	FIXED #
35mm x 102mm (1 3/8" x 4")	SL45607A	45607A	—	—	—
35mm x 127mm (1 3/8" x 5")	SL45608A	45608A	—	—	—
35mm x 152mm (1 3/8" x 6")	SL45609A	45609A	—	—	—
60mm x 102mm (2 3/8" x 4")	SL45604A	45604A	—	—	—
60mm x 127mm (2 3/8" x 5")	SL45605A	45605A	—	—	—
60mm x 152mm (2 3/8" x 6")	SL45606A	45606A	—	—	—

** Peress Blades are available, but may require a longer lead time.

Unique System Components

Anterior Lumbar Ring 61

Laparoscopic 62-63

Cardiothoracic 64-65

 Bolling Mitral Valve 64-65

 Drake Hook..... 65

Hip 66-67

MIS PLA® 68-71


CCA 72

Body Contouring 73



EXTENSION ARM (CURRENT)

ITEM DESCRIPTION	PART #
Extension Arm with Fixed Hinge Joint	44615FH



DUAL ANGLE HANDLE (CURRENT)

ITEM DESCRIPTION	PART #
Dual Angle Handle	SL42101



RAIL CLAMP WITH EXTENSION ARM (CURRENT)

ITEM DESCRIPTION	PART #
Elite III Rail Clamp with Extension Arm Fixed Hinge Joint	44616EFH



FREE ANGLE HANDLE (CURRENT)

ITEM DESCRIPTION	PART #
Free Angle Handle	SL42102



RING FRAME HALF (CURRENT)

ITEM DESCRIPTION	PART #
Ring Frame Half 13 1/2" (2 required for use)	44224



CIRCLE JOINT (CURRENT)


ITEM DESCRIPTION	PART #
Circle Joint 1/2" x 1/4"	42100



HEX TOOL WITH T-HANDLE (CURRENT)

ITEM DESCRIPTION	PART #
Hex Tool 3/16" with T-Handle	60021

NOTICE: This page shows parts unique to the Anterior Lumbar Ring System. Additional components common to other systems are shown with standard system components and blades.



RIGID LAP HOLDERS FOR PRC (CURRENT)

REF	ITEM DESCRIPTION	PART #
A	Universal Rigid Lap Holder for Power Rail Clamp	42139U
B	Rigid Laparoscopic Holder for Power Rail Clamp	42139N
C	S-Lock Rigid Laparoscopic Holder for PRC	SL42139N

OLDER GEN. RIGID LAPAROSCOPIC HOLDERS


ITEM DESCRIPTION	PART #
Laparoscopic Instrument Holder	42139NC
Retractor/Instrument Holder	42139H
Instrument Holder	42134N
Laparoscopic Instrument Holder	42139HC



RIGID LAP HOLDERS FOR ELITE (CURRENT)


REF	ITEM DESCRIPTION	PART #
D	Universal Rigid Lap Holder for Elite (or IHRC)	42139EU
E	Rigid Laparoscopic Holder for Elite (or IHRC)	42139E
F	S-Lock Rigid Laparoscopic Holder for Elite (or IHRC)	SL42139E

Rail clamps shown on pages 6 - 7.




LAPAROSCOPIC ARTICULATING ARMS (CURRENT)

REF	ITEM DESCRIPTION	PART #
G	Laparoscopic Articulating Arm for Power Rail Clamp	42138L
H	Laparoscopic Articulating Arm for Elite (or IHRC)	42138LE



FLEXIBLE LAPAROSCOPIC HOLDERS (CURRENT)

REF	ITEM DESCRIPTION	PART #
I	Flexible Laparoscopic Holder for Power Rail Clamp	42133DN
J	Flexible Laparoscopic Holder for Elite (or IHRC)	42133DEN



OLDER GENERATION SCOPE HOLDERS

REF	ITEM DESCRIPTION	PART #
	Scope Holder "Stay Put Arm"	42137H
K	Scope Holder "Stay Put Arm" for Elite	42137HE



ADJUSTABLE SCOPE GRIP (CURRENT)

REF	ITEM DESCRIPTION	PART #
L	Adjustable Scope Grip	33181




DUO HEAD WITH TOOL GRASP KNOB (CURRENT)

REF	ITEM DESCRIPTION	PART #
P	Duo Head with Tool Grasp Knob 3/5/10mm	42160N

OLDER GENERATION DUO HEAD


ITEM DESCRIPTION	PART #
Duo Head 3/5/10mm	42160



SCOPE CLIP * (CURRENT)

REF	ITEM DESCRIPTION	PART #
M	Scope Clip 5mm	42137F
	Scope Clip 7mm	42137S
	Scope Clip 10mm	42137T
	Scope Clip 11mm	42137E

* Scope Clips may be custom ordered to fit any size scope diameter.



S-LOCK® NATHANSON HOOKS (CURRENT)

REF	ITEM DESCRIPTION	PART #
	S-Lock Small 5mm (169mm x 125mm x 12mm)	SL42200
Q	S-Lock Medium 5mm (203mm x 167mm x 16mm)	SL42201
	S-Lock Large 5mm (240mm x 222mm x 21mm)	SL42202
	S-Lock Small 6 ½mm (169mm x 125mm x 12mm)	SL42203
	S-Lock Medium 6 ½mm (203mm x 167mm x 16mm)	SL42204
	S-Lock Large 6 ½mm (240mm x 222mm x 21mm)	SL42205
R	S-Lock Nathanson Adapter	SL42161


NATHANSON HOOKS (CURRENT)

REF	ITEM DESCRIPTION	PART #
	Small 5mm (169mm x 125mm x 12mm)	42200
S	Medium 5mm (203mm x 167mm x 16mm)	42201
	Large 5mm (240mm x 222mm x 21mm)	42202
	Small 6 ½mm (169mm x 125mm x 12mm)	42203
	Medium 6 ½mm (203mm x 167mm x 16mm)	42204
	Large 6 ½mm (240mm x 222mm x 21mm)	42205

OLDER GENERATION SCOPE GRIP AND BUSHINGS

REF	ITEM DESCRIPTION	PART #
N	Scope Grip	33154
O	Scope Grip Bushing 5mm **	42136F
	Scope Grip Bushing 10mm	42136T

** Additional sizes from 5.1mm - 5.8mm also available.




L-BAR ATTACHMENTS (CURRENT)

REF	ITEM DESCRIPTION	PART #
A	CV Clamp L-Bar 6" x 8" (152mm x 203mm)	47202N
B	CV Clamp L-Bar Reversed 6" x 8" (152mm x 203mm)	47202NR

OLDER GENERATION L-BAR ATTACHMENT

REF	ITEM DESCRIPTION	PART #
C	CV Clamp L-Bar 6" x 8" (152mm x 203m)	47202




L-BAR ATTACHMENTS (CURRENT)

REF	ITEM DESCRIPTION	PART #
D	CV Clamp L-Bar 2" (51mm)	47200N
E	CV Clamp Mini L-Bar Reversed 2" (51mm)	47200NR
F	CV Clamp Mini L-Bar Offset 2" (51mm)	47200NF

OLDER GENERATION L-BAR ATTACHMENT

REF	ITEM DESCRIPTION	PART #
G	CV Clamp L-Bar 2" (51mm)	47200



T-BAR ATTACHMENTS (CURRENT)

REF	ITEM DESCRIPTION	PART #
H	CV Clamp T-Bar 8" (203mm)	47201N
I	CV Clamp T-Bar Offset 8" (203mm)	47201NF

OLDER GENERATION T-BAR ATTACHMENT

REF	ITEM DESCRIPTION	PART #
J	CV Clamp T-Bar 8" (203mm)	47201



STERNAL SPREADERS (CURRENT)

ITEM DESCRIPTION	PART #
Morse Sternal Spreader Blades 28mm (w) x 22mm (d) (1 1/8" x 7/8")	300251
Morse Sternal Spreader-Large Blades 28mm (w) x 38mm (d) (1 1/8" x 1 1/2")	300253




MALLEABLE PADDLE RETRACTORS (CURRENT)

ITEM DESCRIPTION	PART #
5" (127mm)	43100
6 ½" (165mm)	43101
8" (203mm)	43102



CUSHING RETRACTORS (CURRENT)

ITEM DESCRIPTION	PART #
5" (127mm)	43103
6 ½" (165mm)	43104
8" (203mm)	43105



BASKET RETRACTORS (CURRENT)

ITEM DESCRIPTION	PART #
Small 5" (127mm)	43106
Small 6 ½" (165mm)	43107
Small 8" (203mm)	43108
Medium 5" (127mm)	43109
Medium 6 ½" (165mm)	43110
Medium 8" (203mm)	43111
Large 5" (127mm)	43112
Large 6 ½" (165mm)	43113
Large 8" (203mm)	43114




UNIVERSAL JOINT (CURRENT)

ITEM DESCRIPTION	PART #
Universal Joint	47204



DRAKE HOOK (CURRENT)

ITEM DESCRIPTION	PART #
Drake Hook 258mm (10")	51251




HIP CROSSBAR (CURRENT)

REF	ITEM DESCRIPTION	PART #
A	Crossbar w/2 Cam Joints 7" x 19 ½"	43990


OLDER GENERATION HIP CROSSBAR

REF	ITEM DESCRIPTION	PART #
B	Crossbar w/2 Cam Joints 7" x 19 ½"	41990



CURVED ARMS (CURRENT)

REF	ITEM DESCRIPTION	PART #
C	20" Curved Arm 1	44201
D	20" Curved Arm 2	44202
E	20" Curved Arm 1 with Cleats	44203
F	20" Curved Arm 2 with Cleats	44204



HIP TETHER (CURRENT)

REF	ITEM DESCRIPTION	PART #
G	Hip Tethers (5 Pack)	42105




S-LOCK ARTICULATING ARM (CURRENT)

REF	ITEM DESCRIPTION	PART #
H	S-Lock Articulating Arm	SL42190




S-LOCK HIBBS BLADES (CURRENT)

REF	ITEM DESCRIPTION	PART #
I	Hibbs with Teeth Long 25mm x 75mm	SL46411T
	Hibbs with Teeth 25mm x 75mm (1" x 3")	SL46410T
J	Hibbs with Teeth 25mm x 115mm (1" x 4 ½")	SL46414T
	Hibbs with Teeth 25mm x 165mm (1" x 6 ½")	SL46418T



FEMUR ELEVATORS (CURRENT)

REF	ITEM DESCRIPTION	PART #
K	20" Straight Arm with Elevator Handle	44300
	25" Straight Arm with Elevator Handle	44301
L	Femur Hook 9"	44309
M	Femur Hook 12"	44312




HIP HANDLES (CURRENT)

REF	ITEM DESCRIPTION	PART #
N	Cam II Clip-On Blade Adapter 8 1/2"	42126BP
O	Orthopedic Cam II Clip-On 8 1/2" (60025 required for use)	42126OP
P	Orthopedic Retractor Handle Knob	60025
Q	Orthopedic Cam II Clip-On Hook 8 1/2"	42126OHP

OLDER GENERATION HIP HANDLES

ITEM DESCRIPTION	PART #
Cam II Clip-On Blade Adapter 8 1/2"	42126B
Orthopedic Retractor Handle Knob	42126O
Orthopedic Cam II Clip-On Hook 8 1/2"	42126OH




HOHMANN BLADES (CURRENT)

REF	ITEM DESCRIPTION	PART #
R	Hohmann Hip Blade 11"	45305
S	Curved Hohmann Narrow Hip Blade 12" [Coated]	45306C
	Curved Hohmann Narrow Hip Blade 12"	45306
T	Curved Hohmann Wide Hip Blade 12"	45307



NOTICE:

Only parts unique to our Hip Systems shown on these pages.
Additional parts shown with standard system components.




COBRA BLADES (CURRENT)

REF	ITEM DESCRIPTION	PART #
U	Cobra Hip Blade 11"	45301
V	Curved Cobra Hip Blade 12 1/2" [Coated]	45302C
	Curved Cobra Hip Blade 12 1/2"	45302




FEMORAL NECK ELEVATOR (CURRENT)

REF	ITEM DESCRIPTION	PART #
W	Femoral Neck Elevator 14"	45300



DUAL PRONG BLADES (CURRENT)

REF	ITEM DESCRIPTION	PART #
X	Dual Prong Hip Blade 11"	45309
Y	Curved Dual Prong Hip Blade 12"	45308




MIS SPINE ARTICULATING ARM (CURRENT)

REF	ITEM DESCRIPTION	PART #
A	MIS Spine Articulating Arm	42138N

OLDER GENERATION MIS SPINE ARTICULATING ARM


REF	ITEM DESCRIPTION	PART #
B	MIS Spine Articulating Arm	42138S



MIS SPINE FRAME RACKS (CURRENT)

REF	ITEM DESCRIPTION	PART #
C	MIS Spine Frame Rack 200mm (7 7/8")	60010
	MIS Spine Frame Rack 150mm (5 7/8")	60011
	MIS Spine Frame Rack 250mm (9 7/8")	60012

NOTICE:
This page shows parts unique to the MIS PLA® System only.
Rail Clamp is shown on page 6.



MIS SPINE FRAME ARMS (CURRENT)

REF	ITEM DESCRIPTION	PART #
D	MIS Spine Frame Right Arm	60008
E	MIS Spine Frame Left Arm	60009




MIS SPINE FRAME ADD-ONS (CURRENT)

REF	ITEM DESCRIPTION	PART #
F	MIS Spine Frame U-Bar Add-On	47212
G	MIS Spine Frame L-Bar Add-On	47213

OLDER GENERATION MIS SPINE FRAME ADD-ONS

REF	ITEM DESCRIPTION	PART #
H	MIS Spine Frame U-Bar Add-On	47210
I	MIS Spine Frame L-Bar Add-On	47211



MIS SPINE FRAME ATTACHMENTS (CURRENT)

REF	ITEM DESCRIPTION	PART #
J	S-Lock Quick Angle Handle 6"	SL45003G
K	¼" Side Clip x ¼" Slide with Thumb Screw	47205T
L	S-Lock Interchangeable Head with ¼" Scissor	SL45012S

OLDER GEN. MIS SPINE FRAME ATTACHMENT

REF	ITEM DESCRIPTION	PART #
M	S-Lock Angling Handle 6"	SL45002AS



MIS PLA SEQUENTIAL DILATORS (CURRENT)

REF	ITEM DESCRIPTION	PART #
N	Dilator #1 - 7mm (¼")	45831
O	Dilator #2 - 12mm (½")	45832
P	Dilator #3 - 17mm (⅝")	45833
Q	Dilator #4 - 21mm (⅞")	45834

MIS PLA® Blades shown on next pages...



CAUDAL / CEPHALAD BLADES

ITEM DESCRIPTION	LEFT PART #
Radiolucent R 22mm x 40mm ($\frac{7}{8}$ " x 1 $\frac{1}{2}$ ")	45924R
Radiolucent L 22mm x 40mm ($\frac{7}{8}$ " x 1 $\frac{1}{2}$ ")	45924L
Radiolucent R 22mm x 50mm ($\frac{7}{8}$ " x 2")	45925R
Radiolucent L 22mm x 50mm ($\frac{7}{8}$ " x 2")	45925L
Radiolucent R 22mm x 60mm ($\frac{7}{8}$ " x 2 $\frac{3}{8}$ ")	45926R
Radiolucent L 22mm x 60mm ($\frac{7}{8}$ " x 2 $\frac{3}{8}$ ")	45926L
Radiolucent R 22mm x 70mm ($\frac{7}{8}$ " x 2 $\frac{3}{4}$ ")	45927R
Radiolucent L 22mm x 70mm ($\frac{7}{8}$ " x 2 $\frac{3}{4}$ ")	45927L
Radiolucent R 22mm x 80mm ($\frac{7}{8}$ " x 3 $\frac{1}{8}$ ")	45928R
Radiolucent L 22mm x 80mm ($\frac{7}{8}$ " x 3 $\frac{1}{8}$ ")	45928L
Radiolucent R 22mm x 90mm ($\frac{7}{8}$ " x 3 $\frac{1}{2}$ ")	45929R
Radiolucent L 22mm x 90mm ($\frac{7}{8}$ " x 3 $\frac{1}{2}$ ")	45929L
Radiolucent R 22mm x 100mm ($\frac{7}{8}$ " x 4")	45930R
Radiolucent L 22mm x 100mm ($\frac{7}{8}$ " x 4")	45930L
Radiolucent R 22mm x 110mm ($\frac{7}{8}$ " x 4 $\frac{3}{8}$ ")	45931R
Radiolucent L 22mm x 110mm ($\frac{7}{8}$ " x 4 $\frac{3}{8}$ ")	45931L
Radiolucent R 22mm x 120mm ($\frac{7}{8}$ " x 4 $\frac{3}{4}$ ")	45932R
Radiolucent L 22mm x 120mm ($\frac{7}{8}$ " x 4 $\frac{3}{4}$ ")	45932L



MEDIAL / LATERAL BLADES W/LIP

ITEM DESCRIPTION	PART #
Radiolucent 13.5mm x 40mm ($\frac{1}{2}$ " x 1 $\frac{1}{2}$ ")	SL46704AT
Radiolucent 13.5mm x 50mm ($\frac{1}{2}$ " x 2")	SL46705AT
Radiolucent 13.5mm x 60mm ($\frac{1}{2}$ " x 2 $\frac{3}{8}$ ")	SL46706AT
Radiolucent 13.5mm x 70mm ($\frac{1}{2}$ " x 2 $\frac{3}{4}$ ")	SL46707AT
Radiolucent 13.5mm x 80mm ($\frac{1}{2}$ " x 3 $\frac{1}{8}$ ")	SL46708AT
Radiolucent 13.5mm x 90mm ($\frac{1}{2}$ " x 3 $\frac{1}{2}$ ")	SL46709AT
Radiolucent 13.5mm x 100mm ($\frac{1}{2}$ " x 4")	SL46710AT
Radiolucent 13.5mm x 110mm ($\frac{1}{2}$ " x 4 $\frac{3}{8}$ ")	SL46711AT
Radiolucent 13.5mm x 120mm ($\frac{1}{2}$ " x 4 $\frac{3}{4}$ ")	SL46712AT

MEDIAL / LATERAL BLADES W/FINE TEETH

ITEM DESCRIPTION	PART #
Radiolucent 13.5mm x 40mm ($\frac{1}{2}$ " x 1 $\frac{1}{2}$ ")	SL46604ATT
Radiolucent 13.5mm x 50mm ($\frac{1}{2}$ " x 2")	SL46605ATT
Radiolucent 13.5mm x 60mm ($\frac{1}{2}$ " x 2 $\frac{3}{8}$ ")	SL46606ATT
Radiolucent 13.5mm x 70mm ($\frac{1}{2}$ " x 2 $\frac{3}{4}$ ")	SL46607ATT
Radiolucent 13.5mm x 80mm ($\frac{1}{2}$ " x 3 $\frac{1}{8}$ ")	SL46608ATT
Radiolucent 13.5mm x 90mm ($\frac{1}{2}$ " x 3 $\frac{1}{2}$ ")	SL46609ATT
Radiolucent 13.5mm x 100mm ($\frac{1}{2}$ " x 4")	SL46610ATT
Radiolucent 13.5mm x 110mm ($\frac{1}{2}$ " x 4 $\frac{3}{8}$ ")	SL46611ATT
Radiolucent 13.5mm x 120mm ($\frac{1}{2}$ " x 4 $\frac{3}{4}$ ")	SL46612ATT



HOOK BLADES

ITEM DESCRIPTION	PART #
Radiolucent 40mm (1 1/2")	SL45764
Radiolucent 50mm (2")	SL45765
Radiolucent 60mm (2 3/8")	SL45766
Radiolucent 70mm (2 3/4")	SL45767
Radiolucent 80mm (3 1/8")	SL45768
Radiolucent 90mm (3 1/2")	SL45769
Radiolucent 100mm (4")	SL45770
Radiolucent 110mm (4 3/8")	SL45771
Radiolucent 120mm (4 3/4")	SL45772




PLA® MIDLINE BLADES

ITEM DESCRIPTION	PART #
Radiolucent 22mm x 30mm (7/8" x 1 1/8")	45973
Radiolucent 22mm x 40mm (7/8" x 1 1/2")	45974
Radiolucent 22mm x 50mm (7/8" x 2")	45975
Radiolucent 22mm x 60mm (7/8" x 2 3/8")	45976
Radiolucent 22mm x 70mm (7/8" x 2 3/4")	45977
Radiolucent 22mm x 80mm (7/8" x 3 1/8")	45978
Radiolucent 22mm x 90mm (7/8" x 3 1/2")	45979
Radiolucent 22mm x 100mm (7/8" x 4")	45980
Radiolucent 22mm x 110mm (7/8" x 4 3/8")	45981
Radiolucent 22mm x 120mm (7/8" x 4 3/4")	45982
Radiolucent 22mm x 130mm (7/8" x 5 1/8")	45983
Radiolucent 27mm x 30mm (1" x 1 1/8")	45984
Radiolucent 27mm x 40mm (1" x 1 1/2")	45985
Radiolucent 27mm x 50mm (1" x 2")	45986
Radiolucent 27mm x 60mm (1" x 2 3/8")	45987
Radiolucent 27mm x 70mm (1" x 2 3/4")	45988
Radiolucent 27mm x 80mm (1" x 3 1/8")	45989
Radiolucent 27mm x 90mm (1" x 3 1/2")	45990
Radiolucent 27mm x 100mm (1" x 4")	45991
Radiolucent 27mm x 110mm (1" x 4 3/8")	45992
Radiolucent 27mm x 120mm (1" x 4 3/4")	45993
Radiolucent 27mm x 130mm (1" x 5 1/8")	45994



PLA® DUAL BLADE ADAPTER

ITEM DESCRIPTION	PART #
MIS Spine Frame Dual Blade Adapter (2 required for use)	60015



CCA HANDLES (CURRENT)

REF	ITEM DESCRIPTION	PART #
A	S-Lock Hinged Handle 8"	SL45003LH
B	S-Lock Clip-On Hinged Handle 8" (to use CCA Bone Lever Blades with ½" frame components)	SL42126LH



CLIP-CLIP JOINT (CURRENT)

ITEM DESCRIPTION	PART #
------------------	--------

Clip-Clip Joint 3⁄8" x 1⁄4"	42123
-----------------------------	-------



MANIPULATOR HANDLE (CURRENT)

ITEM DESCRIPTION	PART #
------------------	--------

Manipulator Handle	42128AC
--------------------	---------



SILICONE SLEEVES (CURRENT)

ITEM DESCRIPTION	PART #
Silicone Sleeve Large (single use)	64049
Silicone Sleeve Small (single use)	64047




CCA BONE LEVER BLADES (CURRENT)

ITEM DESCRIPTION	PART #
Bone Lever 14mm x 40mm (½" x 1 ½")	SL46982AC
Bone Lever 14mm x 50mm (½" x 2")	SL46983AC
Bone Lever 14mm x 60mm (½" x 2 ¾")	SL46984AC
Bone Lever 14mm x 70mm (½" x 2 ¾")	SL46985AC
Bone Lever 14mm x 80mm (½" x 3 ⅛")	SL46986AC



NOTICE:
This page shows parts unique to the CCA System. Additional components common to other systems are shown with standard system components and blades.



TOWEL CLIP MICRO-ADJUSTABLE (CURRENT)

ITEM DESCRIPTION	PART #
Towel Clip Micro-Adjustable	35028



SALDANHA ARCHED BLADES (CURRENT)

ITEM DESCRIPTION	S-LOCK #	INTERCH. #
76mm x 140mm (3" x 5 1/2")	SO46103	46103



TUMMY TUCK BLADES (CURRENT)

ITEM DESCRIPTION	S-LOCK #	INTERCH. #
76mm x 140mm (3" x 5 1/2")	SO46400	46400
102mm x 140mm (4" x 5 1/2")	SO46401	46401



CAT'S PAW BLADES (CURRENT)

ITEM DESCRIPTION	S-LOCK #	INTERCH. #
6mm x 38mm (1/4" x 1 1/2")	SO46402	46402

 **NOTICE:**

This page shows parts unique to the Body Contouring System. Additional components common to other systems are shown with standard system components and blades.



TOWEL CLIP HOOKS (CURRENT)

ITEM DESCRIPTION	S-LOCK #	INTERCH. #
45mm x 64mm (1 3/4" x 2 1/2")	SO48001	48001

Product Index

Index by Part Number 75-87

Current Systems by Number 87

Product Measurement Guide 88-89

Service and Warranty..... 90

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
5844	Rail Extender 15" Single Clamp	7	42125N	Slide-On Joint ½" x ¼"	14, 17
5844B	Rail Extender 15" Single Clamp (reverse).....	7	42125NM	Universal Slide-On Joint Reversed ¼" x ½"	14
33154	Scope Grip	63	42125OC	Cam I Slide-on Joint ½" x ¾"	15
33181	Adjustable Scope Grip 5mm-11mm.....	63	42125S	Cam II Slide-on Joint ½" x ¼"	15
35028	Towel Clip Micro-Adjustable	73	42125WA	Cam II Slide-on Angling 8"	17
40001RX	Retractor Lite Xe with ACMI Cable	20	42125WH	Cam II Slide-on (non-angling) 8"	17
40001X	Lite Wand Xe with ACMI Cable.....	20	42126B	Cam II Clip-On Blade Adapter 8 ½"	67
40002MIS	MIS Retractor Lite LED	20	42126BP	Cam II Clip-On Blade Adapter 8 ½"	67
40002MRL	Mini Retractor Lite LED	20	42126C	Cam II Clip-On Joint ½" x ¼"	14
40002RX	Retractor Lite Xe with Storz Cable.....	20	42126H	1 Piece Universal Hinge Joint ½" x ¼"	14
40002X	Lite Wand Xe with Storz Cable.....	20	42126LC	Cam Latch Clip-On Joint ½" x ¼"	15
40004RX	Retractor Lite Xe with Wolf Cable.....	20	42126LN	Screw Down Latch Universal Clip Joint ½" x ¼"	14
40004X	Lite Wand Xe with Wolf Cable.....	20	42126NCP	Clip-On Joint ½" x ¼"	14, 17
40017	Lite Wand II Olympus Adapter.....	20	42126O	Orthopedic Cam II Clip-On 8 ½"	67
40018	Lite Wand II Wolf Adapter.....	20	42126OH	Orthopedic Cam II Clip-On Hook 8 ½"	67
40019	Lite Wand II Storz Adapter.....	20	42126OHP	Orthopedic Cam II Clip-On Hook 8 ½"	67
40019RX	Retractor Lite Xe with Olympus Cable	20	42126OP	Orthopedic Cam II Clip-On 8 ½"	67
40019X	Lite Wand Xe with Olympus Cable.....	20	42126WA	Cam II Clip-on Angling Handle 8"	17
40020RX	Retractor Lite Xe.....	20	42126WAL	Cam II Clip On Angling Handle 10"	17
40020X	Lite Wand Xe	20	42126WGLP	Cam II Clip-on Quick Angle 10"	17
40026	Lite Clip	20	42126WGP	Cam II Clip-on Quick Angle 8"	17
40026MRL	Mini Retractor Lite LED - Lite Clip.....	20	42126WH	Cam II Clip-on 8" (non-angling)	17
40732	90° Rigid Rod 10" x 20"	11	42126WHP	Cam II Clip-on 8" (non-angling)	17
40733	90° Rigid Rod 14" x 19"	11	42127	Universal (Screw Down) Joint Wrench	19
40733A	90° Rigid Rod 17" x 19"	11	42127WA	Cam II Clip-on Angling 8" (¾")	17
40733B	90° Rigid Rod 21" x 21"	11	42127WGP	Cam II Clip-on Quick Angle 8" (¾").....	17
41900BAC	Crossbar with 2 Adjustable Cam Joints 24" (11" x 13")	8	42128	Hand Held Adapter (for Fixed or Swivel retractors)	19
41900BC	Crossbar with 2 Cam Joints 24" (11" x 13").....	8	42128AC	Manipulator Handle (CCA)	72
41900BCW	Crossbar with 2 Cam Joints 24" (8" x 16").....	8	42128C	Locking Cam Hand Held Adapter	19
41900F	Bilateral Crossbar 7 ¼" x 13 ½" x 7 ¼"	8	42133DEN	Laparoscopic Flexible Holder for Elite	62
41900FC	Bilateral Crossbar with 2 Joints	8	42133DN	Laparoscopic Flexible Holder for Power Rail Clamp	62
41900H	Bilateral Crossbar Hinged 34" (11" x 8 ½" x 11").....	8	42134N	Instrument Holder	62
41900HL	Bilateral Crossbar Hinged 41 ½" (14" x 10" x 14").....	8	42136F	Scope Grip Bushing 5mm	63
41900PH	Bilateral Crossbar Hinged 34" (11 ½" x 8 ½" x 11 ½").....	8	42136T	Scope Grip Bushing 10mm	63
41900PHL	Bilateral Crossbar Hinged 41 ½" (14 ½" x 10" x 14 ½").....	8	42137E	Scope Clip 11mm.....	63
41902AC	Elite II Rail Clamp with 1 Cam Joint 18"	6	42137F	Scope Clip 5mm.....	63
41902ACL	Elite II Rail Clamp with 1 Cam Joint 16"	6	42137H	Scope Holder "Stay Put Arm"	62
41902ACS	Elite II Rail Clamp with 1 Cam Joint 11"	6	42137HE	Scope Holder "Stay Put Arm" for Elite	62
41902IRC	Infinite Height Rail Clamp with 1 Cam Joint	6	42137S	Scope Clip 7mm.....	63
41902S	Elite II Rail Clamp with 1 (¾") Cam Joint 11"	6	42137T	Scope Clip 10mm.....	63
41903AC	Elite II Rail Clamp with 2 Cam Joints 18"	6	42138L	Laparoscopic Articulating Arm for Power Rail Clamp	62
41903ACL	Elite II Rail Clamp with 2 Cam Joints 16"	6	42138LE	Laparoscopic Articulating Arm for Elite.....	62
41903IRC	Infinite Height Rail Clamp with 2 Cam Joints.....	6	42138N	MIS Spine Frame Articulating Arm	68
41904AC	Elite II Rail Clamp with 1 Cam Joint 22"	6	42138S	MIS Spine Frame Articulating Arm	68
41905AC	Elite II Rail Clamp with 2 Cam Joints 22"	6	42139E	Laparoscopic Rigid Holder for Elite.....	62
41906AC	Elite II Rail Clamp with 1 Cam Joint 26"	6	42139EU	Universal Laparoscopic Rigid Holder for Elite	62
41907AC	Elite II Rail Clamp with 2 Cam joints 26"	6	42139H	Laparoscopic Retractor Holder on Rod	62
41917	Rail Extender 14" Long with 2 ¼" Offset Dual Clamp	7	42139HC	Laparoscopic Instrument Holder	62
41927	Jackson Frame Adapter.....	7	42139N	Laparoscopic Rigid Holder for Power Rail Clamp	62
41929	Rail Extender 22" Dual Clamp.....	7	42139NC	Lap Retractor Holder on Rod	62
41933	Hysterectomy Tray	21	42139U	Universal Laparoscopic Rigid Holder for Power Rail Clamp	62
41938	Rail Extender 20" Single Clamp	7	42140C	Flexible Blade Holder Cam	19
41950	Thompson Power Rail Clamp	7	42160	Duo Head 3mm/5mm/10mm	63
41990	Crossbar with 2 Cam Joints 7" x 19 ½"	8, 66	42160N	Duo Head 3mm/5mm/10mm with Tool Grasp Knob	63
42100	Circle Joint ½" x ¼"	61	42200	Nathanson Hook Sm (5mm) 169mm x 125mm x 12mm (6 ⅝" x 5" x ½")	63
42105	Hip Tethers 24" - 5 Pack.....	66	42201	Nathanson Hook Med (5mm) 203mm x 167mm x 16mm (8" x 6 ½" x ⅝")	63
42110C	Cam I Slide-on Joint ½" x ½"	12	42202	Nathanson Hook Lg (5mm) 240mm x 222mm x 21mm (9 ½" x 8 ¾" x ⅞").....	63
42110N	Slide-on Joint ½" x ½"	12	42203	Nathanson Hook Sm (6.5mm) 169mm x 125mm x 12mm (6 ⅝" x 5" x ½").....	63
42110S	Cam II Slide-on Joint ½" x ½"	12	42204	Nathanson Hook Med (6.5mm) 203mm x 167mm x 16mm (8" x 6 ½" x ⅝").....	63
42110WH	Cam II Slide-on (non-angling) with ½" Handle Rod	17	42205	Nathanson Hook Lg (6.5mm) 240mm x 222mm x 21mm (9 ½" x 8 ¾" x ⅞").....	63
42111H	1 Piece Universal Hinge Joint ½" x ½"	13	42228	Manipulator Handle.....	19
42111LC	Latch Cam Joint ½" x ½"	12	42520	Fixed ½" Handle Weinberg 102mm x 152mm (4" x 6")	28
42111LN	Latch Joint ½" x ½"	13	42521	Fixed ½" Handle Weinberg 83mm x 133mm (3 ¼" x 5 ¼")	28
42112C	Cam I Serrated Slide-on Joint ½" x ½"	12	42522	Fixed ½" Handle Weinberg 102mm x 203mm (4" x 8")	28
42112N	Universal Serrated Slide-on Joint ½" x ½"	13	43032	Angling Head Wrench	19
42123	Clip-Clip Joint ⅝" x ¼"	72	43100	Malleable Paddle 5" (127mm).....	65
42124C	Cam I Slide-on Joint ¼" x ¼"	15	43101	Malleable Paddle 6 ½" (165mm).....	65
42125C	Cam I Slide-On Joint ½" x ¼"	15	43102	Malleable Paddle 8" (203mm).....	65

PRODUCT INDEX

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
43103	Cushing 5" (127mm).....	65	43900BC	Crossbar with 2 Cam III Joints 24" (11" x 13").....	8
43104	Cushing 6 ½" (165mm).....	65	43900BCW	Crossbar with 2 Cam III Joints 24" (8" x 16").....	8
43105	Cushing 8" (203mm).....	65	43902AC	Elite III Rail Clamp with 1 Joint 18".....	6
43106	Small Basket 5" (127mm).....	65	43902ACL	Elite III Rail Clamp with 1 Joint 16".....	6
43107	Small Basket 6 ½" (165mm).....	65	43902ACS	Elite III Rail Clamp with 1 Joint 11".....	6
43108	Small Basket 8" (203mm).....	65	43902IRC	Infinite Height Rail Clamp with 1 Cam III Joint.....	6
43109	Medium Basket 5" (127mm).....	65	43902S	Elite III Rail Clamp with 1 ¾" Joint 11".....	6
43110	Medium Basket 6 ½" (165mm).....	65	43903AC	Elite III Rail Clamp with 2 Joints 18".....	6
43111	Medium Basket 8" (203mm).....	65	43903ACL	Elite III Rail Clamp with 2 Joints 16".....	6
43112	Large Basket 5" (127mm).....	65	43903IRC	Infinite Height Rail Clamp with 2 Cam III Joints.....	6
43112C	Cam III Serrated Slide-on Joint.....	12	43904AC	Elite III Rail Clamp with 1 Joint 22".....	6
43113	Large Basket 6 ½" (165mm).....	65	43905AC	Elite III Rail Clamp with 2 Joints 22".....	6
43114	Large Basket 8" (203mm).....	65	43906AC	Elite III Rail Clamp with 1 Joint 26".....	6
43118	Fixed Malleable 25mm x 203mm (1" x 8").....	30	43907AC	Elite III Rail Clamp with 2 Joints 26".....	6
43118R	Fixed Malleable Renal Vein 25mm x 203mm (1" x 8").....	38	43990	Crossbar with 2 Cam III Joints 7" x 19 ½".....	8, 66
43119C	Fixed Rigid Renal Vein 25mm x 191mm (1" x 7 ½").....	37	44012N	Straight Arm 12".....	10
43122	Fixed Malleable 51mm x 203mm (2" x 8").....	30	44014	Angled Arm 12".....	9
43122A	Fixed Malleable 51mm x 254mm (2" x 10").....	30	44014C	Extension Arm with Clip Joint 14".....	10
43124	Fixed Malleable 64mm x 254 (2 ½" x 10").....	30	44014H	Extension Arm 12" with Hinge Joint (¾").....	10
43124A	Fixed Malleable 76mm x 254mm (3" x 10").....	30	44014HS	Extension Arm 12" with Hinge Joint.....	10
43124B	Fixed Malleable 102 x 254mm (4" x 10").....	30	44014LC	Extension Arm 12" with Cam Latch Joint.....	10
43125	Fixed Richardson 19mm x 25mm (¾" x 1").....	35	44014N	Angled Arm 12".....	9
43125A	Fixed Richardson 19mm x 51mm (¾" x 2").....	35	44014S	Angled Arm 12" (6" x 6").....	9
43126	Fixed Richardson 25mm x 32mm (1" x 1 ¼").....	35	44016C	Extension Arm with Clip Joint 16".....	10
43127	Fixed Kelly 38mm x 51mm (1 ½" x 2").....	28	44020C	Extension Arm with Clip Joint 20".....	10
43128	Fixed Kelly 51mm x 64mm (2" x 2 ½").....	28	44109B	Angled Closure Arm 9" (4 ½" x 4 ½").....	10
43129	Fixed Kelly 64mm x 76mm (2 ½" x 3").....	28	44110SP	Lower Abdominal Bar 16" (10" x 6").....	9
43130	Fixed Kelly 76mm x 89mm (3" x 3 ½").....	28	44111SP	Lower Abdominal Bar 18" (10" x 8").....	9
43134	Fixed Fence 102mm x 102mm (4" x 4").....	29	44112SP	Lower Abdominal Bar 20" (14" x 6").....	9
43135	Fixed Fence 102mm x 127mm (4" x 5").....	29	44113B	Angled Closure Arm 13" (6 ½" x 6 ½").....	10
43136	Fixed Fence 102mm x 152mm (4" x 6").....	29	44115	Hinged Extension Arm 12" x 12".....	10
43136A	Fixed Fence 102mm x 203mm (4" x 8").....	29	44117	Angled Arm 17" (10" x 7") @ 45° (¾").....	9
43137	Fixed Holtgrewe 60mm x 152mm (2 ¾" x 6") [Malleable].....	43	44117N	Angled Arm 17" (10" x 7") @ 45°.....	9
43139	Fixed "Y" Bladder 81mm x 216mm (3 ½" x 8 ½").....	43	44118B	Angled Closure Arm 18" (9" x 9").....	10
43140	Fixed Balfour [with lips] 70mm x 73mm (2 ¾" x 3").....	25	44118C	Straight Arm with 1 Cam Joint 18".....	10
43140C	Fixed Balfour [small lips] 70mm x 73mm (2 ¾" x 3").....	25	44118N	Angled Arm 18" (8" x 10") @ 45°.....	9
43140P	Fixed Balfour [with lips] 76mm x 51mm (3" x 2").....	25	44119	Angled Arm 19" (12" x 7") @ 45° (¾").....	9
43141	Fixed Balfour [with lips] 83mm x 73mm (3 ¼" x 3").....	25	44120	Angled Arm 20" (10" x 10") @ 45°.....	9
43141C	Fixed Balfour [small lips] 83mm x 73mm (3 ¼" x 3").....	25	44120C	Straight Arm with 1 Cam Joint 20".....	10
43141D	Fixed Balfour [with lips] 83mm x 123mm (3 ¼" x 5").....	25	44120N	Straight Arm 20".....	10
43141F	Fixed Balfour [with lips] 83mm x 175mm (3 ½" x 6 ¾").....	25	44121B	Angled Closure Arm 21" (10 ½" x 10 ½").....	10
43142	Fixed Balfour [with lips] 32mm x 32mm (1 ¼" x 1 ¼").....	25	44122DB	Double Angled Arm 22" (7" x 9" x 7").....	9
43142P	Fixed Balfour [with lips] 51mm x 54mm (2" x 2").....	25	44124	Angled Arm 24" (6" x 18") @ 45°.....	9
43143	Fixed Fenestrated 44mm x 89mm (1 ¾" x 3 ½").....	38	44124N	Angled Arm 24" (8" x 16") @ 45°.....	9
43150	Fixed Richardson 51mm x 102mm (2" x 4").....	34	44124W	Angled Arm 24" (11" x 13") @ 45°.....	9
43151	Fixed Richardson 51mm x 127mm (2" x 5").....	34	44126	Angled Arm 26" (6" x 20") @ 45°.....	9
43152	Fixed Richardson 51mm x 152mm (2" x 6").....	34	44128DB	Double Angled Arm 28" (7 ¼" x 13 ½" x 7 ¼").....	9
43160	Fixed Harrington 64mm x 152mm (2 ½" x 6").....	27	44200	Extension Arm 7 ½" x 8 ½" x 7 ½".....	10
43161	Fixed Harrington 114mm x 254mm (4 ½" x 10").....	27	44201	20" Curved Arm 1.....	66
43162P	Fixed Harrington 38mm x 127mm (1 ½" x 5").....	27	44202	20" Curved Arm 2.....	66
43163	Fixed Harrington 64mm x 203mm (2 ½" x 8").....	27	44203	20" Curved Arm with Cleats 1.....	66
43164	Fixed Harrington 64mm x 254mm (2 ½" x 10").....	27	44204	20" Curved Arm with Cleats 2.....	66
43165	Fixed Harrington 64mm x 102mm (2 ½" x 4").....	27	44218	Ring Half Hinged 14 ½".....	9
43170	Fixed Deaver 25mm x 127mm (1" x 5").....	36	44220	Ring Extension Hinged 3".....	9
43172	Fixed Deaver 38mm x 127mm (1 ½" x 5").....	36	44224	Ring Frame Half 13 ½".....	61
43174	Fixed Deaver 51mm x 127mm (2" x 5").....	36	44300	Straight Arm with Elevator Handle 20".....	66
43175	Fixed Deaver 64mm x 127mm (2 ½" x 5").....	36	44301	Straight Arm with Elevator Handle 25".....	66
43176	Fixed Deaver 76mm x 127mm (3" x 5").....	36	44309	Femur Hook 9".....	66
43180	Fixed Heaney 25mm x 102mm (1" x 4").....	42	44312	Femur Hook 12".....	66
43182	Fixed Heaney 25mm x 127mm (1" x 5").....	42	44614LC	Extension Arm 14" with Cam Latch Mount.....	9, 10
43188	Fixed St Marks [with lip] 64mm x 178mm (2 ½" x 7").....	35	44615FH	Extension Arm with Fixed Hinge Joint.....	61
43190	Fixed Double Bladed 76mm x 114mm (3" x 4 ½").....	42	44616EFH	Elite III Rail Clamp with Extension Arm Fixed Hinge Joint.....	61
43191	Fixed Malleable Finger 152mm (6").....	31	44643	Extension Arm 90° 10" x 6".....	11
43214	Fixed Murphy Sharp 25mm x 19mm (1" x ¾").....	39	44646C	Upright Arm with Serrated Cam Joint 20".....	11
43241A	Fixed Wide Balfour [no lips] 83mm x 57mm (3 ¼" x 2 ¼").....	26	44646N	Upright Arm with Universal Joint 20".....	11
43241B	Fixed Wide Balfour [no lips] 100mm x 61mm (4" x 2 ¾").....	26	44733	90° Collapsible Rod 14" x 19".....	11
43241C	Fixed Wide Balfour [no lips] 114mm x 63mm (4 ½" x 2 ½").....	26	44733A	Collapsible 90° Rod 17" x 19".....	11
43800	Crossbar with 1 Cam III Joint Mini-Incision.....	8	44733AC	Collapsible 90° Rod with 1 Cam Joint 17" x 19".....	11

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
44733B	Collapsible 90° Rod 21" x 21"	11	45018A	Rotatable Malleable 25mm x 254mm (1" x 10")	30
44733BC	Collapsible 90° Rod with 1 Cam Joint 21" x 21"	11	45018B	Rotatable Malleable 38mm x 254mm (1 ½" x 10")	30
44733C	Collapsible 90° Rod with 1 Cam Joint 14" x 19"	11	45018CE	Rotatable Radiolucent Rigid Renal Vein 25mm x 152mm (1" x 6")	37
44734	90° Collapsible Rod with Grooves	11	45018R	Rotatable Renal Vein 25mm x 203mm (1" x 8")	38
44800	Crossbar with 1 Cam Joint Mini-Incision	8	45018RET	Rotatable Radiolucent Renal Vein 25mm x 203mm (1" x 8")	37
45000	Rotatable Handle 6"	17	45019	Rotatable Malleable 25mm x 127mm (1" x 5")	30
45001	Rotatable Handle 8"	17	45019A	Rotatable Malleable 25mm x 152mm (1" x 6")	31
45001L	Rotatable Handle 10"	17	45019B	Rotatable Rigid Concave Renal Vein 25mm x 127mm (1" x 5 ½")	37
45003	Interchangeable Handle 7 ¾" (non-angling)	17	45019BE	Rotatable Radiolucent Malleable Renal Vein 25mm x 140mm (1" x 5 ½")	37
45003A	Angling Handle 8"	17	45019BET	Rotatable Radiolucent Concave Renal 25mm x 140mm (1" x 5 ½")	37
45003AS	Angling II Handle 8"	17	45019C	Rotatable Rigid Concave Renal Vein 25mm x 191mm (1" x 7 ½")	37
45003ASL	Angling II Handle 10"	17	45019CE	Rotatable Radio. Malleable Renal Vein 25mm x 191mm (1" x 7 ½") [Aluminum]	37
45003GS	Quick Angle Handle 8"	17	45019RE	Rotatable Radiolucent Rigid Renal Vein 25mm x 140mm (1" x 5 ½")	37
45003GSL	Quick Angle Handle 10"	17	45020	Rotatable Malleable 38mm x 102mm (1 ½" x 4")	31
45003L	Interchangeable Handle 8" (non-angling)	17	45020A	Rotatable Malleable 38mm x 127mm (1 ½" x 5")	31
45003LM	Interchangeable Handle 10" (non-angling)	17	45020B	Rotatable Malleable 38mm x 152mm (1 ½" x 6")	31
45004	Rotatable Thompson Adapter Handle for Omni-Tract® 8 ¾"	18	45020RE	Rotatable Radiolucent Rigid Renal Vein 25mm x 178mm (1" x 7")	37
45004AL	Angling Thompson Adapter Handle for Omni-Tract® 8 ¾"	18	45021RE	Rotatable Radiolucent Rigid Renal Vein 25mm x 254mm (1" x 10")	37
45004AL	Quick Angle Thompson Adapter Handle for Omni-Tract® 8 ¾"	18	45022	Rotatable Malleable 51mm x 203mm (2" x 8")	30
45004L	Thompson Adapter Handle for Omni-Tract® 8 ¾" (non-angling)	18	45022A	Rotatable Malleable 51mm x 254mm (2" x 10")	30
45005	Ratchet Retractor Handle Rotatable	18	45022E	Rotatable Radiolucent Malleable 51mm x 203mm (2" x 8")	31
45005AL	Ratchet Angling Retractor Handle 10 ¼"	18	45022RE	Rotatable Radiolucent Rigid Renal Vein 25mm x 305mm (1" x 12")	37
45005GL	Ratchet Quick Angle Retractor Handle 10 ¼"	18	45023	Rotatable Malleable 38mm x 203mm (1 ½" x 8")	30
45005L	Ratchet Retractor Handle 10 ¼" (non-angling)	18	45024	Rotatable Malleable 64mm x 254mm (2 ½" x 10")	30
45006BS	Micro-Adjustable II Slide-on 10 ¾" Bariatric (non-angling)	16	45024A	Rotatable Malleable 76mm x 254mm (3" x 10")	30
45006BT	Micro-Adjustable II Slide-on 8 ¾" Bariatric (non-angling)	16	45024B	Rotatable Malleable 102mm x 254mm (4" x 10")	30
45006BTC	Micro-Adjustable II Clip-on 10 ¾" Bariatric (non-angling)	16	45024C	Rotatable Aortic Aneurysm Large 457mm x 229mm (18" x 9")	40
45006BTCP	Micro-Adjustable II Clip-on 10 ¾" Bariatric (non-angling)	16	45024D	Rotatable Aortic Aneurysm Small 305mm x 216mm (12" x 8 ½")	40
45006CA	Micro-Adjustable II Clip-on Angling 10"	16	45027	Rotatable Kelly 38mm x 51mm (1 ½" x 2")	28
45006CAR	Micro-Adjustable II Clip-on Angling 15"	16	45028	Rotatable Kelly 51mm x 64mm (2" x 2 ½")	28
45006CGP	Micro-Adjustable II Clip-on Quick Angle 10"	16	45029	Rotatable Kelly 64mm x 76mm (2 ½" x 3")	28
45006CGRP	Micro-Adjustable II Clip-on Quick Angle 15"	16	45030	Rotatable Kelly 76mm x 89mm (3" x 3 ½")	28
45006CL	Micro-Adjustable Latch 10 ¾" (non-angling)	16	45034	Rotatable Fence 102mm x 102mm (4" x 4")	29
45006CR	Micro-Adjustable II Clip-on 15" Handle (non-angling)	16	45035	Rotatable Fence 102mm x 127mm (4" x 5")	29
45006CS	Micro-Adjustable II Slide-on 10 ¾" (non-angling)	16	45036	Rotatable Fence 102mm x 152mm (4" x 6")	29
45006LH	Micro-Adjustable Slide-on 10 ¾" (non-angling)	16	45036A	Rotatable Fence 102mm x 203mm (4" x 8")	29
45006LHB	Micro-Adjustable Slide-on 10 ¾" (non-angling)	16	45036E	Rotatable Radiolucent Fence 102mm x 152mm (4" x 6")	29
45006ST	Micro-Adjustable II Slide-on 8 ¾" (non-angling)	16	45037	Rotatable Holtgrewe Blade 60mm x 152mm (2 ¾" x 6") [Malleable]	43
45006STA	Micro-Adjustable II Slide-on Angling 10"	16	45039	Rotatable "Y" Bladder 81mm x 216mm (3 ½" x 8 ½")	43
45006STC	Micro-Adjustable II Clip-on 10 ¾" (non-angling)	16	45040	Rotatable Balfour [with lips] 70mm x 73mm (2 ¾" x 3")	25
45006STCP	Micro-Adjustable II Clip-on 10 ¾" (non-angling)	16	45040C	Rotatable Balfour [small lips] 70mm x 73mm (2 ¾" x 3")	25
45007AH	Micro-Adjustable Universal Slide-on Angling 10"	16	45040CE	Rotatable Radiolucent Balfour [small lips] 65mm x 72mm (2 ¾" x 3")	25
45007BTC	Micro-Adjustable II Clip-on 8 ¾" Bariatric (non-angling)	16	45040E	Rotatable Radiolucent Balfour [with lips] 65mm x 72mm (2 ¾" x 3")	25
45007BTCP	Micro-Adjustable II Clip-on 8 ¾" Bariatric (non-angling)	16	45040P	Rotatable Balfour [with lips] 76mm x 51mm (2 ¾" x 2")	25
45007CA	Micro-Adjustable II Clip-on Angling 7 ¾"	16	45041	Rotatable Balfour [with lips] 83mm x 73mm (3 ¼" x 3")	25
45007CGP	Micro-Adjustable II Clip-on Quick Angle 7 ¾"	16	45041C	Rotatable Balfour [small lips] 83mm x 73mm (3 ¼" x 3")	25
45007H	Micro-Adjustable Universal Slide-on 8 ¾" (non-angling)	16	45041D	Rotatable Balfour [with lips] 83mm x 123mm (3 ¼" x 5")	25
45007P	Micro-Adjustable II Clip-on Angling 20 ½"	16	45041E	Rotatable Radiolucent Balfour [with lips] 83mm x 70mm (3 ¼" x 3")	25
45007PGP	Micro-Adjustable II Clip-on Quick Angle 20"	16	45041F	Rotatable Balfour [with lips] 83mm x 175mm (3 ¼" x 6 ⅞")	25
45007STA	Micro-Adjustable II Slide-on Angling 7 ¾"	16	45042	Rotatable Balfour [with lips] 32mm x 32mm (1 ¼" x 1 ¼")	25
45007STC	Micro-Adjustable II Clip-on 8 ¾" (non-angling)	16	45042P	Rotatable Balfour [with lips] 51mm x 54mm (2" x 2")	25
45007STCP	Micro-Adjustable II Clip-on 8 ¾" (non-angling)	16	45045E	Rotatable Radiolucent Splanchnic 51mm x 114mm (2" x 4 ½")	29
45007STG	Micro-Adjustable II Slide-on Quick Angle 7 ¾"	16	45046E	Rotatable Radiolucent Splanchnic 51mm x 152mm (2" x 6")	29
45008C	Micro-Adjustable II Rotatable Latch 10"	16	45047E	Rotatable Radiolucent Splanchnic 51mm x 203mm (2" x 8")	29
45008CH	Micro-Adjustable II Slide-on Rotatable 10"	16	45048E	Rotatable Radiolucent Splanchnic 25mm x 152mm (1" x 6")	29
45008CS	Micro-Adjustable II Slide-on Rotatable 10"	16	45049E	Rotatable Radiolucent Splanchnic 25mm x 203mm (1" x 8")	29
45008STC	Micro-Adjustable II Clip-on Rotatable 10"	16	45050	Rotatable Richardson 51mm x 102mm (2" x 4")	34
45008STCP	Micro-Adjustable II Clip-on Rotatable 10"	16	45051	Rotatable Richardson 51mm x 127mm (2" x 5")	34
45009CL	Micro-Adjustable Latch Angling 10"	16	45052	Rotatable Richardson 51mm x 152mm (2" x 6")	34
45009CS	Micro-Adjustable Slide-on Angling 10"	16	45053	Rotatable Richardson 51mm x 178mm (2" x 7")	34
45010A	Mini Micro-Adjustable Clip-On Angling 8" (cervical use only)	16	45055	Rotatable Mini-Open Richardson 51mm x 76mm (2" x 3")	33
45014AL	Angling Thompson Adapter Handle for SynFrame® 8 ¾"	18	45060	Rotatable Harrington 64mm x 152mm (2 ½" x 6")	27
45014GL	Quick Angle Thompson Adapter Handle for SynFrame® 8 ¾"	18	45060E	Rotatable Radiolucent Harrington 64mm x 152mm (2 ½" x 6")	27
45014L	Thompson Adapter Handle for SynFrame® 8 ¾"	18	45061	Rotatable Harrington 114mm x 254mm (4 ½" x 10")	27
45015	Rotatable Malleable 13mm x 64mm (½" x 2 ½")	30	45063	Rotatable Harrington 64mm x 203mm (2 ½" x 8")	27
45016	Rotatable Malleable 13mm x 152mm (½" x 6")	30	45064	Rotatable Harrington 64mm x 254mm (2 ½" x 10")	27
45017	Rotatable Malleable 13mm x 203mm (½" x 8")	30	45070	Rotatable Deaver 25mm x 127mm (1" x 5")	36
45018	Rotatable Malleable 25mm x 203mm (1" x 8")	30	45072	Rotatable Deaver 38mm x 127mm (1 ½" x 5")	36

PRODUCT INDEX

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
45074	Rotatable Deaver 51mm x 127mm (2" x 5")	36	45136A	Swivel Fence 102mm x 203mm (4" x 8")	29
45075	Rotatable Deaver 64mm x 127mm (2 ½" x 5")	36	45137	Swivel Holtgrewe 60mm x 152mm (2 ¾" x 6")	43
45076	Rotatable Deaver 76mm x 127mm (3" x 5")	36	45139	Swivel "Y" Bladder 81mm x 216mm (3 ⅜" x 8 ½")	43
45076E	Rotatable Radiolucent Deaver 76mm x 127mm (3" x 5")	36	45140	Swivel Balfour [with lips] 70mm x 73mm (2 ¾" x 3")	25
45080	Rotatable Heaney 25mm x 102mm (1" x 4")	42	45140C	Swivel Balfour [small lips] 70mm x 73mm (2 ¾" x 3")	25
45082	Rotatable Heaney 25mm x 127mm (1" x 5")	42	45140P	Swivel Balfour [with lips] 76mm x 51mm (2 ¾" x 2")	25
45083	Rotatable Allison Lung 229mm (9")	40	45141	Swivel Balfour [with lips] 83mm x 73mm (3 ¼" x 3")	25
45087A	Rotatable St Marks [no lip] 64mm x 178mm (2 ½" x 7")	35	45141C	Swivel Balfour [small lips] 83mm x 73mm (3 ¼" x 3")	25
45088	Rotatable St Marks [with lip] 64mm x 178mm (2 ½" x 7")	35	45141D	Swivel Balfour [with lips] 83mm x 123mm (3 ¼" x 5")	25
45089	Rotatable St Marks [with lip] 64mm x 140mm (2 ½" x 5 ½")	35	45141F	Swivel Balfour [with lips] 83mm x 175mm (3 ¼" x 6 ⅞")	25
45090	Rotatable Double Bladed 76mm x 114mm (3" x 4 ½")	42	45142	Swivel Balfour [with lips] 32mm x 32mm (1 ¼" x 1 ¼")	25
45091	Rotatable Malleable Finger 152mm (6")	31	45142P	Swivel Balfour [with lips] 51mm x 54mm (2" x 2")	25
45092	Rotatable Bladder 35mm x 89mm (1 ⅜" x 3 ½")	43	45143	Swivel Fenestrated 44mm x 89mm (1 ¾" x 3 ½")	38
45092E	Rotatable Reverse Lip Malleable 25mm x 203mm (1" x 8")	47	45144	Swivel Fenestrated 48mm x 127mm (1 ⅞" x 5")	38
45093E	Rotatable Reverse Lip Malleable 32mm x 203mm (1 ¼" x 8")	47	45145	Swivel Splanchnic 51mm x 114mm (2" x 4 ½")	29
45094E	Rotatable Reverse Lip Rigid 25mm x 102mm (1" x 4")	48	45145E	Swivel Radiolucent Splanchnic 51mm x 152mm (2" x 6")	29
45094RE	Rotatable Reverse Lip Rigid 25mm x 102mm (1" x 4")	48	45150	Swivel Richardson 51mm x 102mm (2" x 4")	34
45095E	Rotatable Reverse Lip Rigid 25mm x 152mm (1" x 6")	48	45146	Swivel Splanchnic 51mm x 152mm (2" x 6")	29
45095RE	Rotatable Reverse Lip Rigid 25mm x 152mm (1" x 6")	48	45150A	Concave Lateral 19mm x 55mm (⅞" x 2 ⅞")	50
45095TE	Rotatable Reverse Lip Tapered 25mm x 152mm (1" x 6")	48	45150AT	Radiolucent Concave Lateral 22mm x 60mm (⅞" x 2 ¾")	51
45096E	Rotatable Reverse Lip Rigid 25mm x 203mm (1" x 8")	48	45150ATT	Radiolucent Concave Lateral with Teeth 22mm x 60mm (⅞" x 2 ¾")	51
45096RE	Rotatable Reverse Lip Rigid 25mm x 203mm (1" x 8")	48	45151	Swivel Richardson 51mm x 127mm (2" x 5")	34
45096TE	Rotatable Reverse Lip Tapered 25mm x 203mm (1" x 8")	48	45151A	Concave Lateral 25mm x 50mm (1" x 2")	50
45097E	Rotatable Reverse Lip Rigid 32mm x 102mm (1 ¼" x 4")	48	45151AT	Radiolucent Concave Lateral 22mm x 70mm (⅞" x 2 ¾")	51
45097RE	Rotatable Reverse Lip Rigid 32mm x 102mm (1 ¼" x 4")	48	45151ATT	Radiolucent Concave Lateral with Teeth 22mm x 70mm (⅞" x 2 ¾")	51
45098E	Rotatable Reverse Lip Rigid 32mm x 152mm (1 ¼" x 6")	48	45152	Swivel Richardson 51mm x 152mm (2" x 6")	34
45098RE	Rotatable Reverse Lip Rigid 32mm x 152mm (1 ¼" x 6")	48	45152A	Concave Lateral 25mm x 85mm (1" x 3 ⅜")	50
45099E	Rotatable Reverse Lip Rigid 32mm x 203mm (1 ¼" x 8")	48	45152AT	Radiolucent Concave Lateral 22mm x 80mm (1" x 3 ⅞")	51
45099RE	Rotatable Reverse Lip Rigid 32mm x 203mm (1 ¼" x 8")	48	45152ATT	Radiolucent Concave Lateral with Teeth 22mm x 80mm (1" x 3 ⅞")	51
45112	Swivel Malleable 51mm x 165mm (2" x 6 ½")	31	45153	Swivel Richardson 51mm x 178mm (2" x 7")	34
45113	Swivel Malleable 76mm x 165mm (3" x 6 ½")	31	45153A	Concave Medial 32mm x 70mm (1 ¼" x 2 ¾")	56
45114	Swivel Malleable 102mm x 165mm (4" x 6 ½")	31	45154A	Concave Medial 38mm x 70mm (1 ½" x 2 ¾")	56
45115	Swivel Malleable 13mm x 64mm (½" x 2 ½")	30	45155A	Concave Medial 38mm x 85mm (1 ½" x 3 ⅜")	56
45116	Swivel Malleable 13mm x 152mm (½" x 6")	30	45156A	Concave Lateral 25mm x 30mm (1" x 1 ¼")	50
45117	Swivel Malleable 13mm x 203mm (½" x 8")	30	45157A	Concave Lateral 25mm x 40mm (1" x 1 ½")	50
45118	Swivel Malleable 25mm x 203mm (1" x 8")	30	45158	Swivel Mini-Open Richardson 51mm x 152mm (2" x 6")	33
45118A	Swivel Malleable 25mm x 254mm (1" x 10")	30	45158A	Concave Lateral 25mm x 60mm (1" x 2 ¾")	50
45118B	Swivel Malleable 38mm x 254mm (1 ½" x 10")	30	45159A	Concave Lateral 25mm x 70mm (1" x 2 ¾")	50
45118P	Swivel Malleable 25mm x 102mm (1" x 4")	30	45160	Swivel Harrington 64mm x 152mm (2 ½" x 6")	27
45118R	Swivel Malleable Renal Vein 25mm x 203mm (1" x 8")	38	45160A	Concave Longitudinal 23mm x 30mm (1" x 1 ¼")	53
45119	Swivel Malleable 25mm x 127mm (1" x 5")	30	45160AT	Radiolucent Concave Longitudinal 23mm x 30mm (1" x 1 ¼")	53
45119A	Swivel Malleable 25mm x 152mm (1" x 6")	31	45161	Swivel Harrington 114mm x 254mm (4 ½" x 10")	27
45119B	Swivel Concave Renal Vein 25mm x 140mm (1" x 5 ½")	37	45161A	Concave Longitudinal 23mm x 50mm (1" x 2")	53
45119RE	Swivel Radiolucent Renal Vein 25mm x 140mm (1" x 5 ½")	37	45161AT	Radiolucent Concave Longitudinal 23mm x 40mm (1" x 1 ½")	53
45120	Swivel Malleable 38mm x 102mm (1 ½" x 4")	31	45162A	Concave Longitudinal 23mm x 70mm (1" x 2 ¾")	53
45120A	Swivel Malleable 38mm x 127mm (1 ½" x 5")	31	45162AT	Radiolucent Concave Longitudinal 23mm x 50mm (1" x 2")	53
45120B	Swivel Malleable 38mm x 152mm (1 ½" x 6")	31	45162P	Swivel Harrington 38mm x 127mm (1 ½" x 5")	27
45120RE	Swivel Radiolucent Renal Vein 25mm x 191mm (1" x 7 ½")	37	45163	Swivel Harrington 64mm x 203mm (2 ½" x 8")	27
45121RE	Swivel Radiolucent Renal Vein 25mm x 254mm (1" x 10")	37	45163AT	Radiolucent Concave Longitudinal 23mm x 60mm (1" x 2 ¾")	53
45122	Swivel Malleable 51mm x 203mm (2" x 8")	30	45164	Swivel Harrington 64mm x 254mm (2 ½" x 10")	27
45122A	Swivel Malleable 51mm x 254mm (2" x 10")	30	45164AT	Radiolucent Concave Longitudinal 23mm x 70mm (1" x 2 ¾")	53
45122P	Swivel Malleable 51mm x 102mm (2" x 4")	30	45165AT	Radiolucent Concave Longitudinal 23mm x 80mm (1" x 3 ⅜")	53
45123	Swivel Malleable 38mm x 203mm (1 ½" x 8")	30	45166AT	Radiolucent Concave Longitudinal 23mm x 90mm (1" x 3 ½")	53
45124	Swivel Malleable 57mm x 254mm (2 ½" x 10")	30	45167AT	Radiolucent Concave Longitudinal 23mm x 100mm (1" x 3 ⅞")	53
45124A	Swivel Malleable 76mm x 254mm (3" x 10")	30	45168AT	Radiolucent Concave Longitudinal 23mm x 110mm (1" x 4 ¼")	53
45124B	Swivel Malleable 102mm x 254mm (4" x 10")	30	45170	Swivel Deaver 25mm x 127mm (1" x 5")	36
45124C	Swivel Aortic Aneurysm Large 457mm x 229mm (18" x 9")	40	45172	Swivel Deaver 38mm x 127mm (1 ½" x 5")	36
45124D	Swivel Aortic Aneurysm Small 305mm x 216mm (12" x 8 ½")	40	45174	Swivel Deaver 51mm x 127mm (2" x 5")	36
45125	Swivel Richardson 19mm x 25mm (¾" x 1")	35	45175	Swivel Deaver 64mm x 127mm (2 ½" x 5")	36
45125A	Swivel Richardson 19mm x 51mm (¾" x 2")	35	45176	Swivel Deaver 76mm x 127mm (3" x 5")	36
45126	Swivel Richardson 25mm x 32mm (1" x 1 ¼")	35	45180	Swivel Heany 25mm x 102mm (1" x 4")	42
45127	Swivel Kelly 38mm x 51mm (1 ½" x 2")	28	45182	Swivel Heany 25mm x 127mm (1" x 5")	42
45128	Swivel Kelly 51mm x 51mm (2" x 2")	28	45183	Swivel Allison Lung 229mm (9")	40
45129	Swivel Kelly 64mm x 76mm (2 ½" x 3")	28	45185	Swivel Hill-Ferguson 32mm x 76mm (1 ¼" x 3")	45
45130	Swivel Kelly 76mm x 89mm (3" x 3 ½")	28	45186	Swivel Hill-Ferguson 22mm x 64mm (⅞" x 2 ½")	45
45134	Swivel Fence 102mm x 102mm (4" x 4")	29	45187	Swivel Rectal 13mm x 117mm (½" x 4 ½")	45
45135	Swivel Fence 102mm x 127mm (4" x 5")	29	45187A	Swivel St Marks [no lip] 64mm x 178mm (2 ½" x 7")	35
45136	Swivel Fence 102mm x 152mm (4" x 6")	29	45188	Swivel St Marks [with lip] 64mm x 178mm (2 ½" x 7")	35

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
45189	Swivel St Marks [with lip] 64mm x 140mm (2 1/2" x 5 1/2").....	35	45616A	Meyerding 51mm x 125mm (2" x 5").....	57
45190	Swivel Double Bladed 76mm x 114mm (3" x 4 1/2").....	42	45616AT	Meyerding with Teeth 51mm x 125mm (2" x 5")	57
45191	Swivel Malleable Finger 152mm (6").....	31	45620A	Army Navy 13mm x 55mm (1/2" x 2 1/8").....	56
45192	Swivel Bladder 35mm x 89mm (1 3/8" x 3 1/2").....	43	45620Q	Fixed Army Navy 13mm x 55mm (1/2" x 2 1/8").....	56
45194RE	Swivel Reverse Lip Rigid 25mm x 102mm (1" x 4")	48	45621A	Army Navy 13mm x 65mm (1/2" x 2 3/8").....	56
45195E	Swivel Reverse Lip Rigid 25mm x 152mm (1" x 6")	48	45621Q	Fixed Army Navy 13mm x 65mm (1/2" x 2 1/2").....	56
45195RE	Swivel Reverse Lip Rigid 25mm x 152mm (1" x 6")	48	45622A	Army Navy 13mm x 75mm (1/2" x 3 1/8").....	56
45196RE	Swivel Reverse Lip Rigid 25mm x 203mm (1" x 8")	48	45622Q	Fixed Army Navy 13mm x 75mm (1/2" x 3 1/8").....	56
45197RE	Swivel Reverse Lip Rigid 32mm x 102mm (1 1/4" x 4")	48	45624AT	Radiolucent MIS Tapered 15mm x 23mm x 35mm (5/8" x 1" x 1 3/8")	55
45198RE	Swivel Reverse Lip Rigid 32mm x 152mm (1 1/4" x 6")	48	45625AT	Radiolucent MIS Tapered 15mm x 23mm x 40mm (5/8" x 1" x 1 1/2")	55
45199RE	Swivel Reverse Lip Rigid 32mm x 203mm (1 1/4" x 8")	48	45626AT	Radiolucent MIS Tapered 15mm x 23mm x 45mm (5/8" x 1" x 1 3/4")	55
45200RE	Swivel Reverse Lip Rigid 25mm x 254mm (1" x 10")	48	45627AT	Radiolucent MIS Tapered 15mm x 23mm x 23mm x 50mm (5/8" x 1" x 2")	55
45201RE	Swivel Reverse Lip Rigid 32mm x 254mm (1 1/4" x 10")	48	45628AT	Radiolucent MIS Tapered 15mm x 23mm x 55mm (5/8" x 1" x 2 1/8").....	55
45210	Swivel Volkman Sharp 29mm x 19mm (1 1/8" x 3/4).....	39	45629AT	Radiolucent MIS Tapered 15mm x 23mm x 60mm (5/8" x 1" x 2 3/8").....	55
45211	Swivel Volkman Blunt 29mm x 19mm (1 1/8" x 3/4).....	39	45630AT	Radiolucent MIS Tapered 15mm x 23mm x 65mm (5/8" x 1" x 2 1/2")	55
45212	Swivel Volkman Sharp 19mm x 19mm (3/4" x 3/4").....	39	45631AT	Radiolucent MIS Tapered 15mm x 23mm x 70mm (5/8" x 1" x 2 3/4")	55
45212A	Concave Medial 13mm x 50mm (1/2" x 2")	56	45632AT	Radiolucent MIS Tapered 15mm x 23mm x 75mm (5/8" x 1" x 3")	55
45213	Swivel Volkman Blunt 16mm x 19mm (5/8" x 3/4").....	39	45640A	Concave Lateral 18mm x 70mm (5/8" x 2 3/4").....	50
45214	Swivel Murphy Sharp 25mm x 19mm (1" x 3/4").....	39	45641A	Concave Tapered 45mm x 50mm (1 3/4" x 2").....	52
45214A	Concave Medial 13mm x 75mm (1/2" x 3")	56	45641AT	Radiolucent Concave Tapered 51mm x 50mm (2" x 2").....	52
45215	Swivel Murphy Blunt 25mm x 22mm (1" x 7/8")	39	45641ATT	Radiolucent Concave Tapered with Teeth 51mm x 50mm (2" x 2").....	53
45216	Swivel Murphy Sharp 19mm x 22mm (3/4" x 7/8")	39	45642A	Concave Tapered 45mm x 65mm (1 3/4" x 2 1/2").....	52
45216A	Concave Medial 13mm x 100mm (1/2" x 4")	56	45642AT	Radiolucent Concave Tapered 51mm x 65mm (2" x 2 1/2")	52
45217	Swivel Murphy Blunt 19mm x 22mm (3/4" x 7/8").....	39	45642ATT	Radiolucent Concave Tapered with Teeth 51mm x 65mm (2" x 2 1/2").....	53
45222A	Concave Medial 19mm x 50mm (3/4" x 2").....	56	45643A	Concave Tapered 19mm x 35mm (3/4" x 1 3/8").....	52
45224A	Concave Medial 19mm x 75mm (3/4" x 3")	56	45643AT	Radiolucent Concave Tapered 23mm x 35mm (1" x 1 3/8").....	52
45226A	Concave Medial 19mm x 100mm (3/4" x 4")	56	45643ATT	Radiolucent Concave Tapered with Teeth 23mm x 35mm (1" x 1 3/8")	52
45232A	Concave Medial 25mm x 50mm (1" x 2").....	56	45644AT	Radiolucent Concave Tapered 23mm x 40mm (1" x 1 1/2").....	52
45234A	Concave Medial 25mm x 75mm (1" x 3")	56	45644ATT	Radiolucent Concave Tapered with Teeth 23mm x 40mm (1" x 1 1/2").....	52
45236A	Concave Medial 25mm x 100mm (1" x 4")	56	45645A	Concave Tapered 19mm x 45mm (3/4" x 1 3/8").....	52
45300	Femoral Neck Elevator 14"	67	45645AT	Radiolucent Concave Tapered 23mm x 45mm (1" x 1 3/8").....	52
45301	Cobra Hip Blade 11"	67	45645ATT	Radiolucent Concave Tapered with Teeth 23mm x 45mm (1" x 1 3/8").....	52
45302	Curved Cobra Hip Blade 12 1/2"	67	45646A	Concave Tapered 19mm x 50mm (3/4" x 2").....	52
45302C	Curved Cobra Hip Blade 12 1/2" [Coated]	67	45646AT	Radiolucent Concave Tapered 23mm x 50mm (1" x 2")	52
45305	Hohmann Hip Blade 11"	67	45646ATT	Radiolucent Concave Tapered with Teeth 23mm x 50mm (1" x 2").....	52
45306	Curved Hohmann Narrow Hip Blade 12"	67	45647AT	Radiolucent Concave Tapered 23mm x 55mm (1" x 2 1/8").....	52
45306C	Curved Hohmann Narrow Hip Blade 12" [Coated].....	67	45647ATT	Radiolucent Concave Tapered with Teeth 23mm x 55mm (1" x 2 1/8")	52
45307	Curved Hohmann Wide Hip Blade 12"	67	45648AT	Radiolucent Concave Tapered 23mm x 60mm (1" x 2 3/8").....	52
45308	Curved Dual Prong Hip Blade 12"	67	45648ATT	Radiolucent Concave Tapered with Teeth 23mm x 60mm (1" x 2 3/8")	52
45309	Dual Prong Hip Blade 11"	67	45649AT	Radiolucent Concave Tapered 23mm x 65mm (1" x 2 1/2").....	52
45600A	Meyerding 16mm x 75mm (5/8" x 3").....	57	45649ATT	Radiolucent Concave Tapered with Teeth 23mm x 65mm (1" x 2 1/2").....	52
45600AT	Meyerding with Teeth 16mm x 75mm (5/8" x 3")	57	45650A	Nerve Root Retractor 9mm x 114mm (3/8" x 4 1/2").....	59
45601A	Meyerding 16mm x 100mm (5/8" x 4").....	57	45650AT	Radiolucent Concave Tapered 23mm x 70mm (1" x 2 3/4").....	52
45601AT	Meyerding with Teeth 16mm x 100mm (5/8" x 4")	57	45650ATT	Radiolucent Concave Tapered with Teeth 23mm x 70mm (1" x 2 3/4").....	52
45602A	Meyerding 25mm x 100mm (1" x 4")	57	45651A	Nerve Root Retractor 13mm x 114mm (1/2" x 4 1/2")	59
45602AT	Meyerding with Teeth 25mm x 100mm (1" x 4")	57	45656A	Malleable 25mm x 110mm (1" x 4 1/2")	30
45603A	Meyerding 25mm x 125mm (1" x 5")	57	45656AT	Radiolucent Malleable 25mm x 110mm (1" x 4 1/2")	30
45603AT	Meyerding with Teeth 25mm x 125mm (1" x 5")	57	45657A	Malleable 13mm x 90mm (1/2" x 3 1/2").....	30
45604A	Peress 60mm x 102mm (2 3/8" x 4").....	59	45657AT	Radiolucent Malleable 13mm x 90mm (1/2" x 3 1/2")	30
45605A	Peress 60mm x 127mm (2 3/8" x 5").....	59	45660	Rotatable Double Angled 25mm x 76mm (1" x 3").....	41
45606A	Peress 60mm x 152mm (2 3/8" x 6").....	59	45660Q	Swivel Double Angled 25mm x 76mm (1" x 3").....	41
45607A	Peress 35mm x 102mm (1 3/8" x 4").....	59	45661A	Rotatable Double Angled 25mm x 64mm (1" x 2 1/2").....	41
45608A	Peress 35mm x 127mm (1 3/8" x 5").....	59	45661Q	Swivel Double Angled 25mm x 64mm (1" x 2 1/2").....	41
45609A	Peress 35mm x 152mm (1 3/8" x 6").....	59	45662Q	Swivel Double Angled 13mm x 76mm (1/2" x 3").....	41
45610A	Meyerding 25mm x 60mm (1" x 2 3/8").....	57	45663Q	Swivel Double Angled 13mm x 64mm (1/2" x 2 1/2").....	41
45610AT	Meyerding with Teeth 25mm x 60mm (1" x 2 3/8").....	57	45664Q	Swivel Double Angled Bulb Retractor 25mm x 64mm (1" x 2 1/2").....	41
45611A	Meyerding 25mm x 75mm (1" x 3").....	57	45665Q	Swivel Prostate Retractor 8mm (3/8").....	44
45611AT	Meyerding with Teeth 25mm x 75mm (1" x 3")	57	45666AT	Radiolucent Odontoid 20mm x 90mm (3/4" x 3 1/2")	54
45612A	Meyerding 51mm x 60mm (2" x 2 3/8").....	57	45670ATTC	Radiolucent Concave Posterior Lumbar 23mm x 80mm (1" x 3 1/8").....	58
45612AET	Radiolucent Meyerding 51mm x 60mm (2" x 2 3/8").....	57	45671ATTC	Radiolucent Concave Posterior Lumbar 45mm x 80mm (1 3/4" x 3 1/8").....	58
45612AETT	Radiolucent Meyerding with Teeth 51mm x 60mm (2" x 2 3/8").....	57	45672ATTC	Radiolucent Concave Posterior Lumbar 60mm x 80mm (2 3/8" x 3 1/8").....	58
45612AT	Meyerding with Teeth 51mm x 60mm (2" x 2 3/8").....	57	45673ATTC	Radiolucent Concave Posterior Lumbar 23mm x 100mm (1" x 4").....	58
45613A	Meyerding 51mm x 75mm (2" x 3").....	57	45674ATTC	Radiolucent Concave Posterior Lumbar 45mm x 100mm (1 3/4" x 4").....	58
45613AT	Meyerding with Teeth 51mm x 75mm (2" x 3")	57	45675ATTC	Radiolucent Concave Posterior Lumbar 60mm x 100mm (2 3/8" x 4").....	58
45614A	Meyerding 16mm x 51mm (5/8" x 2").....	57	45676ATTC	Radiolucent Concave Posterior Lumbar 23mm x 120mm (1" x 4 3/8").....	58
45614AT	Meyerding with Teeth 16mm x 51mm (5/8" x 2").....	57	45677ATTC	Radiolucent Concave Posterior Lumbar 45mm x 120mm (1 3/4" x 4 3/8").....	58
45615A	Meyerding 51mm x 100mm (2" x 4").....	57	45678ATTC	Radiolucent Concave Posterior Lumbar 60mm x 120mm (2 3/8" x 4 3/8").....	58
45615AT	Meyerding with Teeth 51mm x 100mm (2" x 4")	57	45679ATTC	Radiolucent Concave Posterior Lumbar 23mm x 150mm (1" x 5 7/8").....	58
			45680ATTC	Radiolucent Concave Posterior Lumbar 45mm x 150mm (1 3/4" x 5 7/8").....	58

PRODUCT INDEX

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
45681ATTC	Radiolucent Concave Posterior Lumbar 60mm x 150mm (2 3/8" x 5 7/8")	58	46117ET	Radiolucent Malleable 13mm x 203mm (1/2" x 8")	30
45682ATTC	Radiolucent Concave Posterior Lumbar 23mm x 180mm (1" x 7")	58	46118	Malleable 25mm x 203mm (1" x 8")	30
45683ATTC	Radiolucent Concave Posterior Lumbar 45mm x 180mm (1 3/4" x 7")	58	46118A	Malleable 25mm x 254mm (1" x 10")	30
45684ATTC	Radiolucent Concave Posterior Lumbar 60mm x 180mm (2 3/8" x 7")	58	46118B	Malleable 38mm x 254mm (1 1/2" x 10")	30
45750A	Midline R 70mm x 83mm (2 3/4" x 3 1/4")	55	46118C	Malleable 25mm x 76mm (1" x 3")	31
45750AT	Radiolucent Midline R 70mm x 83mm (2 3/4" x 3 1/4")	55	46118CET	Radiolucent Renal Vein 25mm x 152mm (1" x 6")	37
45751A	Midline L 70mm x 83mm (2 3/4" x 3 1/4")	55	46118DET	Radiolucent Renal Vein 51mm x 152mm (2" x 6")	37
45751AT	Radiolucent Midline L 70mm x 83mm (2 3/4" x 3 1/4")	55	46118EET	Radiolucent Rigid Renal Vein 51mm x 203mm (2" x 8")	37
45752A	Nerve Root Retractor 90° 5mm x 123mm (1/4" x 4 7/8")	59	46118ET	Radiolucent Malleable 25mm x 203mm (1" x 8")	30
45753A	Nerve Root Retractor 45° 8mm x 145mm (3/8" x 5 3/4")	59	46118P	Malleable 25mm x 102mm (1" x 4")	30
45753AM	Prostate Retractor 8mm (3/8")	44	46118PET	Radiolucent Malleable 25mm x 102mm (1" x 4")	30
45754A	Midline R 57mm x 83mm (2 1/4" x 3 1/4")	55	46118R	Malleable Renal Vein 25mm x 203mm (1" x 8")	38
45755A	Midline L 57mm x 83mm (2 1/4" x 3 1/4")	55	46118RET	Radiolucent Rigid Renal Vein 25mm x 203mm (1" x 8")	37
45760A	Prong Retractor 75mm (3")	58	46119	Malleable 25mm x 127mm (1" x 5")	30
45761A	Prong Retractor 100mm (4")	58	46119A	Malleable 25mm x 152mm (1" x 6")	31
45831	Dilator #1 - 7mm	69	46119B	Concave Renal Vein 25mm x 140mm (1" x 5 1/2")	37
45832	Dilator #2 - 12mm	69	46119BET	Radiolucent Malleable Renal Vein 25mm x 140mm (1" x 5 1/2")	37
45833	Dilator #3 - 17mm	69	46119C	Concave Renal Vein 25mm x 191mm (1" x 7 1/2")	37
45834	Dilator #4 - 21mm	69	46119CET	Radiolucent Malleable Renal Vein 25mm x 190mm (1" x 7 1/2")	37
45924L	Caudal/Cephalad Left 22mm x 40mm (7/8" x 1 1/2")	70	46119RET	Radiolucent Rigid Renal Vein 25mm x 140mm (1" x 5 1/2")	37
45924R	Caudal/Cephalad Right 22mm x 40mm (7/8" x 1 1/2")	70	46120	Malleable 38mm x 102mm (1 1/2" x 4")	31
45925L	Caudal/Cephalad Left 22mm x 50mm (7/8" x 2")	70	46120A	Malleable 38mm x 127mm (1 1/2" x 5")	31
45925R	Caudal/Cephalad Right 22mm x 50mm (7/8" x 2")	70	46120B	Malleable 38mm x 152mm (1 1/2" x 6")	31
45926L	Caudal/Cephalad Left 22mm x 60mm (7/8" x 2 3/8")	70	46120RET	Radiolucent Rigid Renal Vein 25mm x 190mm (1" x 7 1/2")	37
45926R	Caudal/Cephalad Right 22mm x 60mm (7/8" x 2 3/8")	70	46121RET	Radiolucent Rigid Renal Vein 25mm x 255mm (1" x 10")	37
45927L	Caudal/Cephalad Left 22mm x 70mm (7/8" x 2 3/4")	70	46122	Malleable 51mm x 203mm (2" x 8")	30
45927R	Caudal/Cephalad Right 22mm x 70mm (7/8" x 2 3/4")	70	46122A	Malleable 51mm x 254mm (2" x 10")	30
45928L	Caudal/Cephalad Left 22mm x 80mm (7/8" x 3 1/8")	70	46122AET	Radiolucent Malleable 51mm x 254mm (2" x 10")	30
45928R	Caudal/Cephalad Right 22mm x 80mm (7/8" x 3 1/8")	70	46122ET	Radiolucent Malleable 51mm x 203mm (2" x 8")	31
45929L	Caudal/Cephalad Left 22mm x 90mm (7/8" x 3 1/2")	70	46122P	Malleable 51mm x 102mm (2" x 4")	30
45929R	Caudal/Cephalad Right 22mm x 90mm (7/8" x 3 1/2")	70	46122RET	Radiolucent Rigid Renal Vein 25mm x 305mm (1" x 12")	37
45930L	Caudal/Cephalad Left 22mm x 100mm (7/8" x 4")	70	46122V	Malleable Slotted 51mm x 203mm (2" x 8")	32
45930R	Caudal/Cephalad Right 22mm x 100mm (7/8" x 4")	70	46123	Malleable 38mm x 203mm (1 1/2" x 8")	30
45931L	Caudal/Cephalad Left 22mm x 110mm (7/8" x 4 3/8")	70	46124	Malleable 64mm x 254mm (2 1/2" x 10")	30
45931R	Caudal/Cephalad Right 22mm x 110mm (7/8" x 4 3/8")	70	46124A	Malleable 76mm x 254mm (3" x 10")	30
45932L	Caudal/Cephalad Left 22mm x 120mm (7/8" x 4 3/4")	70	46124AET	Radiolucent Malleable 76mm x 254mm (3" x 10")	31
45932R	Caudal/Cephalad Right 22mm x 120mm (7/8" x 4 3/4")	70	46124AV	Malleable Slotted 76mm x 203mm (3" x 8")	32
45973	PLA Midline 22mm x 30mm (7/8" x 1 1/8")	71	46124B	Malleable 102mm x 254mm (4" x 10")	30
45974	PLA Midline 22mm x 40mm (7/8" x 1 1/2")	71	46124BET	Radiolucent Malleable 102mm x 254mm (4" x 10")	31
45975	PLA Midline 22mm x 50mm (7/8" x 2")	71	46124C	Aortic Aneurysm Large 457mm x 229mm (18" x 9")	40
45976	PLA Midline 22mm x 60mm (7/8" x 2 3/8")	71	46124D	Aortic Aneurysm Small 305mm x 216mm (12" x 8 1/2")	40
45977	PLA Midline 22mm x 70mm (7/8" x 2 3/4")	71	46124ET	Radiolucent Malleable 64mm x 254mm (2 1/2" x 10")	31
45978	PLA Midline 22mm x 80mm (7/8" x 3 1/8")	71	46125	Richardson 19mm x 25mm (3/4" x 1")	35
45979	PLA Midline 22mm x 90mm (7/8" x 3 1/2")	71	46125A	Richardson 19mm x 50mm (3/4" x 2")	35
45980	PLA Midline 22mm x 100mm (7/8" x 4")	71	46126	Richardson 25mm x 32mm (1" x 1 1/4")	35
45981	PLA Midline 22mm x 110mm (7/8" x 4 3/8")	71	46127	Kelly 38mm x 51mm (1 1/2" x 2")	28
45982	PLA Midline 22mm x 120mm (7/8" x 4 3/4")	71	46127ET	Radiolucent Kelly 38mm x 51mm (1 1/2" x 2")	28
45983	PLA Midline 22mm x 130mm (7/8" x 5 1/8")	71	46128	Kelly 51mm x 64mm (2" x 2 1/2")	28
45984	PLA Midline 27mm x 30mm (1" x 1 1/8")	71	46128ET	Radiolucent Kelly 51mm x 64mm (2" x 2 1/2")	28
45985	PLA Midline 27mm x 40mm (1" x 1 1/2")	71	46129	Kelly 64mm x 76mm (2 1/2" x 3")	28
45986	PLA Midline 27mm x 50mm (1" x 2")	71	46129ET	Radiolucent Kelly 64mm x 76mm (2 1/2" x 3")	28
45987	PLA Midline 27mm x 60mm (1" x 2 3/8")	71	46130	Kelly 76mm x 89mm (3" x 3 1/2")	28
45988	PLA Midline 27mm x 70mm (1" x 2 3/4")	71	46130ET	Radiolucent Kelly 76mm x 89mm (3" x 3 1/2")	28
45989	PLA Midline 27mm x 80mm (1" x 3 1/8")	71	46134	Fence 102mm x 102mm (4" x 4")	29
45990	PLA Midline 27mm x 90mm (1" x 3 1/2")	71	46134ET	Radiolucent Fence 102mm x 102mm (4" x 4")	29
45991	PLA Midline 27mm x 100mm (1" x 4")	71	46135	Fence 102mm x 127mm (4" x 5")	29
45992	PLA Midline 27mm x 110mm (1" x 4 3/8")	71	46135ET	Radiolucent Fence 102mm x 127mm (4" x 5")	29
45993	PLA Midline 27mm x 120mm (1" x 4 3/4")	71	46136	Fence 102mm x 152mm (4" x 6")	29
45994	PLA Midline 27mm x 130mm (1" x 5 1/8")	71	46136A	Fence 102mm x 203mm (4" x 8")	29
46103	Saldanha Arched 76mm x 140mm (3" x 5 1/2")	73	46136AET	Radiolucent Fence 102mm x 203mm (4" x 8")	29
46112	Malleable 51mm x 165mm (2" x 6 1/2")	31	46136ET	Radiolucent Fence 102mm x 152mm (4" x 6")	29
46113	Malleable 76mm x 165mm (3" x 6 1/2")	31	46137	Holtgrewe 60mm x 152mm (2 3/8" x 6")	43
46114	Malleable 102mm x 165mm (4" x 6 1/2")	31	46139	"Y" Bladder 81mm x 216mm (3 1/8" x 8 1/2")	43
46115	Malleable 13mm x 64mm (1/2" x 2 1/2")	30	46140	Balfour [with lips] 70mm x 73mm (2 3/4" x 3")	25
46115ET	Radiolucent Malleable 13mm x 64mm (1/2" x 2 1/2")	30	46140C	Balfour [small lips] 70mm x 73mm (2 3/4" x 3")	25
46116	Malleable 13mm x 152mm (1/2" x 6")	30	46140CEB	Radiolucent Balfour [small lips] 65mm x 72mm (2 1/2" x 2 3/4") [Anodized]	25
46116ET	Radiolucent Malleable 13mm x 152mm (1/2" x 6")	30	46140CET	Radiolucent Balfour [small lips] 65mm x 72mm (2 1/2" x 2 3/4")	25
46117	Malleable 13mm x 203mm (1/2" x 8")	30	46140EB	Radiolucent Balfour [with lips] 65mm x 72mm (2 1/2" x 2 3/4") [Anodized]	25
			46140ET	Radiolucent Balfour [with lips] 65mm x 72mm (2 1/2" x 2 3/4")	25

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
46140P	Balfour [with lips] 76mm x 51mm (3" x 2").....	25	46187A	St Marks [no lip] 64mm x 178mm (2 ½" x 7").....	35
46141	Balfour [with lips] 83mm x 73mm (3 ¼" x 3").....	25	46188	St Marks [with lip] 64mm x 178mm (2 ½" x 7")	35
46141C	Balfour [small lips] 83mm x 73mm (3 ¼" x 3")	25	46188ET	Radiolucent St Marks [with lip] 64mm x 178mm (2 ½" x 7")	35
46141CEB	Radiolucent Balfour [small lips] 83mm x 70mm (3 ¼" x 2 ¾") [Anodized].....	25	46189	St Marks [with lip] 64mm x 140mm (2 ½" x 5 ½")	35
46141CET	Radiolucent Balfour [small lips] 83mm x 70mm (3 ¼" x 2 ¾")	25	46189ET	Radiolucent St Marks [with lip] 64mm x 140mm (2 ½" x 5 ½")	35
46141D	Balfour [with lips] 83mm x 123mm (3 ¼" x 5").....	25	46190	Double Bladed 76mm x 114mm (3" x 4 ½").....	42
46141DEB	Radiolucent Balfour [with lips] 83mm x 127mm (3 ¼" x 5") [Anodized].....	25	46191	Malleable Finger 152mm (6").....	31
46141EB	Radiolucent Balfour [with lips] 83mm x 70mm (3 ¼" x 2 ¾") [Anodized].....	25	46191ET	Radiolucent Malleable Finger 152mm (6").....	31
46141ET	Radiolucent Balfour [with lips] 83mm x 70mm (3 ¼" x 2 ¾").....	25	46191L	Malleable Finger 203mm (8").....	31
46141F	Balfour [with lips] 83mm x 175mm (3 ¼" x 6 ¾").....	25	46191P	Malleable Finger 102mm (4").....	31
46142	Balfour [with lips] 32mm x 32mm (1 ¼" x 1 ¼").....	25	46191SP	Malleable Finger 76mm (3").....	31
46142EB	Radiolucent Balfour [with lips] 32mm x 32mm (1 ¼" x 1 ¼") [Anodized].....	25	46192	Bladder Blade 35mm x 89mm (1 ¾" x 3 ½").....	43
46142P	Balfour [with lips] 51mm x 54mm (2" x 2").....	25	46192ET	Reverse Lip Malleable 25mm x 200mm (1" x 7 ¾").....	47
46142PEB	Radiolucent Balfour [with lips] 51mm x 51mm (2" x 2") [Anodized].....	25	46193ET	Reverse Lip Malleable 32mm x 200mm (1 ¼" x 7 ¾").....	47
46143	Fenestrated 44mm x 89mm (1 ¾" x 3 ½").....	38	46194ET	Reverse Lip Rigid 25mm x 100mm (1" x 4")	48
46144	Fenestrated 48mm x 127mm (1 ¾" x 5").....	38	46194RET	Reverse Lip Rigid 25mm x 100mm (1" x 4")	48
46145	Splanchnich 51mm x 114mm (2" x 4 ½")	29	46195ET	Reverse Lip Rigid 25mm x 150mm (1" x 5 ¾")	48
46145ET	Radiolucent Splanchnich 51mm x 115mm (2" x 4 ½").....	29	46195RET	Reverse Lip Rigid 25mm x 150mm (1" x 5 ¾")	48
46146	Splanchnich 51mm x 152mm (2" x 6").....	29	46195TET	Reverse Lip Tapered 25mm x 150mm (1" x 5 ¾")	48
46146ET	Radiolucent Splanchnich 51mm x 152mm (2" x 6")	29	46196CET	Reverse Lip Concave 25mm x 200mm (1" x 7 ¾").....	48
46147ET	Radiolucent Splanchnich 51mm x 203mm (2" x 8")	29	46196ET	Reverse Lip Rigid 25mm x 200mm (1" x 7 ¾")	48
46148ET	Radiolucent Splanchnich 25mm x 152mm (1" x 6")	29	46196RET	Reverse Lip Rigid 25mm x 200mm (1" x 7 ¾")	48
46149	Richardson 44mm x 51mm (1 ¾" x 2").....	34	46196TET	Reverse Lip Tapered 25mm x 200mm (1" x 7 ¾")	48
46149ET	Radiolucent Splanchnich 25mm x 203mm (1" x 8")	29	46197ET	Reverse Lip Rigid 32mm x 100mm (1 ¼" x 4")	48
46150	Richardson 51mm x 102mm (2" x 4").....	34	46197RET	Reverse Lip Rigid 32mm x 100mm (1 ¼" x 4")	48
46151	Richardson 51mm x 127mm (2" x 5").....	34	46198ET	Reverse Lip Rigid 32mm x 150mm (1 ¼" x 5 ¾")	48
46152	Richardson 51mm x 152mm (2" x 6").....	34	46198RET	Reverse Lip Rigid 32mm x 150mm (1 ¼" x 5 ¾")	48
46153	Richardson 51mm x 178mm (2" x 7").....	34	46199CET	Reverse Lip Concave 32mm x 200mm (1 ¼" x 7 ¾").....	48
46153ET	Radiolucent Richardson 51mm x 178mm (2" x 7").....	34	46199ET	Reverse Lip Rigid 32mm x 200mm (1 ¼" x 7 ¾")	48
46154	Richardson 51mm x 76mm (2" x 3").....	34	46199RET	Reverse Lip Rigid 32mm x 200mm (1 ¼" x 7 ¾")	48
46155	Mini-Open Richardson 51mm x 76mm (2" x 3").....	33	46200CET	Reverse Lip Concave 25mm x 255mm (1" x 10")	48
46156	Mini-Open Richardson 51mm x 102mm (2" x 4")	33	46200ET	Reverse Lip Rigid 25mm x 255mm (1" x 10")	48
46157	Mini-Open Richardson 51mm x 127mm (2" x 5")	33	46200RET	Reverse Lip Rigid 25mm x 255mm (1" x 10")	48
46158	Mini-Open Richardson 51mm x 152mm (2" x 6")	33	46201CET	Reverse Lip Concave 32mm x 255mm (1 ¼" x 10")	48
46160	Harrington 64mm x 152mm (2 ½" x 6").....	27	46201ET	Reverse Lip Rigid 32mm x 255mm (1 ¼" x 10")	48
46160E	Radiolucent Harrington 64mm x 152mm (2 ½" x 6").....	27	46201RET	Reverse Lip Rigid 32mm x 255mm (1 ¼" x 10")	48
46160ET	Radiolucent Harrington 64mm x 152mm (2 ½" x 6").....	27	46202	Malleable Slotted 51mm x 165mm (2" x 6 ½").....	32
46160S	Malleable Harrington 64mm x 152mm (2 ½" x 6")	27	46203	Malleable Slotted 76mm x 165mm (3" x 6 ½").....	32
46161	Harrington 114mm x 254mm (4 ½" x 10").....	27	46210	Volkman Sharp 6 Prong 25mm x 13mm (1" x ½").....	39
46161M	Harrington 114mm x 203mm (4 ½" x 8").....	27	46211	Volkman Blunt 6 Prong 25mm x 13mm (1" x ½").....	39
46162P	Harrington 38mm x 127mm (1 ½" x 5").....	27	46212	Volkman Sharp 4 Prong 16mm x 13mm (5/8" x ½").....	39
46163	Harrington 64mm x 203mm (2 ½" x 8").....	27	46213	Volkman Blunt 4 Prong 16mm x 13mm (5/8" x ½")	39
46164	Harrington 64mm x 254mm (2 ½" x 10").....	27	46214	Murphy Sharp 4 Prong 25mm x 13mm (1" x ½")	39
46165	Harrington 64mm x 102mm (2 ½" x 4").....	27	46215	Murphy Blunt 4 Prong 25mm x 13mm (1" x ½")	39
46166	Harrington 38mm x 203mm (1 ½" x 8").....	27	46216	Murphy Sharp 3 Prong 19mm x 13mm (¾" x ½").....	39
46168B	Non-Slip Balfour L 77mm x 66mm (3" x 2 ½").....	26	46217	Murphy Blunt 3 Prong 19mm x 13mm (¾" x ½")	39
46168C	Non-Slip Balfour L 77mm x 94mm (3" x 3 ¾").....	26	46220	Mini-Open Malleable 25mm x 75mm (1" x 3").....	32
46169B	Non-Slip Balfour R 77mm x 66mm (3" x 2 ½").....	26	46221	Mini-Open Malleable 25mm x 100mm (1" x 4")	32
46169C	Non-Slip Balfour R 77mm x 94mm (3" x 3 ¾").....	26	46222	Mini-Open Malleable 25mm x 125mm (1" x 5")	32
46170	Deaver 25mm x 127mm (1" x 5").....	36	46223	Prostate Malleable 280mm x 40mm (11" x 1 ½")	44
46170ET	Radiolucent Deaver 25mm x 127mm (1" x 5").....	36	46224	Mini-Open Double Angled 51mm x 127mm (2" x 5")	33
46172	Deaver 38mm x 127mm (1 ½" x 5").....	36	46225	Mini-Open Double Angled 51mm x 152mm (2" x 6")	33
46172ET	Radiolucent Deaver 38mm x 127mm (1 ½" x 5").....	36	46226	Mini-Open Double Angled 51mm x 165mm (2" x 6 ½")	33
46174	Deaver 51mm x 127mm (2" x 5").....	36	46227	Mini-Open Deep Double Angled 51mm x 127mm (2" x 5")	33
46174ET	Radiolucent Deaver 51mm x 127mm (2" x 5").....	36	46241A	Wide Balfour [no lips] 83mm x 57mm (3 ¼" x 2 ¼").....	26
46175	Deaver 64mm x 127mm (2 ½" x 5").....	36	46241AEB	Radiolucent Wide Balfour [no lips] 83mm x 57mm (3 ¼" x 2 ¼") [Anodized].....	26
46175ET	Radiolucent Deaver 64mm x 127mm (2 ½" x 5").....	36	46241B	Wide Balfour [no lips] 100mm x 61mm (4" x 2 ¾").....	26
46176	Deaver 76mm x 127mm (3" x 5").....	36	46241BEB	Radiolucent Wide Balfour [no lips] 102mm x 57mm (4" x 2 ¾") [Anodized]	26
46176ET	Radiolucent Deaver 76mm x 127mm (3" x 5").....	36	46241C	Wide Balfour [no lips] 114mm x 63mm (4 ½" x 2 ½").....	26
46178	Deaver 102mm x 127mm (4" x 5").....	36	46241CEB	Radiolucent Wide Balfour [no lips] 114mm x 57mm (4 ½" x 2 ½") [Anodized].....	26
46178ET	Radiolucent Deaver 100mm x 127mm (4" x 5").....	36	46260RET	Reverse Lip Rigid 25mm x 100mm (1" x 4")	46
46180	Heany 25mm x 102mm (1" x 4").....	42	46260TET	Reverse Lip Tapered 25mm x 100mm (1" x 4").....	47
46182	Heany 25mm x 127mm (1" x 5").....	42	46261RET	Reverse Lip Rigid 25mm x 110mm (1" x 4 ¾")	46
46183	Allison Lung Retractor 229mm (9")	40	46262RET	Reverse Lip Rigid 25mm x 120mm (1" x 4 ¾")	46
46185	Hill-Ferguson 32mm x 76mm (1 ¼" x 3").....	45	46263RET	Reverse Lip Rigid 25mm x 130mm (1" x 5 ¼")	46
46186	Hill-Ferguson 25mm x 64mm (1" x 2 ½").....	45	46264RET	Reverse Lip Rigid 25mm x 140mm (1" x 5 ½")	46
46187	Rectal 13mm x 117mm (½" x 4 ½")	45	46265RET	Reverse Lip Rigid 25mm x 150mm (1" x 5 ¾")	46

PRODUCT INDEX

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
46265TET	Reverse Lip Tapered 25mm x 150mm (1" x 5 7/8")	47	46402	Cat's Paw 6mm x 38mm (1/4" x 1 1/2")	73
46266RET	Reverse Lip Rigid 25mm x 160mm (1" x 6 1/4")	46	46510AT	Radiolucent Concave Micro 14mm x 35mm (1/2" x 1 3/8")	54
46267RET	Reverse Lip Rigid 25mm x 170mm (1" x 6 3/4")	46	46510ATT	Radiolucent Concave Micro with Teeth 14mm x 35mm (1/2" x 1 3/8")	54
46268RET	Reverse Lip Rigid 25mm x 180mm (1" x 7")	46	46512AT	Radiolucent Concave Micro 14mm x 40mm (1/2" x 1 1/2")	54
46269RET	Reverse Lip Rigid 25mm x 190mm (1" x 7 1/2")	46	46512ATT	Radiolucent Concave Micro with Teeth 14mm x 40mm (1/2" x 1 1/2")	54
46270RET	Reverse Lip Rigid 25mm x 200mm (1" x 7 7/8")	46	46514AT	Radiolucent Concave Micro 14mm x 45mm (1/2" x 1 3/4")	54
46270TET	Reverse Lip Tapered 25mm x 200mm (1" x 7 7/8")	47	46514ATT	Radiolucent Concave Micro with Teeth 14mm x 45mm (1/2" x 1 3/4")	54
46271RET	Reverse Lip Rigid 25mm x 210mm (1" x 8 1/4")	46	46516AT	Radiolucent Concave Micro 14mm x 50mm (1/2" x 2")	54
46272RET	Reverse Lip Rigid 25mm x 220mm (1" x 8 5/8")	46	46516ATT	Radiolucent Concave Micro with Teeth 14mm x 50mm (1/2" x 2")	54
46273RET	Reverse Lip Rigid 25mm x 230mm (1" x 9")	46	46518AT	Radiolucent Concave Micro 14mm x 55mm (1/2" x 2 1/8")	54
46274RET	Reverse Lip Rigid 25mm x 240mm (1" x 9 1/2")	46	46518ATT	Radiolucent Concave Micro with Teeth 14mm x 55mm (1/2" x 2 1/8")	54
46275RET	Reverse Lip Rigid 25mm x 250mm (1" x 9 3/4")	46	46520AT	Radiolucent Concave Micro 14mm x 60mm (1/2" x 2 3/4")	54
46275TET	Reverse Lip Tapered 25mm x 250mm (1" x 9 3/4")	47	46520ATT	Radiolucent Concave Micro with Teeth 14mm x 60mm (1/2" x 2 3/4")	54
46280RET	Reverse Lip Rigid 32mm x 100mm (1 1/4" x 4")	46	46521	Depth Gauge Anterior Cervical 157.5mm (6 1/4")	21
46281RET	Reverse Lip Rigid 32mm x 110mm (1 3/4" x 4 3/8")	46	46523AT	Radiolucent Concave Micro 14mm x 65mm (1/2" x 2 1/2")	54
46282RET	Reverse Lip Rigid 32mm x 120mm (1 3/4" x 4 3/4")	46	46523ATT	Radiolucent Concave Micro with Teeth 14mm x 65mm (1/2" x 2 1/2")	54
46283RET	Reverse Lip Rigid 32mm x 130mm (1 3/4" x 5 1/8")	46	46524AT	Radiolucent Concave Micro 14mm x 70mm (1/2" x 2 3/4")	54
46284RET	Reverse Lip Rigid 32mm x 140mm (1 3/4" x 5 1/2")	46	46524ATT	Radiolucent Concave Micro with Teeth 14mm x 70mm (1/2" x 2 3/4")	54
46285RET	Reverse Lip Rigid 32mm x 150mm (1 3/4" x 5 7/8")	46	46526AT	Radiolucent Concave Micro 14mm x 75mm (1/2" x 3")	54
46286RET	Reverse Lip Rigid 32mm x 160mm (1 3/4" x 6 1/4")	46	46526ATT	Radiolucent Concave Micro with Teeth 14mm x 75mm (1/2" x 3")	54
46287RET	Reverse Lip Rigid 32mm x 170mm (1 3/4" x 6 5/8")	46	46528AT	Radiolucent Concave Micro 14mm x 80mm (1/2" x 3 1/8")	54
46288RET	Reverse Lip Rigid 32mm x 180mm (1 3/4" x 7")	46	46528ATT	Radiolucent Concave Micro with Teeth 14mm x 80mm (1/2" x 3 1/8")	54
46289RET	Reverse Lip Rigid 32mm x 190mm (1 3/4" x 7 1/2")	47	46558	Radiolucent Concave Renal Vein 25mm x 80mm (1" x 3 1/8")	49
46290RET	Reverse Lip Rigid 32mm x 200mm (1 3/4" x 7 7/8")	47	46559	Radiolucent Concave Renal Vein 25mm x 90mm (1" x 3 1/2")	49
46291RET	Reverse Lip Rigid 32mm x 210mm (1 1/4" x 8 1/4")	47	46560	Radiolucent Concave Renal Vein 25mm x 100mm (1" x 4")	49
46292RET	Reverse Lip Rigid 32mm x 220mm (1 1/4" x 8 5/8")	47	46561	Radiolucent Concave Renal Vein 25mm x 110mm (1" x 4 3/8")	49
46293RET	Reverse Lip Rigid 32mm x 230mm (1 1/4" x 9")	47	46562	Radiolucent Concave Renal Vein 25mm x 120mm (1" x 4 3/4")	49
46294RET	Reverse Lip Rigid 32mm x 240mm (1 1/4" x 9 1/2")	47	46563	Radiolucent Concave Renal Vein 25mm x 130mm (1" x 5 1/8")	49
46295RET	Reverse Lip Rigid 32mm x 250mm (1 1/4" x 9 3/4")	47	46564	Radiolucent Concave Renal Vein 25mm x 140mm (1" x 5 1/2")	49
46312AT	Radiolucent Concave Lateral 23mm x 35mm (1" x 1 3/8")	50	46565	Radiolucent Concave Renal Vein 25mm x 150mm (1" x 5 5/8")	49
46312ATT	Radiolucent Concave Lateral with Teeth 23mm x 35mm (1" x 1 3/8")	51	46566	Radiolucent Concave Renal Vein 25mm x 160mm (1" x 6 1/4")	49
46314AT	Radiolucent Concave Lateral 23mm x 40mm (1" x 1 1/2")	50	46567	Radiolucent Concave Renal Vein 25mm x 170mm (1" x 6 3/4")	49
46314ATT	Radiolucent Concave Lateral with Teeth 23mm x 40mm (1" x 1 1/2")	51	46568	Radiolucent Concave Renal Vein 25mm x 180mm (1" x 7")	49
46316AT	Radiolucent Concave Lateral 23mm x 45mm (1" x 1 3/4")	50	46569	Radiolucent Concave Renal Vein 25mm x 190mm (1" x 7 1/2")	49
46316ATT	Radiolucent Concave Lateral with Teeth 23mm x 45mm (1" x 1 3/4")	51	46570	Radiolucent Concave Renal Vein 25mm x 200mm (1" x 7 7/8")	49
46318AT	Radiolucent Concave Lateral 23mm x 50mm (1" x 2")	50	46571	Radiolucent Concave Renal Vein 25mm x 210mm (1" x 8 1/4")	49
46318ATT	Radiolucent Concave Lateral with Teeth 23mm x 50mm (1" x 2")	51	46572	Radiolucent Concave Renal Vein 25mm x 220mm (1" x 8 5/8")	49
46320AT	Radiolucent Concave Lateral 23mm x 55mm (1" x 2 1/8")	50	46573	Radiolucent Concave Renal Vein 25mm x 230mm (1" x 9")	49
46320ATT	Radiolucent Concave Lateral with Teeth 23mm x 55mm (1" x 2 1/8")	51	46574	Radiolucent Concave Renal Vein 25mm x 240mm (1" x 9 1/2")	49
46322AT	Radiolucent Concave Lateral 23mm x 60mm (1" x 2 3/8")	50	46575	Radiolucent Concave Renal Vein 25mm x 250mm (1" x 9 7/8")	49
46322ATT	Radiolucent Concave Lateral with Teeth 23mm x 60mm (1" x 2 3/8")	51	46660	Double Angled 25mm x 76mm (1" x 3")	41
46324AT	Radiolucent Concave Lateral 23mm x 65mm (1" x 2 1/2")	50	46661	Double Angled 25mm x 64mm (1" x 2 1/2")	41
46324ATT	Radiolucent Concave Lateral with Teeth 23mm x 65mm (1" x 2 1/2")	51	46662	Double Angled 13mm x 76mm (1/2" x 3")	41
46326AT	Radiolucent Concave Lateral 23mm x 70mm (1" x 2 3/4")	50	46663	Double Angled 13mm x 70mm (1/2" x 2 3/4")	41
46326ATT	Radiolucent Concave Lateral with Teeth 23mm x 70mm (1" x 2 3/4")	51	46664	Double Angled Bulb 25mm x 64mm (1" x 2 1/2")	41
46328AT	Radiolucent Concave Lateral 23mm x 75mm (1" x 3")	50	47200	CV - Clamp Assembly 2" L-Bar	64
46328ATT	Radiolucent Concave Lateral with Teeth 23mm x 75mm (1" x 3")	51	47200N	CV Clamp Mini L-Bar	64
46330AT	Radiolucent Concave Lateral 23mm x 80mm (1" x 3 1/8")	50	47200NF	CV Clamp Mini L-Bar Offset	64
46330ATT	Radiolucent Concave Lateral with Teeth 23mm x 80mm (1" x 3 1/8")	51	47200NR	CV Clamp Mini L-Bar Reverse	64
46332AT	Radiolucent Concave Lateral 23mm x 85mm (1" x 3 3/8")	50	47201	CV - Clamp Assembly 8" T-Bar	64
46332ATT	Radiolucent Concave Lateral with Teeth 23mm x 85mm (1" x 3 3/8")	51	47201N	CV Clamp T-Bar	64
46334AT	Radiolucent Concave Lateral 23mm x 90mm (1" x 3 1/2")	50	47201NF	CV Clamp 8" T-Bar Offset	64
46334ATT	Radiolucent Concave Lateral with Teeth 23mm x 90mm (1" x 3 1/2")	51	47202	CV - Clamp Assembly 6" x 8" L-Bar	64
46336AT	Radiolucent Concave Lateral 23mm x 95mm (1" x 3 3/4")	50	47202N	CV Clamp L-Bar	64
46336ATT	Radiolucent Concave Lateral with Teeth 23mm x 95mm (1" x 3 3/4")	51	47202NR	CV Clamp L-Bar - Reverse	64
46338AT	Radiolucent Concave Lateral 23mm x 100mm (1" x 4")	50	47204	Universal Joint	65
46338ATT	Radiolucent Concave Lateral with Teeth 23mm x 100mm (1" x 4")	51	47205T	1/4" Side Clip x 1/4" Slide with Thumb Screw	69
46340AT	Radiolucent Concave Lateral 23mm x 105mm (1" x 4 1/8")	50	47210	MIS Spine Frame U-Bar	68
46340ATT	Radiolucent Concave Lateral with Teeth 23mm x 105mm (1" x 4 1/8")	51	47211	MIS Spine Frame L-Bar	68
46342AT	Radiolucent Concave Lateral 23mm x 110mm (1" x 4 3/8")	50	47212	MIS Spine Frame U-Bar Add On	68
46342ATT	Radiolucent Concave Lateral with Teeth 23mm x 110mm (1" x 4 3/8")	51	47213	MIS Spine Frame L-Bar Add on	68
46344AT	Radiolucent Concave Lateral 23mm x 115mm (1" x 4 1/2")	50	48001	Towel Clip Hook 45mm x 64mm (1 3/4" x 2 1/2")	73
46344ATT	Radiolucent Concave Lateral with Teeth 23mm x 115mm (1" x 4 1/2")	51	50000A3	Instrument Case - Accessory 3 Compartments (23" x 10 3/4" x 2 1/2")	22
46346AT	Radiolucent Concave Lateral 23mm x 120mm (1" x 4 3/4")	50	50000ACT	Instrument Case - Anterior Cervical (22" x 11" x 3 1/2")	23
46346ATT	Radiolucent Concave Lateral with Teeth 23mm x 120mm (1" x 4 3/4")	51	50000ALR	Instrument Case - AL Renal Vein Blades (23" x 10 3/4" x 2 1/2")	23
46400	Tummy Tuck 76mm x 140mm (3" x 5 1/2")	73	50000ALT	Instrument Case - Anterior Lumbar (22" x 11" x 3 1/2")	23
46401	Tummy Tuck 102mm x 140mm (4" x 5 1/2")	73	50000BL	Instrument Case (26" x 10" x 5")	22

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
50000BOL	Instrument Case - Bolling Mitral Valve (23" x 10 3/4" x 2 1/2")	22	S45622Q	Swivel Army Navy 13mm x 75mm (1/2" x 3 1/16")	56
50000CCAB	Instrument Case - CCA Blades (23" x 10 3/4" x 2 1/4")	23	S45640Q	Swivel Concave 18mm x 70mm (3/8" x 2 3/4")	50
50000CCAF	Instrument Case - CCA Frame (23" x 10 3/4" x 2 3/4")	23	S45645Q	Swivel Cloward 19mm x 44mm (3/4" x 1 3/4")	52
50000DH	Instrument Case - Drake Hook (12" x 4 1/2" x 2 1/4")	22	SL42101	Dual Angle Handle	61
50000EAC	Instrument Case - Expanded Cervical (23" x 10 3/4" x 2")	23	SL42102	Free Angle Handle	61
50000ERL	Instrument Case - Expanded Reverse Lip (23 1/2" x 11" x 5 1/4")	23	SL42125WA	Cam II Slide-on Angling 8"	17
50000G	Instrument Case with Pin Mat (23" x 11" x 3 1/2")	22	SL42125WH	Cam II Slide-On 8" (non-angling)	17
50000HPB	Instrument Case - Hepatobiliary (22" x 11" x 3 1/2")	22	SL42126LH	Clip-On with Hinge 8"	72
50000LW	Instrument Case (6 1/4" x 14" x 2 1/2")	22	SL42126WA	Cam II Angling Handle 8"	17
50000ML	Instrument Case - PLA Midline Blades (22" x 10" x 2 1/2")	23	SL42126WAL	Cam II Angling Handle 10"	17
50000MSF	Instrument Case - MIS Spine Frame (23 1/4" x 11" x 4 1/2")	23	SL42126WGLP	Cam II Clip-on Quick Angle 10"	17
50000PLA	Instrument Case - PLA (23 1/2" x 11" x 4")	23	SL42126WGP	Cam II Clip-on Quick Angle 8"	17
50000PLAO	Instrument Case - PLA and One Step Dilator (23 1/2" x 11" x 4")	23	SL42126WH	Cam II Clip-on 8" (non-angling)	17
50000RLT	Instrument Case - Reverse Lip (22" x 11" x 3 1/2")	23	SL42126WHP	Cam II Clip-on 8" (non-angling)	17
50000UAC	Instrument Case - Cervical Blade (22" x 10" x 2 1/2")	23	SL42127WA	Cam II Clip-on Angling 8" (3/8")	17
50000UACT	Instrument Case - Cervical Blade Tray (Insert for 50000UAC)	23	SL42127WGP	Cam II Clip-on Quick Angle 8" (3/8")	17
51229	Stylet for SALS Tip	21	SL42128F	Angling Hand Held Adapter	19
51234	Suction for AL Surgery with SALS Tip Stylet	21	SL42128G	Quick Angle Hand Held Adapter	19
51236	Depth Gauge Anterior Lumbar 335mm	21	SL42139E	Laparoscopic Rigid Holder for Elite	62
51251	Drake Hook 258mm (10")	65	SL42139N	Laparoscopic Rigid Holder for Power Rail Clamp	62
60008	MIS Spine Frame Right Arm	68	SL42161	S-Lock Nathanson Adapter	63
60009	MIS Spine Frame Left Arm	68	SL42190	S-Lock Articulating Arm	66
60010	MIS Spine Frame Rack 200mm (7 7/8")	68	SL42200	Nathanson Hook Sm (5mm) 169mm x 125mm x 12mm (6 5/8" x 5" x 1/2")	63
60011	MIS Spine Frame Rack 150mm (5 7/8")	68	SL42201	Nathanson Hook Med (5mm) 203mm x 167mm x 16mm (8" x 6 1/2" x 5/8")	63
60012	MIS Spine Frame Rack 250mm (9 7/8")	68	SL42202	Nathanson Hook Lg (5mm) 240mm x 222mm x 21mm (9 1/2" x 8 3/4" x 7/8")	63
60015	MIS Spine Frame Dual Blade Adapter	71	SL42203	Nathanson Hook Sm (6.5mm) 169mm x 125mm x 12mm (6 5/8" x 5" x 1/2")	63
60020	T-Handle	19	SL42204	Nathanson Hook Med (6.5mm) 203mm x 167mm x 16mm (8" x 6 1/2" x 5/8")	63
60021	Hex Tool 3/16" with T-Handle	61	SL42205	Nathanson Hook Lg (6.5mm) 240mm x 222mm x 21mm (9 1/2" x 8 3/4" x 7/8")	63
60025	Orthopedic Retractor Handle Knob	67	SL45002AS	Angling Handle 6" 3 dimple	69
64047	Silicone Sleeve Small	72	SL45003	Interchangeable 8" (non-angling)	17
64049	Silicone Sleeve Large	72	SL45003AS	Angling Handle 8"	17
300251	Morse Sternal Retractor (Blades 28mm x 22mm / 1 1/8" x 7/8")	64	SL45003ASL	Angling Handle 10"	17
300253	Morse Sternal Spreader Large (Blades 28mm x 38mm / 1 1/8" x 1 1/2")	64	SL45003G	Quick Angle Handle 6"	69
LLS-2000-2	LED Light Source	20	SL45003GS	Quick Angle Handle 8"	17
P45018	Rotatable Malleable 25mm x 102mm (1" x 4")	30	SL45003GSL	Quick Angle Handle 10"	17
P45022	Rotatable Malleable 51mm x 102mm (2" x 4")	30	SL45003L	Handle 8" (non-angling)	17
P45062	Rotatable Harrington 38mm x 127mm (1 1/2" x 5")	27	SL45003LH	Hinged Handle 8"	72
P45070	Rotatable Deaver 25mm x 102mm (1" x 4")	36	SL45003LM	Handle 10" (non-angling)	17
P45074	Rotatable Deaver 51mm x 102mm (2" x 4")	36	SL45004AL	Angling Thompson Adapter Handle for Omni-Tract® 8 3/4"	18
S43124	Fixed Malleable 64mm x 254 (2 1/2" x 10") - 6"	30	SL45004GL	Quick Angle Thompson Adapter Handle for Omni-Tract® 8 3/4"	18
S43125	Fixed Richardson 19mm x 25mm (3/4" x 1") - 6"	35	SL45004L	Thompson Adapter Handle for Omni-Tract® 8 3/4" (non-angling)	18
S43126	Fixed Richardson 25mm x 32mm (1" x 1 1/4") - 6"	35	SL45005AL	Ratchet Angling Retractor Handle 10 1/4"	18
S43127	Fixed Kelly 38mm x 51mm (1 1/2" x 2") - 6"	28	SL45005GL	Ratchet Quick Angle Retractor Handle 10 1/4"	18
S43128	Fixed Kelly 51mm x 64mm (2" x 2 1/2") - 6"	28	SL45005L	Ratchet Retractor Handle 10 1/4" (non-angling)	18
S43129	Fixed Kelly 64mm x 76mm (2 1/2" x 3") - 6"	28	SL45006CA	Micro-Adjustable II Clip-On Angling 10"	16
S43130	Fixed Kelly 76mm x 89mm (3" x 3 1/2") - 6"	28	SL45006CAR	Micro-Adjustable II Clip-on Angling 15"	16
S43140	Fixed Balfour [with lips] 70mm x 73mm (2 3/4" x 3") - 6"	25	SL45006CGP	Micro-Adjustable II Clip-on Quick Angle 10"	16
S43141	Fixed Balfour [with lips] 83mm x 73mm (3 1/4" x 3") - 6"	25	SL45006CGRP	Micro-Adjustable II Clip-on Quick Angle 15"	16
S43151	Fixed Richardson 51mm X 127mm (2" x 5") - 6"	34	SL45006STA	Micro-Adjustable II Slide-on Angling 10"	16
S43170	Fixed Deaver 25mm x 127mm (1" x 5") - 6"	36	SL45006STC	Micro-Adjustable II Clip-on 10 3/4" (non-angling)	16
S43172	Fixed Deaver 38mm x 127mm (1 1/2" x 5") - 6"	36	SL45006STCP	Micro-Adjustable II Clip-on 10 3/4" (non-angling)	16
S43176	Fixed Deaver 76mm x 146mm (3" x 5 3/4") - 6"	36	SL45007CA	Micro-Adjustable II Clip-on Angling 7 3/4"	16
S43182	Fixed Handle Heaney 25mm x 127mm (1" x 5") - 6"	42	SL45007CGP	Micro-Adjustable II Clip-on Quick Angle 7 3/4"	16
S43190	Fixed Double Bladed 76mm x 114mm (3" x 4 1/2")	42	SL45007P	Micro-Adjustable II Clip-On Angling 20"	16
S43191	Fixed Malleable Finger 152mm (6") - 6"	31	SL45007PGP	Micro-Adjustable II Clip-on Quick Angle 20"	16
S43192	Fixed Bladder 35mm x 89mm (1 3/8" x 3 1/2")	43	SL45007STA	Micro-Adjustable II Slide-on Angling 7 3/4"	16
S45125	Swivel Richardson 19mm x 25mm (3/4" x 1") - 6"	35	SL45007STC	Micro-Adjustable II Clip-On 8 3/4" (non-angling)	16
S45125A	Swivel Richardson 19mm x 51mm (3/4" x 2") - 6"	35	SL45007STCP	Micro-Adjustable II Clip-On 8 3/4" (non-angling)	16
S45126	Swivel Richardson 25mm x 32mm (1" x 1 1/4") - 6"	35	SL45007STG	Micro-Adjustable II Slide-on Quick Angle 7 3/4"	16
S45127	Swivel Kelly 38mm x 51mm (1 1/2" x 2") - 6"	28	SL45010A	Mini Micro-Adjustable Clip-On Angling 8" (cervical use only)	16
S45142	Swivel Balfour [with lips] 32mm x 32mm (1 1/4" x 1 1/4") - 6"	25	SL45012S	Interchangeable Head with 1/4" Scissor	69
S45151	Swivel Richardson 51mm x 127mm (2" x 5") 6"	34	SL45014AL	Angling Thompson Adapter Handle for SynFrame® 8 3/4"	18
S45601Q	Swivel Meyerding 16mm x 100mm (5/8" x 4")	57	SL45014GL	Quick Angle Thompson Adapter Handle for SynFrame® 8 3/4"	18
S45611Q	Swivel Meyerding 25mm x 76mm (1" x 3")	57	SL45014L	Thompson Adapter Handle for SynFrame® 8 3/4" (non-angling)	18
S45616Q	Swivel Meyerding 51mm x 125mm (2" x 5")	57	SL45150A	Concave Lateral 19mm x 55mm (3/4" x 2 1/4")	50
S45620Q	Swivel Army Navy 13mm x 55mm (1/2" x 2 1/8")	56	SL45150AT	Radiolucent Concave Lateral 22mm x 60mm (7/8" x 2 3/8")	51
S45621Q	Swivel Army Navy 13mm x 65mm (1/2" x 2 5/8")	56	SL45150ATT	Radiolucent Concave Lateral with Teeth 22mm x 60mm (7/8" x 2 3/8")	51

PRODUCT INDEX

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
SL45151A	Concave Lateral 25mm x 50mm (1" x 2")	50	SL45632AT	Radiolucent MIS Tapered 15mm x 23mm x 75mm (5/8" x 1" x 3")	55
SL45151AT	Radiolucent Concave Lateral 22mm x 70mm (7/8" x 2 3/4")	51	SL45640A	Concave Lateral 18mm x 70mm (5/8" x 2 3/4")	50
SL45151ATT	Radiolucent Concave Lateral with Teeth 22mm x 70mm (7/8" x 2 3/4")	51	SL45641A	Concave Tapered 45mm x 50mm (1 3/4" x 2")	52
SL45152A	Concave Lateral 25mm x 85mm (1" x 3 3/8")	50	SL45641AT	Radiolucent Concave Tapered 51mm x 50mm (2" x 2")	52
SL45152AT	Radiolucent Concave Lateral 22mm x 80mm (7/8" x 3 1/8")	51	SL45641ATT	Radiolucent Concave Tapered with Teeth 51mm x 50mm (2" x 2")	53
SL45152ATT	Radiolucent Concave Lateral with Teeth 22mm x 80mm (7/8" x 3 1/8")	51	SL45642A	Concave Tapered 45mm x 65mm (1 3/4" x 2 1/2")	52
SL45153A	Concave Medial 32mm x 70mm (1 1/4" x 2 3/4")	56	SL45642AT	Radiolucent Concave Tapered 51mm x 65mm (2" x 2 1/2")	52
SL45154A	Concave Medial 38mm x 70mm (1 1/2" x 2 3/4")	56	SL45642ATT	Radiolucent Concave Tapered with Teeth 51mm x 65mm (2" x 2 1/2")	53
SL45155A	Concave Medial 38mm x 85mm (1 1/2" x 3 3/8")	56	SL45643A	Concave Tapered 19mm x 35mm (3/4" x 1 3/8")	52
SL45156A	Concave Lateral 25mm x 30mm (1" x 1 1/4")	50	SL45643AT	Radiolucent Concave Tapered 23mm x 35mm (1" x 1 3/8")	52
SL45157A	Concave Lateral 25mm x 40mm (1" x 1 1/2")	50	SL45643ATT	Radiolucent Concave Tapered with Teeth 23mm x 35mm (1" x 1 3/8")	52
SL45158A	Concave Lateral 25mm x 60mm (1" x 2 3/8")	50	SL45644AT	Radiolucent Concave Tapered 23mm x 40mm (1" x 1 1/2")	52
SL45159A	Concave Lateral 25mm x 70mm (1" x 2 3/4")	50	SL45644ATT	Radiolucent Concave Tapered with Teeth 23mm x 40mm (1" x 1 1/2")	52
SL45160A	Concave Longitudinal 23mm x 30mm (1" x 1 1/4")	53	SL45645A	Concave Tapered 19mm x 45mm (3/4" x 1 3/8")	52
SL45160AT	Radiolucent Concave Longitudinal 23mm x 30mm (1" x 1 1/4")	53	SL45645AT	Radiolucent Concave Tapered 23mm x 45mm (3/4" x 1 3/8")	52
SL45161A	Concave Longitudinal 23mm x 50mm (1" x 2")	53	SL45645ATT	Radiolucent Concave Tapered with Teeth 23mm x 45mm (3/4" x 1 3/8")	52
SL45161AT	Radiolucent Concave Longitudinal 23mm x 40mm (1" x 1 1/2")	53	SL45646A	Concave Tapered 19mm x 50mm (3/4" x 2")	52
SL45162A	Concave Longitudinal 23mm x 70mm (1" x 2 3/4")	53	SL45646AT	Radiolucent Concave Tapered 23mm x 50mm (1" x 2")	52
SL45162AT	Radiolucent Concave Longitudinal 23mm x 50mm (1" x 2")	53	SL45646ATT	Radiolucent Concave Tapered with Teeth 23mm x 50mm (1" x 2")	52
SL45163AT	Radiolucent Concave Longitudinal 23mm x 60mm (1" x 2 3/8")	53	SL45647AT	Radiolucent Concave Tapered 23mm x 55mm (1" x 2 1/8")	52
SL45164AT	Radiolucent Concave Longitudinal 23mm x 70mm (1" x 2 3/4")	53	SL45647ATT	Radiolucent Concave Tapered with Teeth 23mm x 55mm (1" x 2 1/8")	52
SL45165AT	Radiolucent Concave Longitudinal 23mm x 80mm (1" x 3 1/8")	53	SL45648AT	Radiolucent Concave Tapered 23mm x 60mm (1" x 2 3/8")	52
SL45166AT	Radiolucent Concave Longitudinal 23mm x 90mm (1" x 3 1/2")	53	SL45648ATT	Radiolucent Concave Tapered with Teeth 23mm x 60mm (1" x 2 3/8")	52
SL45167AT	Radiolucent Concave Longitudinal 23mm x 100mm (1" x 4")	53	SL45649AT	Radiolucent Concave Tapered 23mm x 65mm (1" x 2 1/2")	52
SL45168AT	Radiolucent Concave Longitudinal 23mm x 110mm (1" x 4 3/8")	53	SL45649ATT	Radiolucent Concave Tapered with Teeth 23mm x 65mm (1" x 2 1/2")	52
SL45212A	Concave Medial 13mm x 50mm (1/2" x 2")	56	SL45650AT	Radiolucent Concave Tapered 23mm x 70mm (1" x 2 3/4")	52
SL45214A	Concave Medial 13mm x 75mm (1/2" x 3")	56	SL45650ATT	Radiolucent Concave Tapered with Teeth 23mm x 70mm (1" x 2 3/4")	52
SL45216A	Concave Medial 13mm x 100mm (1/2" x 4")	56	SL45651AT	Radiolucent Concave Tapered 23mm x 75mm (1" x 3")	52
SL45222A	Concave Medial 19mm x 50mm (3/4" x 2")	56	SL45651ATT	Radiolucent Concave Tapered with Teeth 23mm x 75mm (1" x 3")	53
SL45224A	Concave Medial 19mm x 75mm (3/4" x 3")	56	SL45652AT	Radiolucent Concave Tapered 23mm x 80mm (1" x 3 1/8")	52
SL45226A	Concave Medial 19mm x 100mm (3/4" x 4")	56	SL45652ATT	Radiolucent Concave Tapered with Teeth 23mm x 80mm (1" x 3 1/8")	53
SL45232A	Concave Medial 25mm x 50mm (1" x 2")	56	SL45656A	Malleable 25mm x 110mm (1" x 4 1/2")	30
SL45234A	Concave Medial 25mm x 75mm (1" x 3")	56	SL45656AT	Radiolucent Malleable 25mm x 110mm (1" x 4 1/2")	30
SL45236A	Concave Medial 25mm x 100mm (1" x 4")	56	SL45657A	Malleable 13mm x 90mm (1/2" x 3 1/2")	30
SL45600A	Meyerding 16mm x 75mm (5/8" x 3")	57	SL45657AT	Radiolucent Malleable 13mm x 90mm (1/2" x 3 1/2")	30
SL45600AT	Meyerding with Teeth 16mm x 75mm (5/8" x 3")	57	SL45666AT	Radiolucent Odontoid 20mm x 90mm (3/4" x 3 1/2")	54
SL45601A	Meyerding 16mm x 100mm (5/8" x 4")	57	SL45750A	Midline R 70mm x 83mm (2 3/4" x 3 1/4")	55
SL45601AT	Meyerding with Teeth 16mm x 100mm (5/8" x 4")	57	SL45750AT	Radiolucent Midline R 70mm x 83mm (2 3/4" x 3 1/4")	55
SL45602A	Meyerding 25mm x 100mm (1" x 4")	57	SL45751A	Midline L 70mm x 83mm (2 3/4" x 3 1/4")	55
SL45602AT	Meyerding with Teeth 25mm x 100mm (1" x 4")	57	SL45751AT	Radiolucent Midline L 70mm x 83mm (2 3/4" x 3 1/4")	55
SL45603A	Meyerding 25mm x 125mm (1" x 5")	57	SL45753AM	Prostate Retractor 8mm (3/8")	44
SL45603AT	Meyerding with Teeth 25mm x 125mm (1" x 5")	57	SL45754A	Midline R 57mm x 83mm (2 1/4" x 3 1/4")	55
SL45604A	Peress 60mm x 102mm (2 3/8" x 4")	59	SL45755A	Midline L 57mm x 83mm (2 1/4" x 3 1/4")	55
SL45605A	Peress 60mm x 127mm (2 3/8" x 5")	59	SL45760A	Prong 75mm (3")	58
SL45606A	Peress 60mm x 152mm (2 3/8" x 6")	59	SL45761A	Prong 100mm (4")	58
SL45607A	Peress 35mm x 102mm (1 3/8" x 4")	59	SL45764	PLA Hook 40mm (1 1/2")	71
SL45608A	Peress 35mm x 127mm (1 3/8" x 5")	59	SL45765	PLA Hook 50mm (2")	71
SL45609A	Peress 35mm x 152mm (1 3/8" x 6")	59	SL45766	PLA Hook 60mm (2 3/8")	71
SL45610A	Meyerding 25mm x 60mm (1" x 2 3/8")	57	SL45767	PLA Hook 70mm (2 3/4")	71
SL45610AT	Meyerding with Teeth 25mm x 60mm (1" x 2 3/8")	57	SL45768	PLA Hook 80mm (3 1/8")	71
SL45611A	Meyerding 25mm x 75mm (1" x 3")	57	SL45769	PLA Hook 90mm (3 1/2")	71
SL45611AT	Meyerding with Teeth 25mm x 75mm (1" x 3")	57	SL45770	PLA Hook 100mm (4")	71
SL45614A	Meyerding 16mm x 51mm (5/8" x 2")	57	SL45771	PLA Hook 110mm (4 3/8")	71
SL45614AT	Meyerding with Teeth 16mm x 51mm (5/8" x 2")	57	SL45772	PLA Hook 120mm (4 3/4")	71
SL45618AT	Meyerding with Teeth 25mm x 115mm (1" x 4 1/2")	57	SL46115	Malleable 13mm x 64mm (1/2" x 2 1/2")	30
SL45619AT	Meyerding with Teeth 38mm x 165mm (1 1/2" x 6 1/2")	57	SL46115ET	Radiolucent Malleable 13mm x 64mm (1/2" x 2 1/2")	30
SL45620A	Army Navy 13mm x 55mm (1/2" x 2 1/8")	56	SL46116	Malleable 13mm x 152mm (1/2" x 6")	30
SL45621A	Army Navy 13mm x 65mm (1/2" x 2 3/8")	56	SL46116ET	Radiolucent Malleable 13mm x 152mm (1/2" x 6")	30
SL45622A	Army Navy 13mm x 75mm (1/2" x 3 1/8")	56	SL46117	Malleable 13mm x 203mm (1/2" x 8")	30
SL45624AT	Radiolucent MIS Tapered 15mm x 23mm x 35mm (5/8" x 1" x 1 3/8")	55	SL46117ET	Radiolucent Malleable 13mm x 203mm (1/2" x 8")	30
SL45625AT	Radiolucent MIS Tapered 15mm x 23mm x 40mm (5/8" x 1" x 1 1/2")	55	SL46118	Malleable 25mm x 203mm (1" x 8")	30
SL45626AT	Radiolucent MIS Tapered 15mm x 23mm x 45mm (5/8" x 1" x 1 3/4")	55	SL46118A	Malleable 25mm x 254mm (1" x 10")	30
SL45627AT	Radiolucent MIS Tapered 15mm x 23mm x 50mm (5/8" x 1" x 2")	55	SL46118B	Malleable 38mm x 254mm (1 1/2" x 10")	30
SL45628AT	Radiolucent MIS Tapered 15mm x 23mm x 55mm (5/8" x 1" x 2 1/8")	55	SL46118ET	Radiolucent Malleable 25mm x 203mm (1" x 8")	30
SL45629AT	Radiolucent MIS Tapered 15mm x 23mm x 60mm (5/8" x 1" x 2 3/8")	55	SL46118P	Malleable 25mm x 102mm (1" x 4")	30
SL45630AT	Radiolucent MIS Tapered 15mm x 23mm x 65mm (5/8" x 1" x 2 1/2")	55	SL46118PET	Radiolucent Malleable 25mm x 102mm (1" x 4")	30
SL45631AT	Radiolucent MIS Tapered 15mm x 23mm x 70mm (5/8" x 1" x 2 3/4")	55	SL46118R	Malleable Renal Vein 25mm x 203mm (1" x 8")	38

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
SL46118RET	Radiolucent Renal Vein 25mm x 203mm (1" x 8")	37	SL46172	Deaver 38mm x 127mm (1 1/2" x 5")	36
SL46119	Malleable 25mm x 127mm (1" x 5")	30	SL46172ET	Radiolucent Deaver 38mm x 127mm (1 1/2" x 5")	36
SL46119B	Concave Renal Vein 25mm x 140mm (1" x 5 1/2")	37	SL46174	Deaver 51mm x 127mm (2" x 5")	36
SL46119BET	Radiolucent Malleable Renal Vein 25mm x 140mm (1" x 5 1/2")	37	SL46174ET	Radiolucent Deaver 51mm x 127mm (2" x 5")	36
SL46119C	Concave Renal Vein 25mm x 191mm (1" x 7 1/2")	37	SL46176ET	Radiolucent Deaver 76mm x 127mm (3" x 5")	36
SL46119CET	Radiolucent Malleable Renal Vein 25mm x 190mm (1" x 7 1/2")	37	SL46178ET	Radiolucent Deaver 100mm x 127mm (4" x 5")	36
SL46119RET	Radiolucent Rigid Renal Vein 25mm x 140mm (1" x 5 1/2")	37	SL46180	Heany 25mm x 102mm (1" x 4")	42
SL46120RET	Radiolucent Rigid Renal Vein 25mm x 190mm (1" x 7 1/2")	37	SL46182	Heany 25mm x 127mm (1" x 5")	42
SL46121RET	Radiolucent Rigid Renal Vein 25mm x 255mm (1" x 10")	37	SL46183	Allison Lung Retractor 229mm (9")	40
SL46122	Malleable 51mm x 203mm (2" x 8")	30	SL46185	Hill-Ferguson 32mm x 76mm (1 1/4" x 3")	45
SL46122A	Malleable 51mm x 254mm (2" x 10")	30	SL46186	Hill-Ferguson 25mm x 64mm (1" x 2 1/2")	45
SL46122AET	Radiolucent Malleable 51mm x 254mm (2" x 10")	30	SL46187	Rectal 13mm x 117mm (1/2" x 4 1/2")	45
SL46122ET	Radiolucent Malleable 51mm x 203mm (2" x 8")	31	SL46190	Double Bladed 76mm x 114mm (3" x 4 1/2")	42
SL46122P	Malleable 51mm x 102mm (2" x 4")	30	SL46191	Malleable Finger 152mm (6")	31
SL46122RET	Radiolucent Renal Rigid Vein 25mm x 305mm (1" x 12")	37	SL46191ET	Radiolucent Malleable Finger 152mm (6")	31
SL46122V	Malleable Slotted 51mm x 203mm (2" x 8")	32	SL46191L	Malleable Finger 203mm (8")	31
SL46123	Malleable 38mm x 203mm (1 1/2" x 8")	30	SL46191P	Malleable Finger 102mm (4")	31
SL46124	Malleable 64mm x 254mm (2 1/2" x 10")	30	SL46191SP	Malleable Finger 76mm (3")	31
SL46124ET	Radiolucent Malleable 64mm x 254mm (2 1/2" x 10")	31	SL46192ET	Reverse Lip Malleable 25mm x 200mm (1" x 7 7/8")	47
SL46125	Richardson 19mm x 25mm (3/4" x 1")	35	SL46193ET	Reverse Lip Malleable 32mm x 200mm (1 1/4" x 7 7/8")	47
SL46125A	Richardson 19mm x 51mm (3/4" x 2")	35	SL46194ET	Reverse Lip Rigid 25mm x 100mm (1" x 4")	48
SL46126	Richardson 25mm x 32mm (1" x 1 1/4")	35	SL46194RET	Reverse Lip Rigid 25mm x 100mm (1" x 4")	48
SL46127	Kelly 38mm x 51mm (1 1/2" x 2")	28	SL46195ET	Reverse Lip Rigid 25mm x 150mm (1" x 5 7/8")	48
SL46127ET	Radiolucent Kelly 38mm x 51mm (1 1/2" x 2")	28	SL46195RET	Reverse Lip Rigid 25mm x 150mm (1" x 5 7/8")	48
SL46128	Kelly 51mm x 64mm (2" x 2 1/2")	28	SL46195TET	Reverse Lip Tapered 25mm x 150mm (1" x 5 7/8")	48
SL46128ET	Radiolucent Kelly 51mm x 64mm (2" x 2 1/2")	28	SL46196CET	Reverse Lip Concave 25mm x 200mm (1" x 7 7/8")	48
SL46137	Holtgrewe 60mm x 152mm (2 3/8" x 6")	43	SL46196ET	Reverse Lip Rigid 25mm x 200mm (1" x 7 7/8")	48
SL46139	"Y" Bladder 81mm x 216mm (3 1/8" x 8 1/2")	43	SL46196RET	Reverse Lip Rigid 25mm x 200mm (1" x 7 7/8")	48
SL46140ET	Radiolucent Balfour [with lips] 65mm x 72mm (2 1/2" x 2 3/4")	25	SL46196TET	Reverse Lip Tapered 25mm x 200mm (1" x 7 7/8")	48
SL46140P	Balfour [with lips] 76mm x 51mm (3" x 2")	25	SL46197ET	Reverse Lip Rigid 32mm x 100mm (1 1/4" x 4")	48
SL46141ET	Radiolucent Balfour [with lips] 83mm x 70mm (3 1/4" x 2 3/4")	25	SL46197RET	Reverse Lip Rigid 32mm x 100mm (1 1/4" x 4")	48
SL46142	Balfour [with lips] 32mm x 32mm (1 1/4" x 1 1/4")	25	SL46198ET	Reverse Lip Rigid 32mm x 150mm (1 1/4" x 5 7/8")	48
SL46142EB	Radiolucent Balfour [with lips] 32mm x 32mm (1 1/4" x 1 1/4") [Anodized]	25	SL46198RET	Reverse Lip Rigid 32mm x 150mm (1 1/4" x 5 7/8")	48
SL46142ET	Radiolucent Balfour [with lips] 32mm x 32mm (1 1/4" x 1 1/4")	25	SL46199CET	Reverse Lip Concave 32mm x 200mm (1 1/4" x 7 7/8")	48
SL46142P	Balfour [with lips] 51mm x 54mm (2" x 2")	25	SL46199ET	Reverse Lip Rigid 32mm x 200mm (1 1/4" x 7 7/8")	48
SL46142PEB	Radiolucent Balfour [with lips] 51mm x 51mm (2" x 2") [Anodized]	25	SL46199RET	Reverse Lip Rigid 32mm x 200mm (1 1/4" x 7 7/8")	48
SL46143	Fenestrated 44mm x 89mm (1 3/4" x 3 1/2")	38	SL46200CET	Reverse Lip Concave 25mm x 255mm (1" x 10")	48
SL46144	Fenestrated 48mm x 127mm (1 7/8" x 5")	38	SL46200ET	Reverse Lip Rigid 25mm x 255mm (1" x 10")	48
SL46145	Splanchchnich 51mm x 114mm (2" x 4 1/2")	29	SL46200RET	Reverse Lip Rigid 25mm x 255mm (1" x 10")	48
SL46145ET	Radiolucent Splanchchnich 51mm x 115mm (2" x 4 1/2")	29	SL46201CET	Reverse Lip Concave 32mm x 255mm (1 1/4" x 10")	48
SL46146	Splanchchnich 51mm x 152mm (2" x 6")	29	SL46201ET	Reverse Lip Rigid 32mm x 255mm (1 1/4" x 10")	48
SL46146ET	Radiolucent Splanchchnich 51mm x 152mm (2" x 6")	29	SL46201RET	Reverse Lip Rigid 32mm x 255mm (1 1/4" x 10")	48
SL46147ET	Radiolucent Splanchchnich 51mm x 203mm (2" x 8")	29	SL46210	Volkman Sharp 6 Prong - 25mm x 13mm (1" x 1/2")	39
SL46148ET	Radiolucent Splanchchnich 25mm x 152mm (1" x 6")	29	SL46211	Volkman Blunt 6 Prong - 25mm x 13mm (1" x 1/2")	39
SL46149	Richardson 44mm x 51mm (1 3/4" x 2")	34	SL46212	Volkman Sharp 4 Prong - 16mm x 13mm (5/8" x 1/2")	39
SL46149ET	Radiolucent Splanchchnich 25mm x 203mm (1" x 8")	29	SL46213	Volkman Blunt 4 Prong - 16mm x 13mm (5/8" x 1/2")	39
SL46150	Richardson 51mm x 102mm (2" x 4")	34	SL46214	Murphy Sharp 4 Prong - 25mm x 13mm (1" x 1/2")	39
SL46150ET	Radiolucent Richardson 51mm x 102mm (2" x 4")	34	SL46215	Murphy Blunt 4 Prong - 25mm x 13mm (1" x 1/2")	39
SL46154	Richardson 51mm x 76mm (2" x 3")	34	SL46216	Murphy Sharp 3 Prong - 19mm x 13mm (3/4" x 1/2")	39
SL46154ET	Radiolucent Richardson 51mm x 76mm (2" x 3")	34	SL46217	Murphy Blunt 3 Prong - 19mm x 13mm (3/4" x 1/2")	39
SL46155	Mini-Open Richardson 51mm x 76mm (2" x 3")	33	SL46220	Mini-Open Malleable 25mm x 75mm (1" x 3")	32
SL46156	Mini-Open Richardson 51mm x 102mm (2" x 4")	33	SL46221	Mini-Open Malleable 25mm x 100mm (1" x 4")	32
SL46157	Mini-Open Richardson 51mm x 127mm (2" x 5")	33	SL46222	Mini-Open Malleable 25mm x 125mm (1" x 5")	32
SL46158	Mini-Open Richardson 51mm x 152mm (2" x 6")	33	SL46224	Mini-Open Double Angled 51mm x 127mm (2" x 5")	33
SL46160	Harrington 64mm x 152mm (2 1/2" x 6")	27	SL46225	Mini-Open Double Angled 51mm x 152mm (2" x 6")	33
SL46160E	Radiolucent Harrington 64mm x 152mm (2 1/2" x 6") [Aluminum]	27	SL46226	Mini-Open Double Angled 51mm x 165mm (2" x 6 1/2")	33
SL46160ET	Radiolucent Harrington 64mm x 152mm (2 1/2" x 6")	27	SL46227	Mini-Open Deep Double Angled 51mm x 127mm (2" x 5")	33
SL46161	Harrington 114mm x 254mm (4 1/2" x 10")	27	SL46260RET	Reverse Lip Rigid 25mm x 100mm (1" x 4")	46
SL46162P	Harrington 38mm x 127mm (1 1/2" x 5")	27	SL46260TET	Reverse Lip Tapered 25mm x 100mm (1" x 4")	47
SL46163	Harrington 64mm x 203mm (2 1/2" x 8")	27	SL46261RET	Reverse Lip Rigid 25mm x 110mm (1" x 4 3/8")	46
SL46163ET	Radiolucent Harrington 64mm x 203mm (2 1/2" x 8")	27	SL46262RET	Reverse Lip Rigid 25mm x 120mm (1" x 4 3/8")	46
SL46164	Harrington 64mm x 254mm (2 1/2" x 10")	27	SL46263RET	Reverse Lip Rigid 25mm x 130mm (1" x 5 1/8")	46
SL46164ET	Radiolucent Harrington 64mm x 254mm (2 1/2" x 10")	27	SL46264RET	Reverse Lip Rigid 25mm x 140mm (1" x 5 1/2")	46
SL46165	Harrington 64mm x 102mm (2 1/2" x 4")	27	SL46265RET	Reverse Lip Rigid 25mm x 150mm (1" x 5 7/8")	46
SL46165ET	Radiolucent Harrington 64mm x 102mm (2 1/2" x 4")	27	SL46265TET	Reverse Lip Tapered 25mm x 150mm (1" x 5 7/8")	47
SL46170	Deaver 25mm x 127mm (1" x 5")	36	SL46266RET	Reverse Lip Rigid 25mm x 160mm (1" x 6 1/4")	46
SL46170ET	Radiolucent Deaver 25mm x 127mm (1" x 5")	36	SL46267RET	Reverse Lip Rigid 25mm x 170mm (1" x 6 5/8")	46


PRODUCT INDEX

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
SL46268RET	Reverse Lip Rigid 25mm x 180mm (1" x 7")	46	SL46510ATT	Radiolucent Concave Micro with Teeth 14mm x 35mm (½" x 1 ¾")	54
SL46269RET	Reverse Lip Rigid 25mm x 190mm (1" x 7 ½")	46	SL46512AT	Radiolucent Concave Micro 14mm x 40mm (½" x 1 ½")	54
SL46270RET	Reverse Lip Rigid 25mm x 200mm (1" x 7 ¾")	46	SL46512ATT	Radiolucent Concave Micro with Teeth 14mm x 40mm (½" x 1 ½")	54
SL46270TET	Reverse Lip Tapered 25mm x 200mm (1" x 7 ¾")	47	SL46514AT	Radiolucent Concave Micro 14mm x 45mm (½" x 1 ½")	54
SL46271RET	Reverse Lip Rigid 25mm x 210mm (1" x 8 ¼")	46	SL46514ATT	Radiolucent Concave Micro with Teeth 14mm x 45mm (½" x 1 ½")	54
SL46272RET	Reverse Lip Rigid 25mm x 220mm (1" x 8 ¾")	46	SL46516AT	Radiolucent Concave Micro 14mm x 50mm (½" x 2")	54
SL46273RET	Reverse Lip Rigid 25mm x 230mm (1" x 9")	46	SL46516ATT	Radiolucent Concave Micro with Teeth 14mm x 50mm (½" x 2")	54
SL46274RET	Reverse Lip Rigid 25mm x 240mm (1" x 9 ½")	46	SL46518AT	Radiolucent Concave Micro 14mm x 55mm (½" x 2 ½")	54
SL46275RET	Reverse Lip Rigid 25mm x 250mm (1" x 9 ¾")	46	SL46518ATT	Radiolucent Concave Micro with Teeth 14mm x 55mm (½" x 2 ½")	54
SL46275TET	Reverse Lip Tapered 25mm x 250mm (1" x 9 ¾")	47	SL46520AT	Radiolucent Concave Micro 14mm x 60mm (½" x 2 ¾")	54
SL46280RET	Reverse Lip Rigid 32mm x 100mm (1 ¼" x 4")	46	SL46520ATT	Radiolucent Concave Micro with Teeth 14mm x 60mm (½" x 2 ¾")	54
SL46281RET	Reverse Lip Rigid 32mm x 110mm (1 ¼" x 4 ¾")	46	SL46523AT	Radiolucent Concave Micro 14mm x 65mm (½" x 2 ½")	54
SL46282RET	Reverse Lip Rigid 32mm x 120mm (1 ¼" x 4 ¾")	46	SL46523ATT	Radiolucent Concave Micro with Teeth 14mm x 65mm (½" x 2 ½")	54
SL46283RET	Reverse Lip Rigid 32mm x 130mm (1 ¼" x 5 ½")	46	SL46524AT	Radiolucent Concave Micro 14mm x 70mm (½" x 2 ¾")	54
SL46284RET	Reverse Lip Rigid 32mm x 140mm (1 ¼" x 5 ½")	46	SL46524ATT	Radiolucent Concave Micro with Teeth 14mm x 70mm (½" x 2 ¾")	54
SL46285RET	Reverse Lip Rigid 32mm x 150mm (1 ¼" x 5 ¾")	46	SL46526AT	Radiolucent Concave Micro 14mm x 75mm (½" x 3")	54
SL46286RET	Reverse Lip Rigid 32mm x 160mm (1 ¼" x 6 ¼")	46	SL46526ATT	Radiolucent Concave Micro with Teeth 14mm x 75mm (½" x 3")	54
SL46287RET	Reverse Lip Rigid 32mm x 170mm (1 ¼" x 6 ¾")	46	SL46528AT	Radiolucent Concave Micro 14mm x 80mm (½" x 3 ½")	54
SL46288RET	Reverse Lip Rigid 32mm x 180mm (1 ¼" x 7")	46	SL46528ATT	Radiolucent Concave Micro with Teeth 14mm x 80mm (½" x 3 ½")	54
SL46289RET	Reverse Lip Rigid 32mm x 190mm (1 ¼" x 7 ½")	47	SL46558	Radiolucent Concave 25mm x 80mm (1" x 3 ¾")	49
SL46290RET	Reverse Lip Rigid 32mm x 200mm (1 ¼" x 7 ¾")	47	SL46559	Radiolucent Concave 25mm x 90mm (1" x 3 ½")	49
SL46291RET	Reverse Lip Rigid 32mm x 210mm (1 ¼" x 8 ¼")	47	SL46560	Radiolucent Concave 25mm x 100mm (1" x 4")	49
SL46292RET	Reverse Lip Rigid 32mm x 220mm (1 ¼" x 8 ¾")	47	SL46561	Radiolucent Concave 25mm x 110mm (1" x 4 ¾")	49
SL46293RET	Reverse Lip Rigid 32mm x 230mm (1 ¼" x 9")	47	SL46562	Radiolucent Concave 25mm x 120mm (1" x 4 ¾")	49
SL46294RET	Reverse Lip Rigid 32mm x 240mm (1 ¼" x 9 ½")	47	SL46563	Radiolucent Concave 25mm x 130mm (1" x 5 ½")	49
SL46295RET	Reverse Lip Rigid 32mm x 250mm (1 ¼" x 9 ¾")	47	SL46564	Radiolucent Concave 25mm x 140mm (1" x 5 ½")	49
SL46312AT	Radiolucent Concave Lateral 23mm x 35mm (1" x 1 ¾")	50	SL46565	Radiolucent Concave 25mm x 150mm (1" x 5 ¾")	49
SL46312ATT	Radiolucent Concave Lateral with Teeth 23mm x 35mm (1" x 1 ¾")	51	SL46566	Radiolucent Concave 25mm x 160mm (1" x 6 ¼")	49
SL46314AT	Radiolucent Concave Lateral 23mm x 40mm (1" x 1 ½")	50	SL46567	Radiolucent Concave 25mm x 170mm (1" x 6 ¾")	49
SL46314ATT	Radiolucent Concave Lateral with Teeth 23mm x 40mm (1" x 1 ½")	51	SL46568	Radiolucent Concave 25mm x 180mm (1" x 7")	49
SL46316AT	Radiolucent Concave Lateral 23mm x 45mm (1" x 1 ¾")	50	SL46569	Radiolucent Concave 25mm x 190mm (1" x 7 ½")	49
SL46316ATT	Radiolucent Concave Lateral with Teeth 23mm x 45mm (1" x 1 ¾")	51	SL46570	Radiolucent Concave 25mm x 200mm (1" x 7 ¾")	49
SL46318AT	Radiolucent Concave Lateral 23mm x 50mm (1" x 2")	50	SL46571	Radiolucent Concave 25mm x 210mm (1" x 8 ¼")	49
SL46318ATT	Radiolucent Concave Lateral with Teeth 23mm x 50mm (1" x 2")	51	SL46572	Radiolucent Concave 25mm x 220mm (1" x 8 ¾")	49
SL46320AT	Radiolucent Concave Lateral 23mm x 55mm (1" x 2 ½")	50	SL46573	Radiolucent Concave 25mm x 230mm (1" x 9")	49
SL46320ATT	Radiolucent Concave Lateral with Teeth 23mm x 55mm (1" x 2 ½")	51	SL46574	Radiolucent Concave 25mm x 240mm (1" x 9 ½")	49
SL46322AT	Radiolucent Concave Lateral 23mm x 60mm (1" x 2 ¾")	50	SL46575	Radiolucent Concave 25mm x 250mm (1" x 9 ¾")	49
SL46322ATT	Radiolucent Concave Lateral with Teeth 23mm x 60mm (1" x 2 ¾")	51	SL46604ATT	PLA Medial/Lateral with Fine Teeth 13½mm x 40mm (½" x 1 ½")	70
SL46324AT	Radiolucent Concave Lateral 23mm x 65mm (1" x 2 ½")	50	SL46605ATT	PLA Medial/Lateral with Fine Teeth 13½mm x 50mm (½" x 2")	70
SL46324ATT	Radiolucent Concave Lateral with Teeth 23mm x 65mm (1" x 2 ½")	51	SL46606ATT	PLA Medial/Lateral with Fine Teeth 13½mm x 60mm (½" x 2 ¾")	70
SL46326AT	Radiolucent Concave Lateral 23mm x 70mm (1" x 2 ¾")	50	SL46607ATT	PLA Medial/Lateral with Fine Teeth 13½mm x 70mm (½" x 2 ¾")	70
SL46326ATT	Radiolucent Concave Lateral with Teeth 23mm x 70mm (1" x 2 ¾")	51	SL46608ATT	PLA Medial/Lateral with Fine Teeth 13½mm x 80mm (½" x 3 ½")	70
SL46328AT	Radiolucent Concave Lateral 23mm x 75mm (1" x 3")	50	SL46609ATT	PLA Medial/Lateral with Fine Teeth 13½mm x 90mm (½" x 3 ½")	70
SL46328ATT	Radiolucent Concave Lateral with Teeth 23mm x 75mm (1" x 3")	51	SL46610ATT	PLA Medial/Lateral with Fine Teeth 13½mm x 100mm (½" x 4")	70
SL46330AT	Radiolucent Concave Lateral 23mm x 80mm (1" x 3 ¼")	50	SL46611ATT	PLA Medial/Lateral with Fine Teeth 13½mm x 110mm (½" x 4 ¾")	70
SL46330ATT	Radiolucent Concave Lateral with Teeth 23mm x 80mm (1" x 3 ¼")	51	SL46612ATT	PLA Medial/Lateral with Fine Teeth 13½mm x 120mm (½" x 4 ¾")	70
SL46332AT	Radiolucent Concave Lateral 23mm x 85mm (1" x 3 ¾")	50	SL46660	Double Angled 25mm x 76mm (1" x 3")	41
SL46332ATT	Radiolucent Concave Lateral with Teeth 23mm x 85mm (1" x 3 ¾")	51	SL46661	Double Angled 25mm x 64mm (1" x 2 ½")	41
SL46334AT	Radiolucent Concave Lateral 23mm x 90mm (1" x 3 ½")	50	SL46662	Double Angled 13mm x 76mm (½" x 3")	41
SL46334ATT	Radiolucent Concave Lateral with Teeth 23mm x 90mm (1" x 3 ½")	51	SL46663	Double Angled 13mm x 70mm (½" x 2 ¾")	41
SL46336AT	Radiolucent Concave Lateral 23mm x 95mm (1" x 3 ¾")	50	SL46664	Double Angled Bulb 25mm x64mm (1" x 2½")	41
SL46336ATT	Radiolucent Concave Lateral with Teeth 23mm x 95mm (1" x 3 ¾")	51	SL46704AT	PLA Medial/Lateral with Lip 13½mm x 40mm (½" x 1 ½")	70
SL46338AT	Radiolucent Concave Lateral 23mm x 100mm (1" x 4")	50	SL46705AT	PLA Medial/Lateral with Lip 13½mm x 50mm (½" x 2")	70
SL46338ATT	Radiolucent Concave Lateral with Teeth 23mm x 100mm (1" x 4")	51	SL46706AT	PLA Medial/Lateral with Lip 13½mm x 60mm (½" x 2 ¾")	70
SL46340AT	Radiolucent Concave Lateral 23mm x 105mm (1" x 4 ¼")	50	SL46707AT	PLA Medial/Lateral with Lip 13½mm x 70mm (½" x 2 ¾")	70
SL46340ATT	Radiolucent Concave Lateral with Teeth 23mm x 105mm (1" x 4 ¼")	51	SL46708AT	PLA Medial/Lateral with Lip 13½mm x 80mm (½" x 3 ½")	70
SL46342AT	Radiolucent Concave Lateral 23mm x 110mm (1" x 4 ¾")	50	SL46709AT	PLA Medial/Lateral with Lip 13½mm x 90mm (½" x 3 ½")	70
SL46342ATT	Radiolucent Concave Lateral with Teeth 23mm x 110mm (1" x 4 ¾")	51	SL46710AT	PLA Medial/Lateral with Lip 13½mm x 100mm (½" x 4")	70
SL46344AT	Radiolucent Concave Lateral 23mm x 115mm (1" x 4 ½")	50	SL46711AT	PLA Medial/Lateral with Lip 13½mm x 110mm (½" x 4 ¾")	70
SL46344ATT	Radiolucent Concave Lateral with Teeth 23mm x 115mm (1" x 4 ½")	51	SL46712AT	PLA Medial/Lateral with Lip 13½mm x 120mm (½" x 4 ¾")	70
SL46346AT	Radiolucent Concave Lateral 23mm x 120mm (1" x 4 ¾")	50	SL46982AC	Bone Lever 14mm x 40mm (½" x 1 ½")	72
SL46346ATT	Radiolucent Concave Lateral with Teeth 23mm x 120mm (1" x 4 ¾")	51	SL46983AC	Bone Lever 14mm x 50mm (½" x 2")	72
SL46410T	Hibbs with Teeth 25mm x 75mm (1" x 3")	66	SL46984AC	Bone Lever 14mm x 60mm (½" x 2 ¾")	72
SL46411T	Hibbs with Teeth Long 25mm x 75mm (1" x 3")	66	SL46985AC	Bone Lever 14mm x 70mm (½" x 2 ¾")	72
SL46414T	Hibbs with Teeth 25mm x 115mm (1" x 4 ½")	66	SL46986AC	Bone Lever 14mm x 80mm (½" x 3 ½")	72
SL46418T	Hibbs with Teeth 25mm x 165mm (1" x 6 ½")	66	SO45612A	Meyerding 51mm x 60mm (2" x 2 ¾")	57
SL46510AT	Radiolucent Concave Micro 14mm x 35mm (½" x 1 ¾")	54	SO45612AET	Radiolucent Meyerding 51mm x 60mm (2" x 2 ¾")	57

PART #	DESCRIPTION	PAGE	PART #	DESCRIPTION	PAGE
SO45612AETT	Radiolucent Meyerding with Teeth 51mm x 60mm (2" x 2 3/8")	57	SO46157	Mini-Open Richardson 51mm x 127mm (2" x 5")	33
SO45612AT	Meyerding with Teeth 51mm x 60mm (2" x 2 3/8")	57	SO46158	Mini-Open Richardson 51mm x 152mm (2" x 6")	33
SO45613A	Meyerding 51mm x 75mm (2" x 3")	57	SO46160ET	Radiolucent Harrington 64mm x 152mm (2 1/2" x 6")	27
SO45613AT	Meyerding with Teeth 51mm x 75mm (2" x 3")	57	SO46168B	Non-Slip Balfour L 77mm x 66mm (3" x 2 1/2")	26
SO45615A	Meyerding 51mm x 100mm (2" x 4")	57	SO46168C	Non-Slip Balfour L 77mm x 94mm (3" x 3 3/4")	26
SO45615AT	Meyerding with Teeth 51mm x 100mm (2" x 4")	57	SO46169B	Non-Slip Balfour R 77mm x 66mm (3" x 2 1/2")	26
SO45616A	Meyerding 51mm x 125mm (2" x 5")	57	SO46169C	Non-Slip Balfour R 77mm x 94mm (3" x 3 3/4")	26
SO45616AT	Meyerding with Teeth 51mm x 125mm (2" x 5")	57	SO46175	Deaver 64mm x 127mm (2 1/2" x 5")	36
SO45670ATTC	Radiolucent Posterior Lumbar 23mm x 80mm (1" x 3 1/8")	58	SO46175ET	Radiolucent Deaver 64mm x 127mm (2 1/2" x 5")	36
SO45671ATTC	Radiolucent Posterior Lumbar 45mm x 80mm (1 3/4" x 3 1/8")	58	SO46176	Deaver 76mm x 127mm (3" x 5")	36
SO45672ATTC	Radiolucent Posterior Lumbar 60mm x 80mm (2 3/8" x 3 1/8")	58	SO46176ET	Radiolucent Deaver 76mm x 127mm (3" x 5")	36
SO45673ATTC	Radiolucent Posterior Lumbar 23mm x 100mm (1" x 4")	58	SO46178	Deaver 102mm x 127mm (4" x 5")	36
SO45674ATTC	Radiolucent Posterior Lumbar 45mm x 100mm (1 3/4" x 4")	58	SO46178ET	Radiolucent Deaver 100mm x 127mm (4" x 5")	36
SO45675ATTC	Radiolucent Posterior Lumbar 60mm x 100mm (2 3/8" x 4")	58	SO46187A	St Marks [no lip] 64mm x 178mm (2 1/2" x 7")	35
SO45676ATTC	Radiolucent Posterior Lumbar 23mm x 120mm (1" x 4 3/4")	58	SO46188	St Marks [with lip] 64mm x 178mm (2 1/2" x 7")	35
SO45677ATTC	Radiolucent Posterior Lumbar 45mm x 120mm (1 3/4" x 4 3/4")	58	SO46188ET	Radiolucent St Marks [with lip] 64mm x 178mm (2 1/2" x 7")	35
SO45678ATTC	Radiolucent Posterior Lumbar 60mm x 120mm (2 3/8" x 4 3/4")	58	SO46189	St Marks [with lip] 64mm x 140mm (2 1/2" x 5 1/2")	35
SO45679ATTC	Radiolucent Posterior Lumbar 23mm x 150mm (1" x 5 7/8")	58	SO46189ET	Radiolucent St Marks [with lip] 64mm x 140mm (2 1/2" x 5 1/2")	35
SO45680ATTC	Radiolucent Posterior Lumbar 45mm x 150mm (1 3/4" x 5 7/8")	58	SO46241A	Wide Balfour [no lips] 83mm x 57mm (3 1/4" x 2 1/4")	26
SO45681ATTC	Radiolucent Posterior Lumbar 60mm x 150mm (2 3/8" x 5 7/8")	58	SO46241AEB	Radiolucent Wide Balfour [no lips] 83mm x 57mm (3 1/4" x 2 1/4") [Anodized]	26
SO45682ATTC	Radiolucent Posterior Lumbar 23mm x 180mm (1" x 7")	58	SO46241B	Wide Balfour [no lips] 100mm x 61mm (4" x 2 3/8")	26
SO45683ATTC	Radiolucent Posterior Lumbar 45mm x 180mm (1 3/4" x 7")	58	SO46241BEB	Radiolucent Wide Balfour [no lips] 102mm x 57mm (4" x 2 1/4") [Anodized]	26
SO45684ATTC	Radiolucent Posterior Lumbar 60mm x 180mm (2 3/8" x 7")	58	SO46241C	Wide Balfour [no lips] 114mm x 63mm (4 1/2" x 2 1/2")	26
SO46103	Saldanha Archd 76mm x 140mm (3" x 5 1/2")	73	SO46241CEB	Radiolucent Wide Balfour [no lips] 114mm x 57mm (4 1/2" x 2 1/4") [Anodized]	26
SO46124A	Malleable 76mm x 254mm (3" x 10")	30	SO46400	Tummy Tuck 76mm x 140mm (3" x 5 1/2")	73
SO46124AET	Mini-Open Malleable 76mm x 254mm (3" x 10")	31	SO46401	Tummy Tuck 102mm x 140mm (4" x 5 1/2")	73
SO46124AV	Malleable Slotted 76mm x 203mm (3" x 8")	32	SO46402	Cat's Paw 6mm x 38mm (1/4" x 1 1/2")	73
SO46124B	Malleable 102mm x 254mm (4" x 10")	30	SO48001	Towel Clip Hook 45mm x 64mm (1 3/4" x 2 1/2")	73
SO46124BET	Radiolucent Malleable 102mm x 254mm (4" x 10")	31			
SO46124C	Aortic Aneurysm Large 457mm x 229mm (18" x 9")	40			
SO46124D	Aortic Aneurysm Small 305mm x 216mm (12" x 8 1/2")	40			
SO46129	Kelly 64mm x 76mm (2 1/2" x 3")	28			
SO46129ET	Radiolucent Kelly 64mm x 76mm (2 1/2" x 3")	28			
SO46130	Kelly 76mm x 89mm (3" x 3 1/2")	28			
SO46130ET	Radiolucent Kelly 76mm x 89mm (3" x 3 1/2")	28			
SO46134	Fence 102mm x 102mm (4" x 4")	29			
SO46134ET	Radiolucent Fence 102mm x 102mm (4" x 4")	29			
SO46135	Fence 102mm x 127mm (4" x 5")	29			
SO46135ET	Radiolucent Fence 102mm x 127mm (4" x 5")	29			
SO46136	Fence 102mm x 152mm (4" x 6")	29			
SO46136A	Fence 102mm x 203mm (4" x 8")	29			
SO46136AET	Radiolucent Fence 102mm x 203mm (4" x 8")	29			
SO46136ET	Radiolucent Fence 102mm x 152mm (4" x 6")	29			
SO46140	Balfour [with lips] 70mm x 73mm (2 3/4" x 3")	25			
SO46140C	Balfour [small lips] 70mm x 73mm (2 3/4" x 3")	25			
SO46140CEB	Radiolucent Balfour [small lips] 65mm x 72mm (2 1/2" x 2 3/4") [Anodized]	25			
SO46140CET	Radiolucent Balfour [small lips] 65mm x 72mm (2 1/2" x 2 3/4")	25			
SO46140EB	Radiolucent Balfour [with lips] 65mm x 72mm (2 1/2" x 2 3/4") [Anodized]	25			
SO46140ET	Radiolucent Balfour [with lips] 65mm x 72mm (2 1/2" x 2 3/4")	25			
SO46141	Balfour [with lips] 83mm x 73mm (3 1/4" x 3")	25			
SO46141C	Balfour [small lips] 83mm x 73mm (3 1/4" x 3")	25			
SO46141CEB	Radiolucent Balfour [small lips] 83mm x 70mm (3 1/4" x 2 3/4") [Anodized]	25			
SO46141CET	Radiolucent Balfour [small lips] 83mm x 70mm (3 1/4" x 2 3/4")	25			
SO46141D	Balfour [with lips] 83mm x 123mm (3 1/4" x 5")	25			
SO46141DEB	Radiolucent Balfour [with lips] 83mm x 127mm (3 1/4" x 5") [Anodized]	25			
SO46141DET	Radiolucent Balfour [with lips] 83mm x 127mm (3 1/4" x 5")	25			
SO46141EB	Radiolucent Balfour [with lips] 83mm x 70mm (3 1/4" x 2 3/4") [Anodized]	25			
SO46141ET	Radiolucent Balfour [with lips] 83mm x 70mm (3 1/4" x 2 3/4")	25			
SO46141F	Balfour [with lips] 83mm x 175mm (3 1/4" x 6 7/8")	25			
SO46151	Richardson 51mm x 127mm (2" x 5")	34			
SO46151ET	Radiolucent Richardson 51mm x 127mm (2" x 5")	34			
SO46152	Richardson 51mm x 152mm (2" x 6")	34			
SO46152ET	Radiolucent Richardson 51mm x 152mm (2" x 6")	34			
SO46153	Richardson 51mm x 178mm (2" x 7")	34			
SO46153ET	Radiolucent Richardson 51mm x 178mm (2" x 7")	34			
SO46155	Mini-Open Richardson 51mm x 76mm (2" x 3")	33			
SO46156	Mini-Open Richardson 51mm x 102mm (2" x 4")	33			


CURRENT SYSTEMS REFERENCE:

SL80043A	Accessory Blade Kit (for Body Contouring System)
SL82001	Bariatric System
SL82002	Colon / Rectal System
SL82003	General / Vascular System
SL82004	HIPEC System
SL82005	Kidney System
SL82006	Pediatric System
SL82007	Liver / Oncology System
SL82008	Radiolucent Liver / Oncology System
SL82009	Small Incision System
SL82010	GYN System
SL82011	GYN Oncology System
SL82012	Urology System
SL82013	Urologic Oncology System
SL82014	Hip Tether System
SL82015	Hip System
SL82016	Body Contouring System
SL82017	Cervical System
SL82018	CCA Retractor System
SL82019	Anterior Lumbar System
SL82020	Posterior Lumbar System
SL82021	MIS PLA® System
SL82022	Bilateral Expansion Kit (for Body Contouring System)
SL82023	Anterior Lumbar Ring System




ELITE RAIL CLAMPS

Measured from bottom of the bed rail to the highest point joint can reach.
(NOTE: Elite III is measured with jaw in the FULLY closed position.)




RAIL EXTENDERS

Length of rail space provided, from end to end.




UNILATERAL CROSSBARS + ANGLED ARMS

Total length equals the sum of both sides. (A + B)



BILATERAL CROSSBARS


Length of usable rod on each section (A x B x C).



LOWER ABDOMINAL BARS


All measure 17" across, between joints. Lateral arm sections are measured separately.

Metric	Approx. English
6mm	¼"
10mm	⅜"
13mm	½"
23mm	⅞"
25mm	1"
30mm	1 ⅛"
32mm	1 ¼"
35mm	1 ⅜"
40mm	1 ½"
45mm	1 ¾"
50mm	2"
55mm	2 ⅛"
60mm	2 ⅜"
65mm	2 ½"
70mm	2 ¾"
75mm	3"
80mm	3 ⅛"
85mm	3 ⅜"
90mm	3 ½"
95mm	3 ¾"
100mm	4"
110mm	4 ⅜"
120mm	4 ¾"
130mm	5 ⅛"
140mm	5 ½"
150mm	5 ⅞"
160mm	6 ¼"
170mm	6 ⅝"
180mm	7"
190mm	7 ½"
200mm	7 ⅞"
210mm	8 ¼"
220mm	8 ⅝"
230mm	9"
240mm	9 ½"
250mm	9 ⅞"
260mm	10 ¼"
270mm	10 ⅝"
280mm	11"
290mm	11 ⅜"
300mm	11 ¾"




HANDLES

Center of nipple hole to end of handle, not including end cap.




MALLEABLE BLADES

Width (A) x total length of blade before bending (B).




BLADES (NON-MALLEABLE)

Width (A) x height (B) from underside of blade at top of bend to distal tip.




NATHANSON HOOKS

Handle length, not including nipple for SL versions (A) x hook height (B) x tip length (C).



BOLLING RETRACTOR BLADES

Total handle length including retractor blade.



HIP BLADES

Length from end to end, not including bend.

Service and Warranty

PRODUCT WARRANTY: Thompson Surgical Instruments, Inc. warrants that all instruments are free from defects in material or workmanship for 10 years *. This warranty is void if failure of the product has resulted from normal wear and tear from instrument use, accident, abuse, misapplication, negligence, or if the product has been damaged, altered or repaired outside of Thompson Surgical's facility. Warranty void if purchased from a non-authorized supplier/distributor.

- * Malleable and Cobra blades carry a two (2) year guarantee against defects in both material and workmanship. Please replace Cobra blades after two (2) years of surgical use or 200 uses if this occurs first to avoid possible fracturing or cracking of Cobra blade tips.
- * Hip Tethers carry a one (1) year guarantee against defects in both material and workmanship. Please replace Hip Tethers after 10 reprocessing cycles to avoid possible breakage.
- * Fiberoptic Xenon lighting products carry a 60 day guarantee. Fiberoptic LED lighting products carry a one year guarantee. Fibers are not warranted against breakage. Please be cautious when handling and avoid stepping on or kicking the fiber optic cable to prevent damage.
- * Blade finishes do not carry a warranty.
- * Please note that aluminum blades are sensitive to high and low pH levels. Please follow cleaning and sterilization instructions exactly. (See Thompson Retractor IFU.) Failure to do so may cause aluminum blade staining, pitting, or other corrosion and will void blade warranty.

Visit thompsonsurgical.com for complete Service and Warranty information.

THIS PAGE INTENTIONALLY LEFT BLANK

For a Free Trial Call Today*

1.800.227.7543

Rev C
062818
mg prodcut



Uncompromised Exposure

VISIT US ONLINE: THOMPSONSURGICAL.COM

Thompson
SURGICAL INSTRUMENTS

* Free trial valid for U.S. customers only. International customers please call for more information.

© 2018 Thompson Surgical Instruments, Inc. Traverse City, Michigan. Printed in U.S.A.

® S-Lock®, PLA®, and the "T Circle" logomark  are Registered Trademarks of Thompson Surgical Instruments, Inc.

Patents: US4971038, US5025780, US5888197, US5897087, US5902233, US5984865, US6033363, US6416465, US6511423, US7338442, US7749163, US8257255, US8360971, US8617064

Other patents pending.

The One Step Dilator (01-01004, 01-02015, 01-01006) is manufactured by Thompson MIS and distributed by Thompson Surgical Instruments. The MIS Retractor Lite LED (#40002MIS), Mini Retractor Lite LED (#40002MRL), and LED Light Source (#LLS-2000-2) are manufactured by Sunoptic Technologies® and distributed by Thompson Surgical Instruments.

Made In
USA